**REFRACTÓMETROS DE MANO
MULTISCALE REFRACTOMETER
REFRACTOMÈTRES À MAIN**

Zuzi

Serie 300/Series 300/Série 300

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Thank you for choosing this instrument. We sincerely wish that you enjoy your Zuzi hand refractometer. We highly recommend looking after this equipment according to what is stated in this manual.

Zuzi develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

Please bear in mind the following:

- ◆ This manual is inseparable from the refractometer, so it should be available for all the users of this equipment.
- ◆ You should carefully handle the refractometer avoiding sudden movements, knocks, free fall of heavy / sharp objects on it.
- ◆ Never dismantle the different pieces of the refractometer to repair it yourself, since it could produce a defective use of the instrument
- ◆ If you have any doubt about operation do not hesitate in contacting your wholesaler.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage (only Universal Refractometer)
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.

**ATTENTION!!
IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.**

INDEX OF LANGUAGES

1. Spanish	2-16
2. English.....	17-31
3. French	32-46

INDEX OF CONTENTS

1. USES OF THE INSTRUMENT.....	18
2. DESCRIPTION.....	19
3. TECHNICAL SPECIFICATIONS.....	19
4. THEORETICAL FOUNDATIONS.....	27
5. OPERATION MODE.....	29
6. MAINTENANCE AND CLEANING.....	31
7. ACCESSORIES.....	31

ANNEX I: °BRIX /T^a CORRECTION CHART.....47

ANNEX II: UNIT CONVERSION CHART49

ANNEX: CE CERTIFICATE.....52

1. USES OF THE INSTRUMENT

ZUZI hand refractometers are very handy instruments that allows measuring the concentration of dissolved solids in aqueous solutions. They can be used for measurements on a broad range of samples such as: sugar concentration of fruits, concentrated juices, condensed milk, carbonated beverages, grape wine, cutting oils, jam, egg yolk, proteins in serum, specific weight of urine, salt water, and so on.

Thanks to its small size and easy handling, this instrument can be used by any person (without previous experience), in any place, both in laboratories and outside, for fieldworks or routine inspections, as for example:

- ◆ Fruit producers: to know sugar concentration and determine its grade of ripeness. They are also suitable to find out the content of sugar in sugar plants.
- ◆ Food industries (sugar refineries, canning industries, distilleries, other beverage factories, honey producers, and so on): to carry out a standardization that makes possible to obtain mixtures with exact proportions in order to maintain the quality of the product and a strict quality control.
- ◆ Medical and veterinary applications (clinical laboratories, hospitals, veterinary clinics): to study the content of total proteins in blood, serum, plasma and urine. They are also used to determine the specific weight of urine, as well as to know the content of sugar in medicines such as syrups.


◆ Chemical industries: related to glass, perfumes, ... generally all kind of chemical factories, research institutes and basic research laboratories.

◆ Oceanographic research (salinity refractometers): control of natural and artificial seawater, the cultivation of seafish and the investigation of fishing grounds.

2. DESCRIPTION

Figure 1. General Refractometer

1.1 Prism

1.2 Prism cover plate


1.3 Scale adjustment screw (Calibration)

1.4 Rubber-covered refractometer body

1.5 Diopter adjustment ring

1.6 Eyepiece

Figure 2. Universal refractometer

1.7 Light inlet

1.8 Scale selector

1.9 Polarization ring

1.10 Cover plate

1.11 Diopter adjustment ring / Eyepiece

3. TECHNICAL SPECIFICATIONS

The following information will help you to choose the most adequate refractometer depending on the application required. The approximate concentration of different types of samples are related to the different models of refractometers.

Brix scale refractometers


Model 50301010

It is used for samples with low concentration of dissolved solids (less than 10%), suitable for tomatoes, low concentrated juices and industrial cutting oils. Calibration is made adjusting to zero with distilled water.

Scale: 0 - 10% Brix

Accuracy: 0.1 %

Size: 27x40x210 mm

Model 50301110

It presents the same specifications as the previous model but with Automatic Temperature Compensation (10-30° C).


Model 50301020

It is suitable for most of fruits except for grape juice; calibration is made adjusting to zero with distilled water.

Scale: 0 - 20% Brix

Accuracy: 0.1 %

Size: 30x40x200 mm


Models 50301030/50301032 (ECO Series)

It is commonly used for low concentrations in fruits, fruit juices, beverages, dairy products, industrial oils and so on. Calibration is made adjusting to zero with distilled water.

Scale: 0 - 32% Brix

Accuracy: 0.2 %

Size: 30x40x200 mm

Models 50301130/50301132 (ECO Series)

They present the same specifications as the previous models but with Automatic Temperature Compensation (10-30° C).


Model 50301040

It is used for mid-range-concentration samples such as concentrated fruit juices, jams, canned products with sugar solutions, etc. Calibration is made with a sugar solution in distilled water at 28% (28 g of sugar in 100 mL of solution).

Scale: 28 - 62% Brix

Accuracy: 0.2 %

Size: 27x40x140 mm


Model 50301050

It is used for high concentrations such as condensed milk, liquid caramel, very concentrated juices, jams, etc. Calibration is made with a sugar solution in distilled water at 45% (28 g of sugar in 100 mL of solution).

Scale: 45 - 82% Brix

Accuracy: 0.5 %

Size: 27x40x140 mm

Model 50301060

It is specially developed for measuring the three main indexes of honey: sugar concentration, Baume degrees and water content.

Calibration is made with Monobromonaphthalene and the test piece supplied with the instrument; as follows: Place 1 or 2 drops of oil on the bright surface of the test piece. Open the cover plate, place the glass on the surface of the prism and close the cover plate. Keep it pressed so as the test piece does not slide down; turn the adjustment screw until the white/blue boundary line fits the reference line. The test piece is not used for normal measurements. Just put 1 or 2 drops of honey on the prism, close the cover plate and looking through the eyepiece read the result on the scale.

Scale: 58 - 90% Brix

38 - 43 ° Baume


17 - 27% Water

Accuracy: 0.5 % Brix

0.5 ° Baume

0.5% Water

Size: 27x40x160 mm


Model 50301080

Thanks to its wide scale it is used for measuring all kind of products such as fruits, juices, alcohol-free beverages and even industrial cutting oils. Calibration is made adjusting to zero with distilled water.

Scale: 0 - 80% Brix

Accuracy: 1 %

Size: 27x40x160 mm

Salinity refractometers

Model 50302010


It is specially used for controlling salt concentration in natural or artificial seawater, sea products, etc. It is also provided with specific weight scale. Calibration is made with distilled water.

Scale: 0 - 10% Sal
1000-1070 g/mL

Accuracy: 0.1 %

0.2 g/mL

Size: 27x40x190 mm

Model 50302110

It presents the same specifications as the previous model but with Automatic Temperature Compensation (10-30° C).


Modelo 50302020/50302022 (ECO Series)

It is used to measure salt concentration in water, oceanography, aquariums, fishing grounds, etc. Calibration is made with distilled water.

Scale: 0 - 28% Sal

Accuracy: 0.2 %

Size: 27x40x160 mm

Modelo 50302122 (ECO Series)

It presents the same specifications as the previous model but with Automatic Temperature Compensation (10-30° C).

Clinical refractometers


Model 50303020/50303022 (ECO Series)

This model is provided with three different scales for serum proteins, urine specific gravity and refractive index. The calibration is made with distilled water and the boundary line must fit with the value 1.333 of the refractive index scale.

Scale: 0 - 12 g/dl Serum proteins
1.000 - 1.040 Specific gravity
1.333 - 1.360 nD

Accuracy: 0.2 %
0.002
0.005

Size: 27x40x160 mm

Models 50303120/50303122 (ECO Series)

It presents the same specifications as the previous model but with Automatic Temperature Compensation (10-30° C).

Alcohol and grape juice refractometers


Model 50305010

It is used for measuring the degree of alcohol in an aqueous solution, it can be used for wine and liquors taking into account other ingredients. Calibration is made with distilled water.

Scale: 0 - 80% Alcohol

Accuracy: 1 %

Size: 27x40x160 mm


Model 50305052 (ECO Series)

This instrument measures the probable degree of alcohol in grape juice, as well as Baumé and Brix degrees, thus allowing to estimate the ripeness point of grapes before and during harvest. Calibration is made with distilled water.

Scale: 0 - 25% Alcohol
0-22°Baumé
0-40% Brix

Accuracy: 0.2 % Alcohol
0.2° Baumé
0.2% Brix

Size: 27x40x160 mm

Models 50305150/50305152 (ECO Series)

It presents the same specifications as the previous model but with Automatic Temperature Compensation (10-30° C).

Model 50305162 (ECO Series)

It is used to measure the alcohol degree with 3 different scales, commonly used in different parts of the world. Calibration is made with distilled water and it is provided with Automatic Temperature Compensation (10-30° C).

Scale: 0-140° Oe
0-25° KMW Babo
0-32% mas sacch

Accuracy: 1° Oe
0.2° KMW Babo
0.2% mas sacch

Size: 27x40x160 mm

Milk refractometer

Model 50308000

It is used to measure the water content in milk. Calibration is made with distilled water.

Scale: 0-20%

Accuracy: 0.2%

Size: 27x40x160 mm

Batteries and anti-freezing refractometers


Model 50304010

It is used to measure the freezing point of cooling systems based on ethylene and propylene glycol, as well as for checking the strength of the electrolytic solutions of batteries. Calibration is made with distilled water.

Scale: -60 a 32°F (ethylene glycol)
-50 a 32 °F (propylene glycol)
(anti-freezing)
1.15 - 1.30 s (Battery)

Accuracy: 10 °F
0.01 s

Size: 27x40x160 mm

Model 50304030/50304032 (ECO Series)

It is used to measure the freezing point of cooling systems based on ethylene and propylene glycol and of cleaning liquid, as well as for checking the strength of the electrolytic solutions of batteries. It is also appropriate to determine the urea content in fuel deposits by using reducing agents (AdBlue, EcoBlue). Calibration is made with distilled water.


Scale: -50 a 0 °C
(ethylene / propylene glycol)
(Anti-freezing)
1.15 - 1.30 s (Battery)
-40 a 0 °C (Cleaning)

Accuracy: 5 °C

0.01 s

5°C

Size: 27x40x160 mm

Model 50304132 (ECO Series)

It presents the same specifications as the previous model but with Automatic Temperature Compensation (10-30° C).

Emulsion refractometers


Model 50308110

It is used to measure the concentration in weight percentage of M10 and MDT emulsions, the most commonly used emulsions in machinery. It is provided with Automatic Temperature Compensation (10-30° C). Calibration is made with distilled water.

Scale: 0-15% M-10/MDT

Accuracy: 0.5%

Size: 27x40x160 mm

Universal refractometers (3 scales)

They are very useful to measure the dissolved solids content (%Brix) or refractive index (nD) in aqueous solutions.

Thanks to its wide measuring range divided into 3 scales, they allow measuring a wide range of products.

Model 50308090

Total scale: 1.333 nD-1.515 nD

Division scales: 1.333-1.400 nD
1.400-1.470 nD
1.470-1520 nD

Accuracy: 0.0005 nD

Size: 37x34x200 mm

Model 50301090

Total range: 0-90 %Brix

Division scales: 0-42 %
42-71%
71-90%

Accuracy: 0.2%

Size: 37x34x200 mm

SAMPLE		CONCENTRATION							
		10	20	30	40	50	60	70	80
FRUITS, FRUIT JUICES SOFT DRINKS	orange, pear	■	■						
	tomato	■							
	apple, melon		■	■					
	strawberry, peach	■	■						
	grape		■	■					
	concentrated fruit juices					■	■	■	■
	fruit juices		■	■					
	tomato puree		■	■					
	tomato juice	■	■						
	carbonated beverages	■	■						
	acidophilus beverages		■	■					
	lactic beverages		■	■					
	FOODS	condensed milk						■	■
liquid sugar								■	■
yeast culture solutions		■	■						
grape wine			■	■					
peanut milk		■	■						
canned fruits			■	■	■				
egg yolk							■		
sauce, tomato ketchup					■	■			
milk			■	■					
jam, marmalade							■		
INDUSTRY		cutting/quenching oil	■	■					
	emulsifying solutions	■							

ENGLISH

Automatic Temperature Compensation (ATC) system:

When using models provided with automatic temperature compensation (ATC), it is not necessary to look up the attached charts to correct the measuring according to the working temperature, since this compensation is made automatically by the instrument.


This system makes possible to measure the concentration of a solution without taking into account the temperature.

4. THEORETICAL FOUNDATIONS

Refractometry

When a light beam reaches the surface of separation between two different environments, part of this light is reflected and the rest of it is refracted (it goes in the second environment), propagating with a different angle and speed than in the first environment. These two magnitudes can be characterized by means of the refractive index, which is a characteristic value for each substance.

Figure 3


The refractive index is related to the mass, the load and the number of particles of the substance through which the light radiation is transmitted. There are different kinds of refractometers but most of them are based on the same principle.

An example that explains clearly the phenomenon of refraction is what happens when we introduce a pencil into a glass filled with water: the part of the pencil that is dipped into the water seems to be bent with regards to the other part that is out of the water. This is the phenomenon of light refraction. If we submerge the pencil in hot sugar-saturated water the point of the pencil seems to be even more bent.


Figure 4

The sample of liquid to be analysed is placed on the prism and then it is illuminated with monochromatic or natural light (white light). The light beam goes through the sample and the prism; then in the eyepiece we can observe the prism surface up to find the BOUNDARY LINE (the boundary line is the refracted ray coming from a ray that has incident with an angle of 90° with the normal).

The reading is carried out dividing the eyepiece field in a dark zone and a light zone. The boundary line will determine both zones. All the refracted rays come from incident rays with an incident angle less than 90° .

Figure 5


Once the eyepiece field has been divided, the value of the refractive index can be read directly on the scale.

Before making any measurement, the instrument must be calibrated with distilled water (refractive index 1). For this purpose, a drop of water must be placed on the prism of the refractometer and make the necessary correction by turning the scale adjustment screw.

A refractometer is an instrument that measure the concentration of aqueous solutions by measuring the refractive index. All aqueous solutions can be crossed by light. The path length is proportional to its concentration.

In general, the refractive index of a thicker transparent substance is higher than a less thick liquid, it means, the light speed is smaller in a substance with higher density. Then, if a ray incides in an oblique way on an environment with a higher refractive index, it will be deviated to the normal, whereas if it incides on an environment with a smaller refractive index, it will move away from it. The rays that incide on the direction of the normal are reflected and refracted in the same direction.


Figure 6

5. OPERATION MODE

- ◆ Take out the instrument from its packing and remove all protections
- ◆ Open the prism cover plate (1.2) and clean both the cover plate and the prism (1.1) with alcohol and a fluff-free soft cloth; be very careful not to scratch the prism surface.
- ◆ Direct the prism towards a bright light and adjust the diopter ring (1.5) until the scale is focused.
- ◆ To calibrate the instrument put one or two drops of distilled water or reference sample on the prism (fig. 7.2) and close the plate.
- ◆ By turning the scale adjustment screw (1.3) make the white/blue boundary line match the value 0 in the scale in case of distilled water, or the corresponding value in case of a known concentration reference sample.
- ◆ Open the plate, clean and wipe the prism and the plate, and put one or two drops of sample on the prism; close the plate and looking through the eyepiece, read the concentration result on the scale (fig. 7.3).

Figura 7

◆ To make use easier and avoid Reading mistakes, these instruments present the following features:

- Diopter adjustment ring (1.5), to focus the scale properly.

-High contrast of boundary line (white/blue), what allows an easier and more reliable reading, and at the same time it saves time.

- Rubber hood: it prevents light from entering through the eyepiece during measurements, and increases the luminosity of the scale and the boundary line.

- Rubber grip (1.4): Measuring results are influenced by temperature changes. This rubber grip avoids the user's body heat to be transmitted to the sample.

- Scale adjustment ring (1.3), to calibrate the instrument.

Universal refractometer

Direct the prism towards a bright light and adjust the diopter ring (1.11) until the scale is focused.

◆ By turning the polarization ring (1.9), remove the light chromaticism until a Sharp and focused boundary line between dark and light zones is observed.

◆ Adjust the light inlet (1.7)

◆ Calibration is not necessary since the equipment is calibrated in origin.

◆ Open the cover plate (1.10), clean the prism surface with a soft fluff-free cloth and put one or two drops of the sample to be analyzed. Close the cover, press it slightly and read the result (refractive index or %Brix depending on the model) by looking through the eyepiece.

◆ Select the most adequate scale by turning the scale selector (1.8).

◆ Once measurement has finished, remove the sample with a soft cloth and wipe the prism.

Use the attached charts (annexes) to convert Brix degrees into refractive index value, and to correct the effect of temperature on readings.

6. MAINTENANCE AND CLEANING

◆ Hand refractometers need the same maintenance as any other optical instrument. They must be always clean and they have to be placed in its box to avoid dust settle on them. Take into account the following advices:

◆ Never dismatle the lenses; if there is any dust or dirtiness on the external surface of lenses wipe them with a fluff-free cloth.

◆ To remove dust settled on lenses blow it up with a pear or clean in with a soft brush.

◆ Dry immediately the optical parts in case they get wet.

Never elt the optical parts wet

◆ If they come into contact with acid or basic substances, clean them immediately with water and then dry.

◆ For the cleaning of mechanical pieces use non-corrosive lubricants, being careful not to touch the optical pieces.

◆ Be careful not to knock the instrument or scratch the surface of the prism.

7. ACCESSORIES

Available only for general refractometers:

<i>Code</i>	<i>Description</i>
90300011	Prism cover plate
90300621	BLinker
90300981	Plastic box

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ANNEX II: If your refractometer has not available the automatic compensation temperature use this chart to correct the results you obtain with respect to ambient temperature.

To do it, localize the temperature in which you are working with the sample in the left column, and the most approximate percentage of solution (brix grades) in the upper line, once you have found out the value of correction you will have to add (Temperature higher than 20°C) or subtract (Temperature lower than 20°C) this value of correction to the reading.

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ANNEX III: In this chart you will be able to find the equivalences among density, brix grades, Baumé grades and probable alcohol grades (in grape juice) referring to 20°C.

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INSTRUCTIONS ON ENVIRONMENT PROTECTION

At the end of its life cycle, please, do not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for recycling. It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

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Fdo: Alfonso Ainciburu Sanz
DIRECTOR | GERENTE

A handwritten signature in black ink, appearing to be the name Alfonso Ainciburu Sanz.

BERIAIN a 02 de MARZO de 2010