AUXILAB S.L.

Material de laboratorio  
Laboratory supplies

DECLARACIÓN DE CONFORMIDAD CE  
CE DECLARATION OF CONFORMITY  
DÉCLARATION DE CONFORMITÉ CE

El fabricante | The manufacturer | Le fournisseur:

AUXILAB, S.L.

Declara que el equipo | Declare that the equipment | Declare que l'appareil:

CENTRÍFUGA | CENTRIFUGE | CENTRIFUGEUSES

Código | Code | Code: 52900000

Modelo | Model | Modèle: 2900

Cumple las siguientes directivas | Meet the following directives | Accomplit les directives suivantes:

73/23/CE   Directiva de seguridad eléctrica Directive for electrical safety Directive de sécurité électrique	89/336/CE   Directiva de Compatibilidad electromagnética (CEM) Directive for electromagnetic compatibility (EMC) Directive de compatibilité électromagnétique (CEM)
--------------------------------------------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Cumple las siguientes Normas: | Meet the following Standards | Accomplit les normes suivantes:

EN 61326 | Material eléctrico para medida control y uso en laboratorio  
Requisitos de compatibilidad electromagnética (CEM.)  
Electrical equipment for measurement, control and laboratory use  
EMC requirements.  
Matériel électriques de mesure, de commande et de laboratoire  
Prescriptions relatives à la CEM.

EN 61010-1 | Requisitos de seguridad de equipos eléctricos de medida, control y uso en laboratorio  
Parte 1: Requisitos generales  
Safety requirements for electrical equipments for measurement, control and laboratory use  
Part 1: General requirements.  
Règles de sécurité pour appareils électriques de mesurage, de régulation et de laboratoire  
Partie 1: Prescriptions générales.

EN 61010-2-020 | Requisitos particulares para centrifugadoras de laboratorio.  
Particular requirements for laboratory centrifuges.  
Prescriptions particulières pour centrifugeuses de laboratoire.

Fdo: Alfonso Ainciburu Sanz  
DIRECTOR | GERENTE

BERIAIN a 16 de MARZO de 2009

# CENTRÍFUGA MICROHEMATOCRITO MANUAL MICROHEMATOCRIT CENTRIFUGE CENTRIFUGEUR MICROHÉMATOCRITE

## Mahita

MODELO 2900 / MODEL 2900 / MODÈLE 2900

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- ◆ Balance the centrifuge by using capillary tubes with the same weight, shape and size and put them in geometrically symmetrical positions. If it is necessary, use capillaries filled with water. In a centrifuge it is very important the concept of weight so it is advisable to have a balance to compensate the weight of the different tubes.
- ◆ This is absolutely necessary because an upset of the rotor may produce vibrations that debilitate the centrifuge and rotor increasing the rate of broken tubes and producing a worse sedimentation.
- ◆ Never force the stoppage of the centrifuge, it is a risk of accident and debilitation and it may produce as well a mixing of the separated parts.
- ◆ Change of fuse: the fuse compartment (1.3) is located at the base of the socket. With the help of a sharp object remove the fuse compartment and change the fuse.

### Cleaning

- ◆ Never use scourers or substances that can grate for cleaning metallic parts such as stainless steel, aluminium, coatings, etc. as they damage the centrifuge and produce an early ageing of the equipment.
- ◆ Use a fluff-free cloth dampened with soaped water that does not contain abrasives.
- ◆ The Nahita digital microhematocrit centrifuge 2900 can reach 12000 rpm. If any tube breaks, the interior of the centrifuge must be cleaned and disinfected with a disinfectant agent and the broken parts must be removed; otherwise different components of the centrifuge will suffer important damages and due to the lack of cleaning mistakes in subsequent analysis will take place.

**ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.**

### INSTRUCTIONS ON ENVIRONMENT PROTECTION

At the end of its life cycle, please, do not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

Thank you for choosing this equipment. We sincerely wish that you enjoy your Nahita microhematocrit centrifuge 2900. We highly recommend looking after this equipment according to what is stated in this manual.

Nahita develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

**TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.**

Please bear in mind the following:

- ◆ This manual is inseparable from the Nahita microhematocrit centrifuge 2900, so it should be available for all the users of this equipment.
- ◆ You should carefully handle the centrifuge avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- ◆ Never dismantle the different pieces of the centrifuge to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the equipment.
- ◆ To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting Nahita Technical Assistance Department by email to [asistencia@auxilab.es](mailto:asistencia@auxilab.es) or by telephone: +34 807 117 040 (1.30 Euros/min).
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the equipment warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the centrifuge warranty.
- ◆ Neither fuses nor accessories (including their loss), are covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Nahita Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Please do not forget filling the warranty certificate and send it before 15 days after the date of purchase.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.

**ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.**

**INDEX OF LANGUAGES**

Spanish ..... 2 -10  
 English..... 11 -18  
 French ..... 19-27

**INDEX OF CONTENTS**

1. USES OF THE INSTRUMENT..... 12  
 2. DESCRIPTION ..... 12  
 3. TECHNICAL SPECIFICATIONS ..... 13  
 4. THEORETICAL FOUNDATIONS ..... 13  
 5. INSTALLATION / SETTING UP..... 14  
 6. MAINTENANCE AND CLEANING ..... 17  
 ANNEX 1: CE CERTIFICATE ..... 28

**1. USES OF THE INSTRUMENT**

Nahita microhematocrit centrifuge model 2900 is a silent, vibration-free desktop model, with a robust design, analogic adjustment and capacity for 24 microhematocrit capillaries. It is equipped with electric brake, 10 min timer and speed selector with 2 position thus being ideal for determination of the hematocrit value in research and diagnostic laboratories.

**2. DESCRIPTION**

- |                            |                             |
|----------------------------|-----------------------------|
| 1.1. Lid                   | 1.6. Speed selecting button |
| 1.2. Opening button        | 1.7. Start                  |
| 1.3. Fuse holder           | 1.8. Stop                   |
| 1.4. Socket                | 1.9. Timer                  |
| 1.5. On/Off general switch | 1.10. Brake button          |

- ◆ Do not put the centrifuge near any warm supply (burners, blowlamps...), nor expose it directly to the sun, etc. Avoid vibrations, dust and dry environments.
- ◆ During operation, dangerous material such as flammable liquids or pathological material should be out of the working area.
- ◆ Nahita centrifuges are not suitable for treatment of flammable liquids, explosive samples or samples capable of forming explosive mixtures.
- ◆ He centrifuge is equipped with a mechanical lid locking system and a safety device that stops the rotor in case the lid is open during operation, and that avoids starting operation if the lid is open.
- ◆ When you are not using the centrifuge for a long period of time please make sure it is unplugged in order to avoid possible accidents.
- ◆ It is essential to have the equipment switched off and unplugged from the net before cleaning, checking components or replacing any piece (e.g. replacement of a fuse).
- ◆ Never try to repair the centrifuge by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the equipment (burns, hurts...).
- ◆ Try not to spill any liquid on the control panel, though it is properly insulated. In case you have any doubts do please immediately unplug the centrifuge (see Maintenance).
- ◆ Made under the European regulations for electrical security, electromagnetic compatibility and security on machines.

**6. MAINTENANCE AND CLEANING**

To get the best results and a higher duration of this equipment it is essential to follow the processes of use.

*Note: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.*

- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of this equipment.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ The centrifuge is supplied with a Schuko standard wire. It has to be plugged to an earth connection and the socket should be handy and ready to unplug the equipment in case of emergency.
- ◆ Never try to repair the centrifuge by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the centrifuge (burns, hurts...) or damages in nearby equipments.
- ◆ In the event of breakdown please contact your distributor to overhaul through Nahita Technical Assistance Department.
- ◆ Use the appropriate tubes resistant to the R.C.F. used.

- ◆ Press the general switch (1.5) at the front panel; the stopped centrifuge indicator (red) will light up.
- ◆ Select the centrifugation speed (11000 or 12000 rpm) by the corresponding knob (1.6)
- ◆ Turn the timer (1.9) to select the centrifugation time.
- ◆ Press Start (1.7); the centrifuge begins to operate.

**Very important:** before starting operation make sure that the brake button (1.10) is in off position.

- ◆ Once the selected time has passed, the centrifuge will stop. To stop operation before the selected time has passed, press Stop (1.8).
- ◆ To stop the centrifuge rapidly without affecting the sedimentation result, activate the brake mechanism by pressing brake button (1.10) to the On position.

**Note:** never take the samples out before the rotor is completely stopped. It is very important not to open the lid and let the rotor stop without applying any external force to avoid accidents.

- ◆ Turn the equipment off by pressing the general switch (1.5). Open the lid by pressing the opening button (1.2) and take the capillary tubes out.

#### Hematocrit value determination

For this determination can be used a hematocrit reader that gives a direct result or a ruler:

Hematocrit reader: put the capillary tube so as the beginning of blood (without taking into account the sealing rubber) coincides with the line 0 and move the capillary tube rightwards or leftwards until the line 100 coincides with the end of the plasma. The hematocrit value will be indicated by the line that matches with the separating area between the plasma and the cells.

Ruler: if a ruler is used, the hematocrit value is calculated as follows:

A-----100  
B----- X

Being,


A: total length.

B: length of the corpuscular part.

X: hematocrit value in %.

#### Security

- ◆ The centrifuge must be used by previously qualified staff that knows how the equipment works thanks to the user manual.
- ◆ Put the centrifuge on top of a horizontal, plane and stable table making a free space at least at 30 cm per side.

### 3. TECHNICAL SPECIFICATIONS

<b>Code</b>	52900000
<b>Model</b>	2900
<b>Adjustment</b>	Analogical
<b>Rotor</b>	Microhematocrit with lid
<b>Capacity</b>	24 capillary tubes
<b>Capillary dimensions (Øxlenght)</b>	1.5x75 mm
<b>Speed</b>	11000 or 12000 rpm
<b>Max. RCF</b>	15300 g
<b>Timer</b>	0-10 min
<b>Brake</b>	2 positions (on/off)
<b>Power</b>	220-230 V, 50 Hz
<b>Dimensions</b>	31x30x27 cm
<b>Weight</b>	13 Kg


### 4. TEORICAL FOUNDATIONS

Centrifugation is a separation technique based on the movement of the particles that are displaced to the distal end of the rotation axis depending on their shape and mass.

The Relative Centrifugal Force (R.C.F.) is the required force to separate particles. It is given as some value times that of gravity (g).

The R.C:F. is calculated with the following formula:

$$F.C.R. = 1.118 \cdot 10^{-5} \cdot r \cdot n^2$$

$1,118 \times 10^{-5}$  is a constant.

$r$  = radius, the horizontal distance in centimetres from the rotation axis to the bottom of the sample tube.

$n$  = rotation speed expressed in revolutions per minute (r.p.m.).

### Microhematocrit technique

The microhematocrit technique is very used due to the small quantity of sample required and because it is a very simple technique that allows the analysis of several samples at the same time.

The hematocrit value of a complete blood sample is the relation between the erythrocyte volume and the total volume expressed as a percentage.

If a blood sample treated with an anticoagulant (heparin, EDTA, oxalate) is centrifuged at a high number of r.p.m., the corpuscular component will remain at the bottom of the tube due to its higher density while the plasmatic component will stay above. The hematocrit value is obtained by measuring the quantity and relative proportion of erythrocytes respecting to the total.

**Technique:** The capillary is filled with blood until three quarters of its total capacity. The filling of the tubes is made by capillarity, inclining the capillary downwards while filling up. The end of the capillary is sealed with wax or by flaming with a burner.

## 5. INSTALLATION / SETTING UP

### Preliminary inspection

- ◆ Unwrap the centrifuge, take off the involving plastic and take off the polyspan protection in which it comes fitted.
- ◆ Take off all protections and, without connecting the centrifuge to the net, check manually that the rotor spins freely, without any obstacle. Make sure that the equipment does not present any damage because of the shipment. In case it presents any damage, tell it immediately to your transport agent or dealer so that they can make the claims in the correct time limit.
- ◆ Please keep the original wrapping; you will always need it for returns enclosed with all the accessories supplied.

Please check that all the accessories are enclosed with the equipment:

- Angular rotor for 24 microhematocrit tubes
- Schuko wire
- Instructions manual
- Warranty certificate

**We will only accept any equipment return within 15 days after delivery and provided it comes in its original wrapping.**

### Installation

Before using this instrument, it is convenient for you to familiarize with its components and basic essentials, as well as with its control functions.

**PLEASE READ THOROUGHLY THE INSTRUCTIONS BEFORE CONNECTING AND OPERATING WITH THIS EQUIPMENT WITH THE AIM OF ACHIEVING THE HIGHEST FEATURES AND THE MAXIMUM DURABILITY OF THE INSTRUMENT.**

- ◆ Please put the centrifuge on top of a horizontal, plane and stable table making a free space at least at 30 cm per side. Do not put the equipment near any warm supply (burners, blowlamps...), nor expose it directly to the sun, etc.
  - ◆ During operation, people and dangerous material such as flammable or infectious liquids must be out of the working area.
  - ◆ Make sure that the general switch (1.5) and the brake button (1.10) are in Off position and insert the Schuko wire that feeds the AC electric current in the base of current 220 V, 50 Hz  $\pm$  10% provided with earth wire and to the other end to the centrifuge connector (1.4).
- Neither the manufacturer nor the distributor will assume any responsibility for the damages produced to the equipment during its installation or damages to persons suffered by the improper use of the electric connection. The tension should be 220 V, 50 Hz  $\pm$  10%.*
- ◆ If you are not using the centrifuge for a long period of time please make sure it is disconnected from the net and protected from dust (this way you will avoid accidents and will extend its working-life).

### Setting up (Operation mode)

- ◆ Press the opening button (1.2) and open the lid of the centrifuge (1.1).
- ◆ Unscrew the lid of the rotor by turning the central screw counter-clockwise.
- ◆ Placing the capillary tubes. The centrifuge must be perfectly equilibrated. Capillaries must be placed with the sealed ends outwards, perfectly adjusted to the external rubber rim and fitted into their corresponding grooves. If capillaries are not correctly placed when centrifugation starts they can get broken.
- ◆ Screw again the lid of the rotor by turning the central screw clockwise.

**NOTE: Either the rotor and the cover of the rotor present a black mark; when closing the cover make both marks coincide.**

- ◆ Close the lid of the centrifuge (1.1). As a security measure, the centrifuge will not be operative while the lid stays open.

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- ◆ Press the general switch (1.5) at the front panel; the stopped centrifuge indicator (red) will light up.