DECLARACIÓN DE CONFORMIDAD CE
CE DECLARATION OF CONFORMITY
DÉCLARATION DE CONFORMITÉ CE

El fabricante | The manufacturer | Le fournisseur:

AUXILAB, S.L.

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Note: According to the in force legislation regarding “Non-automatic weighting instruments” in which balances are included, by means of writ dating from 22nd October, 1994 (BOE 1/3rd/95), these balances must not be used for:

- Commercial transactions
- Calculating taxes, tariffs, rates, indemnities and other similar canons
- Judicial surveys
- Pharmaceutical medicine preparations, as well as analysis made in medical or pharmaceutical laboratories
- Determining the price or total amount in sale price or in pre-packaged preparations

10. ACCESSORIES

WEIGHTS (Check the available codes on our web site www.auxilab.es).

- Class E2. Antimagnetic stainless steel.
- Class F1. Antimagnetic stainless steel.
- Class F2. Antimagnetic stainless steel.
- Class M1. Carbon steel.
- Set f weights.

WEIGHING TRAYS.

<i>Code</i>	<i>Description</i>
40350045	Weighing tray 4.5x4.5 cm, B/500
40350080	Weighing tray 8x8 cm, B/500

INSTRUCTIONS ON ENVIRONMENTAL PROTECTION

At the end of its life cycle, please, do not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking. By recycling material or by other forms of reutilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

Thank you for choosing this equipment. We sincerely wish that you enjoy your Nahita precision balances, series 5023. We highly recommend looking after this equipment according to what is stated in this manual.

Nahita develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT

Please bear in mind the following:

- ◆ This manual is inseparable from the Nahita precision balances 5023, so it should be available for all the users of this equipment.
- ◆ You should carefully handle the balance avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- ◆ Never dismantle the different pieces of the balance to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the balance.
- ◆ To prevent fire or electric discharges avoid the contact of the electric circuit with any liquid. In case it may happen unplug the equipment immediately. Avoid dry or dusty environments.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting Nahita Technical Assistance Department by email to asistencia@auxilab.es or by telephone: +34 807 117 040 (0.35 Euros/min).
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the equipment warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the warranty.
- ◆ Neither fuses nor accessories (including their loss), are covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Nahita Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.

9. MAINTENANCE AND CLEANING

Note: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.

To get the best results and a higher duration of the balance it is essential to follow the processes of use:

- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of this equipment.
- ◆ Since it is a precision equipment, prevent it from vibrations, air flows or any hit or sudden movement.
- ◆ Do not put the balance near any warm supply (burners, blowlamps...), nor expose it directly to the sun, etc.
- ◆ Do not put the balance near any warm supply (burners, blowlamps...), nor expose it directly to the sun, etc.
- ◆ Never use the balance in a wedged location as for example a shelf.
- ◆ Never use sharp objects as pens, etc... to press the buttons of the control panel; only use your fingers.
- ◆ Never place on the pan an object heavier than the maximum capacity of the balance, the sensor could be damaged.
- ◆ Neither submerge the balance nor spill liquids on it.
- ◆ If any liquid comes into contact with the electric parts of the balance please do immediately disconnect it from the net and send it to the technical service, as soon as possible, for its checking and adjustment.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ Never try to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty.
- ◆ In the event of a breakdown, contact your wholesaler to repair the equipment through the Nahita Technical Service Department.

Cleaning

- ◆ Never use scourers or substances that can grate for cleaning metallic parts such as stainless steel, aluminium, coatings, etc. as they damage the balance and produce an early ageing of the equipment.
- ◆ Use a fluff-free cloth dampened with soaped water that does not contain abrasives.

ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

INDEX OF LANGUAGES

1. Spanish	2-13
2. English.....	14-25
3. French.....	26-37

INDEX OF CONTENTS

1. USES OF THE INSTRUMENT	15
2. DESCRIPTION	16
3. TECHNICAL SPECIFICATIONS	17
4. INSTALLATION / SETTING UP.....	17
5. CALIBRATION.....	19
6. FUNCTIONS OF THE BALANCE.....	20
7. SOFTWARE	22
8. TROUBLESHOOTING.....	23
9. MAINTENANCE AND CLEANING	24
10. ACCESSORIES	25
ANNEX I: CE CERTIFICATE	38

1. USES OF THE INSTRUMENT

Nahita digital precision balances are basic instruments at any laboratory to perform weighing tasks in which a high precision is required. Developed with the most advanced technology, they are reliable and durable equipments, and its solid, clean and ergonomic design is suitable also for education, research or industrial use.

- "PRT-2" Automatic printing every 1 minute.
- "PRT-3" Automatic printing every 2 minutes.
- "PRT-4" Continuous printing.


8. TROUBLESHOOTING

PROBLEM	CAUSE	SOLUTION
No digits are displayed	Check the balance receives electrical current	Adjust the socket
Value on display changes constantly	1 The balance is in a place with vibrations 2 The door of the balance chamber is not perfectly closed 3 There is an object between the pan and the pan support that makes the balance unstable	1 Put the balance in a place free from vibrations 2 Close all the doors of the balance chamber 3 Check the interior of the balance is clean and there is no object
The value displayed is wrong	1 The balance is not calibrated 2 The balance has not been tared 3 The balance is not well levelled	1 Calibrate the balance 2 Tare the balance before weighing 3 Level the balance until the bubble is well centred
"-----" is shown at the upper part of display	1 The weight exceed the maximum permitted 3 The balance is not calibrated	1 Remove the weight and calibrate the balance 2 Calibrate the balance according to the "Calibration" point of the manual (time of stabilization must not be higher than 8 seconds)
"-----" is shown at the lower part of display	The weighing pan is not well placed	Remove and replace correctly the pan

2. DESCRIPTION


- 1.1 Control panel
- 1.2 Chamber
- 1.3 Pan
- 1.4 RS232 output

- 1.5 Socket
- 1.6 Screwed legs (back legs)
- 1.7 Bubble level

- 2.1 LCD display
- 2.2 ON
- 2.3 OFF
- 2.4 TAR: Tare
- 2.5 CAL: Calibration

- 2.6 COU: Piece counting
- 2.7 UNT: Weighing unit selection
- 2.8 ASD: Weighing mode selection
- 2.9 INT: Weighing mode selection
- 2.10 PRT: Data printing

3. TECHNICAL SPECIFICATIONS

<i>Code</i>	55023210	55023220	55023230	55023250
<i>Capacity</i>	100 g	200 g	300 g	500 g
<i>Readability</i>	0.001 g			
<i>Repeatability</i>		± 0.005 g (0 ≤ m ≤ 50)		
		± 0.010 (50 ≤ m ≤ 200)		
		± 0.015 (m ≥ 200)		
<i>Pan diameter</i>	80 mm	110 mm	110 mm	110 mm
<i>Chamber dimensions</i>		215x195x200 mm		
<i>External dimensions</i>		320x200x290 mm		
<i>Power</i>		220V 50±1HZ		
<i>Weight</i>		7.5 kg		

Among their main features it is worth noting:

- ◆ Quick setting up and stabilization
- ◆ Easy-to-use and easy-to-read LCD display
- ◆ Electromagnetic sensor
- ◆ Panel resistant to liquids splashes, dust, etc.
- ◆ Bubble level and screwed legs for a perfect levelling of the balance
- ◆ Weighing, piece counting and lower weighing functions
- ◆ Weighing unit selection: grams (g), carats (ct), ounces (oz)
- ◆ Substrative tare in the whole weighing range
- ◆ Calibration through external calibration weight (included)
- ◆ RS232 output

4. INSTALLATION/SETTING UP

Preliminary inspection

- ◆ Unwrap the balance, take off the involving plastic and take off the polyspan protection in which it comes fitted.
- ◆ Make sure that it does not present any damage because of the shipment. In case the balance presents any damage tell it immediately to your transport agent or dealer so that they can make the claims in the correct time limit.

We will only accept any equipment return within 15 days after delivery and provided it comes in its original wrapping with all the accessories and documents supplied.

Nahita balances series 5023 are equipped with a system to perform density measurements and weighings of magnetical products that cannot be weight inside the chamber.

Prepare the balance for lower weighing:

- ◆ Put the balance on a table suitable for performing this kind of measurements; the table must be provided with a hole through which a weighing wire or thread will pass.
- ◆ At the bottom of the balance, remove the lid to uncover the weighing hook, and then tie the weighing thread on this hook.
- ◆ An appropriate container (not supplied) must be held to the free end of the thread (Figure 4).
- ◆ Put the sample in this container and proceed as a normal weighing.

6.2 PIECE COUNTING

With this function, the balance can calculate the quantity of pieces placed on the pan.

For this function proceed as follows:

- ◆ Press COU (2.6) successively and select one of the following options:

"-COU-0"	grams
"-COU-1"	5 pieces
"-COU-2"	10 pieces
"-COU-3"	25 pieces
"-COU-4"	50 pieces

- ◆ Once selected, the option blinking on display (e.g.: "COU-10") indicates the number of pieces that must be put on the pan and that will be the reference for the later count.
- ◆ Put the number of pieces on the pan.
- ◆ Press CAL (2.5), "CAL----" will be displayed for a few seconds and finally the reference number of pieces.
- ◆ Remove the pieces from the pan.
- ◆ "0" will be displayed.
- ◆ Put on the pan the sample of pieces to be count and the quantity will be displayed.

7. SOFTWARE

The balance is equipped with universal standard RS232 output and can be connected to printers and computers.

The output data format is:

Bauds 1200.

10 bits (1 start bit (0), 8 data bits (code ASCII) and 1 stop bit (1)).

Data printing.

By pressing PRT (2.10) select the desired data output mode:

- "PRT-0" Immediate printing by pressing PRT (2.10) key
- "PRT-1" Automatic printing every half minute.

◆ If it is necessary to tare the balance because a substance contained in a recipient is needed to be weight, put the empty recipient on the pan and press TAR (2.4). Wait until "0.000 g" is displayed, the balance is ready to weight the substance.

Note:

1. The value of subtractive tare is saved in the balance memory until TAR (2.4) is pressed again.

2. If more than one substance is needed to be weight in the same recipient and their independent weights are also necessary, add the first substance and when the necessary weight is displayed press TAR (2.4). "0.000 g" will be displayed so as the next substance can be weight independently.

◆ If it is necessary to weight a substance that is already contained in a recipient, put the recipient with the substance on the pan. When its weight stabilizes, press TAR (2.4) and "0.000 g" will be displayed. By means of a spatula, transfer the substance to another recipient; the weight of the substance being transferred will be displayed with minus symbol.

◆ If it is necessary to compare the weights of two different substances, put one of them on the pan and press TAR (2.4). When "0.000 g" is displayed, remove the substance and its weight with minus symbol will be displayed. Put the second substance on the pan and the difference of weight between both substances will be automatically displayed.

Lower weighing

Figure 1: Lower part of the balance



Figure 2: Detail of weighing

Please check that all the accessories are enclosed with the equipment:

- Weighing pan
- Standard Schuko wire
- Calibration weight
- User's manual
- Warranty certificate

Installation

Before using the balance, it is convenient for you to familiarize with its components and basic essentials, as well as with its control functions.

PLEASE READ THOROUGHLY THE INSTRUCTIONS BEFORE CONNECTING AND OPERATING WITH THIS EQUIPMENT WITH THE AIM OF ACHIEVING THE HIGHEST FEATURES AND THE MAXIMUM DURABILITY OF THE INSTRUMENT.

- ◆ Put the balance on top of a horizontal, plane and stable table making the enough free space around for a comfortable operation.
- ◆ Since it is a precision equipment, prevent it from vibrations, air flows or any hit or sudden movement.
- ◆ Since it is equipped with an electromagnetic sensor, avoid the proximity of magnetic material that could interfere in its proper operation.
- ◆ Do not put the balance near any warm supply (burners, blowlamps...), nor expose it directly to the sun, etc.
- ◆ Level the balance by turning the screwed legs (1.6) so as the bubble level is perfectly centred (1.7)
- ◆ Insert the Schuko standard wire in the base of current 220 V, 50 Hz $\pm 10\%$ provided with earth wire and the other end to the balance socket (1.5)

Neither the manufacturer nor the distributor will assume any responsibility for the damages produced to the equipment during its installation or damages to people suffered by the improper use of the electric connection. The tension should be 220 V, 50 Hz $\pm 10\%$.

- ◆ Put the stainless steel weighing pan.
- ◆ The balance is ready for its setting up.

Setting up

NOTE: Before using the balance, it is necessary to calibrate it according to the "Calibration" point of this manual. Whenever observing important variations in weighing or not correct weighings, balance should be calibrated.

- ◆ Make sure there is nothing on the pan and press ON to turn the balance on. Display

- (2.1) will turn on and the balance makes an autotest, finally “0.000 g” is displayed.
- ◆ By pressing ASD (2.8) and INT (2.9) select the desired weighing mode among the following possible combinations:
 - INT-1 / ASD-2 Fast weighing (when an absolute precision is not required)
 - INT-3 / ASD-2 Normal weighing (in normal conditions)
 - INT-3 / ASD-3 Slow weighing (when weight is not stable due to unfavourable environmental conditions)
 - INT-0 / ASD-0 For factory adjustments only. User must not use this mode.
 - ◆ By pressing UNT (2.7) select the desired weighing unit: grams (g), carats (ct) (upper part of “g” symbol), ounces (oz) (lower part of “g” symbol). Se podrían dibujar estos dos símbolos en rojo?
 - ◆ The balance is now ready to operate.
 - ◆ Turn the balance off by pressing OFF (2.3) once operation is finished.

Note: By default the balance always turns on with its standard configuration: weighing in grams and mode INT-3/ADS-2.

The following chart shows the unit conversion

ABBREVIATURE	NAME	CONVERSION
g	Gram	1 g
ct	Carat	0.1999694 g
oz	Ounce	28.3495231 g

NOTE: To avoid possible fluctuations in weighing, it is necessary to keep the balance chamber closed as well as avoiding vibrations on the working surface during measurements.

5. CALIBRATION

Before starting using the balance and whenever observing important variations in weighing or not correct weighings, balance should be calibrated. Some factors such as variation of gravity depending on latitude, change of working place as well as an inappropriate use of the equipment lead to the necessity of recalibrate the balance.

¡Important!: Do not move the balance during calibration

Check the correct calibration weight according to the model of balance

Calibration	Weight
55023210	100 g
55023220	200 g
55023230	200 g
55023250	200 g

For calibration, proceed as follows:

- ◆ Make sure that there is nothing on the pan and turn the balance on by pressing ON (2.2).
- ◆ The balance makes an auto-test and finally stabilizes in “0.000 g”.

Note: If “0.000 g” is not displayed, press TAR (2.4) until this value is displayed.

- ◆ Press CAL (2.5).
- ◆ “CAL-XXX” will be displayed blinking; XXX value will show the weight of the calibration weight required.
- ◆ Put the calibration weight on the pan.
- ◆ At the beginning “CAL----” is displayed for a few seconds and finally the weight of the corresponding calibration weight (e.g.: “200.000”) whose reference is automatically registered.
- ◆ Remove the calibration weight from the pan.
- ◆ When “0.000 g” is displayed the calibration has finished correctly. If this not happens, repeat the calibration process.

It is necessary to make a verification of balance calibration:

- ◆ When “0.000 g” is displayed, put again on the pan the calibration weight and check if the value displayed coincides with its weight.

6. BALANCE FUNCTIONS

6.1 WEIGHING

Normal weighing

- ◆ Take into account the room temperature; if the balance has been moved wait at least for an hour to reach the total stabilization of the temperature of the balance.
- ◆ Calibrate the balance if necessary according to “Calibration” point.
- ◆ Make sure there is nothing on the pan and turn the balance on by pressing ON .
- ◆ When “0.000 g” is displayed and the stabilization symbol (o) disappears, the balance is ready to perform measurements (if necessary, select the desired weighing unit by pressing UNT (2.7)). Weighing measurements can be taken into account once the stabilization symbol (o) disappears.