AUXILAB S.L.

## CE DECLARATION OF CONFORMITY

PRECISION MULTIFUNCTION BALANCES of Auxilab,S.L for the Directive of machines (89/392/CEE modified) and the regulations adopted for its transposition

NAME OF THE MANUFACTURER / IMPORTER:

AUXILAB, S.L.

ADDRESS:

Polígono Morea Norte, 8  
31191 Beriáin (Navarra)

WE STATE THAT:

**PRECISION MULTIFUNCTION BALANCES SERIES 5033**  
**Models 55033210, 55033220, 55033230, 55033250**

Are designed and manufactured according to:

- ♦ Directive 89/392/CEE, including the modifications and the national regulations that transpose them.
- ♦ Directive 73/23/CEE, modified over the electric security.
- ♦ Directive 89/336/CEE, modified over the electromagnetic compatibility.
- ♦ And that the following harmonized rules have applied (or part of them):  
UNE 292-1, UNE 292-2, UNE 292-2/A1, UNE 614-1, UNE 1050, UNE 294, UNE 894-1, UNE 894-2, UNE 60204, UNE 61010-1.

BERIAIN 05th November 2006

Signed by: ALFONSO AINCIBURU SANZ  
DIRECTOR/MANAGER

Polígono Morea Norte, 8 31191 Beriain (Navarra) - Spain. Tel. 948 310 513 Fax 948 312 071  
Internet: [www.auxilab.es](http://www.auxilab.es) · Email: [correo@auxilab.es](mailto:correo@auxilab.es)

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*This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.*

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your fingers.

- ◆ Never place on the pan an object heavier than the maximum capacity of the balance, the sensor could be damaged.
- ◆ Neither submerge the balance nor spill liquids on it.
- ◆ If you use batteries, remove them when you are not using the balance for a long period of time.
- ◆ If any liquid comes into contact with the electric parts of the balance please do immediately disconnect it from the net and send it to the technical service, as soon as possible, for its checking and adjustment.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.

#### CLEANING

- ◆ Never use scourers or substances that can grate for cleaning metallic parts such as stainless steel, aluminium, coatings, etc. as they damage the balance and produce an early ageing of the equipment.
- ◆ Use a fluff-free cloth dampened with soaped water that does not contain abrasives.

**ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.**

*NOTE: According to legislation currently in force, regarding "Non-automatic weighing instruments" in which balances are included, by means of writ dating from 22nd December 1994 (BOE 3/1/95), Nahita balances series 5033 must not be used for:*

- Commercial transactions.
- Calculation of taxes, tariffs, duties, remunerations, indemnities and other similar canons.
- Judicial surveys.
- Pharmaceutical medicine preparations, as well as analysis made in medical or pharmaceutical laboratories.
- Determination of the price or total amount in retail sales and in the preparation of pre-packaged products.

#### INSTRUCTIONS ON ENVIRONMENT PROTECTION

At the end of its life cycle, please, do not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

**7. TROUBLESHOOTING**

TROUBLEA	CAUSA	SOLUCIÓN
Digits are not displayed	Check the balance receives electrical current. The adapter can be damaged. The fuse is blown.	Tighten the power plugs. Change the current adapter. Replace the fuse.
The value on display does not stabilizes	The balance is in a place with vibrations. The door of the balance is not completely closed. There is an odd object under the pan causing destabilization.	Move the balance to a place without vibrations. Close all the doors of the balance.. Check that the balance is clean inside and that there are not strange objects.
The displayed value is wrong	Balance is not calibrated. Balance has not been tared before weighing. Balance is not well leveled.	Calibrate the balance. Tare the balance before each weighing. Level the balance until the bubble situated at the back side is centred.

CASTELLANO

**8. MAINTENANCE**

To get the best results and a higher duration of the pH-meter it is essential to follow the processes of use.

*Note: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.*

- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of this equipment.
- ◆ Prevent the balance from sudden movements and knocks, as well as from direct sunlight or air flows. The balance is a precision instrument, you must handle it carefully.
- ◆ The balance is supplied with an adapter. It has to be plugged to an earth connection and the socket should be handy and ready to unplug the equipment in case of emergency.
- ◆ Never unplug the adapter by pulling the wire, do it from the base.
- ◆ Place the balance in a location with adequate ventilation to prevent heat accumulation. Never use the balance in a wedged location as for example a shelf.
- ◆ Never use sharp objects as pens, etc... to press the buttons of the control panel; only use

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

#### COMMAND MODE

Data will be transmitted from the balance to an external device when they receive an order. Balance can be controlled from the external device by the following orders in the same way as by the frontal panel of the balance.

#### ORDERS (Commands)

<P><CR><LF> ON/OFF order (50 0D 0A): the same as ON/OFF function from the frontal panel.


<Q><CR><LF> Printing order (51 0D 0A): the same as PRINT function from the frontal panel.

<R><CR><LF> Tare order (52 0D 0A): the same as TARE function from the frontal panel.


<U><CR><LF> Mode order (55 0D 0A): the same as MODE function from the frontal panel.

#### CONNECTION WITH EXTERNAL SERVICES:

##### 1. With computer

##### 2. With print

**RS-232C OUTPUT. SPECIFICATIONS**

**Transmission way:** double, asynchronous and bi-directional

**Data format:** digital transmission speed: 300, 600, 1200 and 2400.

**Informative bit:** 7

**Bit of parity:** 1

**Bit of stop:** 1

**Code:** ASCII

Data output: Auto-print, "PRINT" key (1.6), data emission

Data output format: ST, +010.000xxg <CR><LF><LF>  
US, +010.000xxg <CR><LF><LF>

Note:

ST = Stable (53H, 54H)

US = Unstable (55H, 53H)

X = Space (20H)

<CR> = Carriage return (ODH)

<LF> = Line Feed <OAH>

**DATA OUTPUT****AUTO-PRINT MODE**

There are three different ways in Auto Print mode:

## 1. AUTO PRINT A

**Note: the pan must be clean and screen must show Zero after each weighing in AUTO PRINT A mode.**

Data transmission starts when balance begins to stabilize (stability indicator appears (1.7)).

## 2. AUTO PRINT B

Data transmission happens whenever balance stabilizes.

## 3. AUTO PRINT C

Data are transmitted continuously when plugging the balance without giving any additional order to neither the balance nor the computer or printer.

**PRINT MODE**

Data will only be transmitted if PRINT key (1.6) is pressed.

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FUNCTION	PARAMETERS (CX-Y)	NAME	MEANING
C4	C4-0*	Digital	2400
	C4-1	transmission	1200
	C4-2	speed	
	C4-3		300
C5	C5-0*	Data exit	A printing
	C5-1		B printing
	C5-2		Print mode
	C5-3		C printing
	C5-4		Command mode
C6	C6-0*	Back illumination	OFF
	C6-1		ON
C7	C7-0*	Stability band's	1
	C7-1	amplitude allowed	2
	C7-2		4
	C7-3		8
C8	C8-0	Average time	1 second
	C8-1		2 second
	C8-2		4 second
	C8-3*		6 second

\*Factory setting

### DATA RECORDING

- ◆ All the above features (table 3) can be modify directly from the screen.
- ◆ Press TARE and ON/OFF keys simultaneously with the balance turned off.
- ◆ Release the key and 888888888 will be displayed.
- ◆ Press "PRINT" key (1.6).
- ◆ "C1-0" (CX-Y) will be displayed.
- ◆ Press "MODE" key (1.4) to change the X digit between values from 1 to 7 (see table 3).
- ◆ By pressing "PRINT" key (1.6) change the Y value. This digit can take different values depending on X digit value (see table 3).
- ◆ Press ON/OFF key to save the new parameters in the memory.
- ◆ The balance is ready to work with the new parameters, that will be recorded when you turn on/off the balance.

**PERCENTAGE MODE**

In percentage mode the balance perform measurements based on a reference value.

- ◆ By pressing "MODE" key (1.4), select percentage mode (%).
- ◆ "0.00" will be displayed.
- ◆ Press "CAL" key (1.3).
- ◆ By default, "--100.00--" will be displayed; this value indicates the weight percentage corresponding to the reference sample.

**NOTE:** To modify this value, press MODE key (to increase) or PRINT key (to decrease). You can choose between 50, 100, 150, 200, 250, 300, 350 and 400 %.

- ◆ Place on the pan the sample whose weight is going to be considered as the 100 % (or the selected value).
- ◆ Press "TARE" key, and after a few seconds, the corresponding percentage will be displayed.
- ◆ Removed the sample from the pan.
- ◆ "0.00 %" will be displayed.
- ◆ Place the sample to be compared on the pan. As a result of the comparison between the reference sample and the problem sample, the deviation in percentage will be displayed.

**Note:** To set a new reference, press CAL and proceed as previously described (from point 4).

**PARAMETERS**

Nahita balances series 5033 are supplied with a Software to adapt their features, as best as possible, to users' necessities:

TABLE 3

FUNCTION	PARAMETERS (CX-Y)	NAME	MEANING
C1	C1-0		8
	C1-1		16
	C1-2		32
C2	C2-0		1000 digits
	C2-1		2000 digits
	C2-3		5000 digits
	C2-3		10000 digits
C3	C3-0*		ON
	C3-1	Zero and tare	OFF

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- ◆ Turn the balance on by pressing ON/OFF key (1.2) with the pan empty.
- ◆ When 0.000 g and the stabilization indicator (1.7) are displayed, the balance is ready to weight.
- ◆ When tare is necessary because a recipient is going to be used in weighing, place the recipient empty on the pan and press "TARE" key (1.5). Wait until 0.000 g is stabilized; the balance is ready to weight.

### *Note:*

*1. The value of subtractive tare is kept in balance's memory until TARE key (1.4) is pressed once again.*

*2. If you need to weight more than one substance in the same recipient and you require having their separated weights, add the first substance and, after its weight has been displayed, press "TARE" key (1.5). 0.000 g will be displayed, so the next substance can be weight independently.*

- ◆ If you need to weight a substance contained in a recipient, place the full recipient on the pan. When its weight stabilizes, press "TARE" key (1.5) and 0.000 g will be displayed. With the aid of a spatula or similar, transfer the substance to be weight to another recipient. When transferring finishes the substance's weight (with negative sign) will be displayed.

### PIECE COUNTING MODE

In counting mode the balance calculates the number of pieces placed on the pan. Proceed as follows:

1. By pressing "MODE" key, select Counting mode (1.4).
2. "0" will be displayed.
3. Press "CAL" key (1.3).
4. By default, "10" will be displayed; this number indicates the number of pieces that should be placed on the pan and that will be taken as a reference for a later count.

*NOTE: To modify this value, press MODE key (to increase) or PRINT key (to decrease). You can choose between 5, 10, 15, 20, 25, 30, 35 and 40 pieces.*

5. Place 10 pieces of those to be counted, or the selected number of pieces, on the pan.
6. Press TARE key and after a few seconds, the corresponding number of pieces will be displayed.
7. Remove the pieces from the pan.
8. "0" will be displayed.
9. Place on the pan the amount of pieces to be counted and the number of pieces will be displayed.

*Note: To set a new reference, press CAL key and proceed as previously described (from point 4).*

Place the calibration weight on the pan when 0.000 g is displayed, and check that the displayed value coincides with the weight value with a maximum deviation of  $\pm 0.001$  g; otherwise, balance is not correctly calibrated, and calibration should be performed again.

*Note: If "CAL no" is displayed during calibration, it means that calibration weights are not the correct ones because their weight is up  $\pm 10\%$  on the maximum capacity of the balance.*

Change the weights for ones below these limits and do the calibration again. If "CAL no" message remains in spite of using the appropriate weights, please ask Nahita authorised technical service.

### **BALANCE ADJUSTMENT**

◆ If you need more accuracy and precision in measures and you also know the uncertainty of calibration weights, you can adjust manually the weighing values of the balance. Taking as an example that the calibration weight value is 200.5 g, calibration should be performed as follows:

- ◆ Turn the balance on by pressing ON/OFF key (1.2) with the pan empty.
- ◆ Press "CAL" key (1.3) when 0.000 g is displayed.
- ◆ "CAL 0" will be displayed.
- ◆ Press "CAL" key (1.3) again.
- ◆ The maximum capacity of the balance will be displayed (200.0 g in the example).
- ◆ Adjust this value until reaching the exact value of the calibration weight to be used, 200.5 g in the example. For this purpose, press "MODE" or "PRINT" key either to increase or decrease the value initially displayed.
- ◆ Once the calibration weight value has been adjusted, press "TARE" (1.5) key.
- ◆ "CAL 0" will be displayed.
- ◆ Press "TARE" key (1.5) again.
- ◆ "CAL F" will be displayed.
- ◆ Place the calibration weight on the pan.
- ◆ Press "TARE" key (1.5).
- ◆ "CAL END" will be displayed.
- ◆ The adjusted value (200.5 g) will be displayed.
- ◆ Remove the calibration weight from the pan and 0.000 g will be displayed; the balance is correctly calibrated and ready to use.

## **6. FUNCTIONS OF THE BALANCE**

Take into account the temperature of the balance; wait at least one hour for a total stabilization of the temperature if balance has been moved.

Calibrate the balance if necessary, according to "Calibration" point.

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**vibrations on the table in which it is situated the balance while weighing.**

5. Switch the balance off by pressing On/Off key once operation has finished.

**5. CALIBRATION**

Whenever observing important variations in weighing or not correct weighings, balance should be calibrated. Some factors such as variation of gravity according to latitude, change of working place as well as an inappropriate use of the equipment lead to the necessity of recalibrate the balance.

**Important!! Do not move the balance during calibration**

Consult the weights you need to calibrate your balance based on its model.

TABLE 2	Code Model	Max. capacity	Calibration weight
	55033210	100 g	100 g
	55033220	200 g	200 g
	55033230	300 g	300 g

Calibration weights should be previously calibrated.

Choose the calibration weight you need according to the model of the balance to be calibrated; you need a weight with the same value as maximum capacity of the balance (see table 2). Proceed as following:

- Turn the balance on by pressing On/Off key (1.2); the pan should be empty.
- The balance will perform an autotest and finally stabilizes on 0.000 g.

**Note: If 0.000 g does not appear on display, press TARE key (-) until it appears and the stabilization indicator (-) lights up.**


- ◆ Press "CAL" key (1.3) to enter calibration mode.
- ◆ CAL 0" will be displayed.
- ◆ Press "TARE" key (1.5).
- ◆ CAL F" will be displayed.
- ◆ Place the calibration weight on the pan.
- ◆ Press "TARE" key (1.5) again.
- ◆ CAL END" will be displayed.
- ◆ The maximum capacity of the balance will appear on display, it is also the weight of the standard calibration weight used.
- ◆ Removed the calibration weight from the pan.
- ◆ 0.000 g will be displayed, indicating that balance calibration has been well performed.

It is necessary to check the results of calibration:

by the improper use of the electric connection. The tension should be 220 V, 50 Hz  $\pm$  10%.

Place the stainless steel weighing pan. The pan is composed of four pieces that must be inserted according to the order shown in figure:

Check that the level is in the right position (centred), if not, centre the bubble with the knobs situated at the frontal legs of the balance.

The balance is ready to use.

#### SETTING UP

**NOTE:** Before using the balance, it must be calibrated according to "Calibration" point of this manual. Recalibration must be performed whenever observing important variations or mistakes in weighings.

1. Make sure that there is nothing on the pan and switch the balance on (-). The display will light up (-) and the balance will make and autotest that finishes when 0.000 g is visualized on display.

2. Select the desired function: weighing, percentage or pieces count.

3. By pressing MODE key (1.4) select the measuring unit of weight: grams (g), milligrams (mg), ounces (oz) or carats (ct). The following table shows the conversion of these units into grams.

#### UNITS CONVERSION

ABBREVIATION	NAME	CONVERSION
kg	Kilogram	1000 g
g	Gram	1 g
mg	Milligram	1000 g
oz	Ounce	28.3495231 g
ct	Carat	0.1999694 g

4. The balance is ready to weigh.

**NOTE:** To avoid possible oscillations in weighing, it is important to avoid air flow and

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## 4. INSTALLATION

### Preliminary inspection

Unwrap the balance, take off the involving plastic and take off the polispán protection in which it comes fitted. Take off all the protective items and, without connecting the balance to the net, make sure that it does not present any damage because of the shipment. In case the balance presents any damage tell it immediately to your transport agent or dealer so that they can make the claims in the correct time limit.

Please keep the original wrapping; you will always need it for returns enclosed with all the accessories supplied.

Please check that all the accessories are enclosed with the equipment:

- ◆ Pan
- ◆ Calibration weight (according to balance capacity)
- ◆ Cleaning cloth
- ◆ Electric current adapter
- ◆ Instructions manual
- ◆ Warranty

*We will only accept any equipment return within 15 days after delivery and provided it comes in its original wrapping.*

### Installation

Before using the balance, it is convenient for you to familiarize with its components and basic essentials, as well as with its control functions.

**PLEASE READ THOROUGHLY THE INSTRUCTIONS BEFORE CONNECTING AND OPERATING WITH THIS EQUIPMENT WITH THE AIM OF ACHIEVING THE HIGHEST FEATURES AND THE MAXIMUM DURABILITY OF THE INSTRUMENT.**

Please put the balance on top of a horizontal, plane and stable table making a free space at least at 30 cm per side. Do not put the balance near any warm supply (burners, blowlamps...), nor expose it directly to the sun, etc

### Placing the AC adapter


Please insert the wire that feeds the AC electric current in the base of current 220 V, 50 Hz  $\pm$  10% provided with earth wire and to the other end to the balance connector.

*Neither the manufacturer nor the distributor will assume any responsibility for the damages produced to the equipment during its installation or damages to people suffered*

- 1.1 Digital display
- 1.2 ON/OFF
- 1.3 CAL key: calibration
- 1.4 MODE key: select the weighing mode
- 1.5 TARE key
- 1.6 PRINT key
- 1.7 Stability indicator
- 1.8 COUNT: piece count function indicator
- 1.9 %: percentage function indicator
- 1.10 ct (carats) unit indicator
- 1.11 oz (ounces) unit indicator
- 1.12 mg (milligrams) unit indicator
- 1.13 g (grams) unit indicator

**BACK PANEL**

- 2.1 Nivel
- 2.2 Toma corriente
- 2.3 Salida RS 232

**3. TECHNICAL SPECIFICATIONS**

CODE	55033210	55033220	55033230
Capacity	100 g	200 g	300 g
Readability		0.001 g	
Reproducibility		± 0.001 g	
Linearity		± 0.002 g	
Pan diameter		90 mm	
Chamber dimensions		210x175x170 mm	
External dimensions		330x190x310 mm	
Weight		5.25 Kg	

Among its most outstanding characteristics are the following ones:

- ◆ Stable and rapid weighing
- ◆ Easy working, simple and light line
- ◆ Subtractive tare
- ◆ Network connection by means of adapter (included)
- ◆ Automatic calibration
- ◆ Weighing in grams and pieces count and percentage functions
- ◆ RS 232 interface for connecting a computer or a printer

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


**1. USES OF THE INSTRUMENT**

Nahita precision electronic balances are the result of many years of research, design and development in the field of analytical balances. They are basic instruments in every laboratory in weighing processes requiring the highest accuracy. Nahita balances, developed with the most advanced technology, are reliable and durable instruments with a solid, clean and ergonomic design suitable for educational, industrial or researching usages.

**2. DESCRIPTION**

Thank you for choosing this equipment. We sincerely wish that you enjoy your Nahita digital balance, series 5033. We highly recommend looking after this equipment according to what is stated in this manual.

Nahita develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

**TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.**

Please bear in mind the following:

- ◆ This manual is inseparable from the Nahita digital balance series 5033, so it should be available for all the users of this equipment.
- ◆ You should carefully handle the balance avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- ◆ Never dismantle the different pieces of the balance to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the balance.
- ◆ To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting Nahita Technical Assistance Department by email to asistencia@auxilab.es or by telephone: +34 807 117 040 (0.30 Euros/min).
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the equipment warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the warranty.
- ◆ Neither fuses (1A) nor accessories (including their loss), are covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Nahita Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Please do not forget filling the warranty certificate and send it before 15 days after the date of purchase.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.

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