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## 7. TROUBLESHOOTING

PROBLEM	CAUSE	SOLUTION
Digits do not appear on the screen.	- Check that the electrical power gets into the balance. The adapter could be damaged. - The fuse does not work.	- Adjust the power plugs. Change the current adapter. - Replace the fuse.
The value on the screen is changing all the time.	- Your balance is in a place with vibrations. - The balance's door is not completely closed. - An odd object can be between the pan and its placing, causing the destabilization. - The underhook where the thread is hitched in hydrostatic measuring, is open.	- Move the balance to a place far from vibrations. - Close every balance's door. - Check the inner part of the balance is clean and there is no odd object in it.
The value on the screen is wrong.	- The balance is not calibrated. - The balance has not been tared before weighing. - The balance is not leveled.	- Calibrate - Tare the balance before each weighing. - Level the balance until the bubble in the back is centered.

*Note:* According to the in force legislation regarding *Non-automatic weighing instruments* in which balances are included, by means of writ dating from 22nd October, 1994 (BOE 3/1/95), **these balances must not be used for:**

- Commercial Transactions.
- Calculating taxes, tariffs, rates, indemnities and other similar canons.
- Judicial surveys.
- Pharmaceutical medicine preparations, as well as analysis made in medical or pharmaceutical laboratories.
- Determining the price or total amount in sale price or in prepackaged preparations.

- ◆ Do not press the display or any key with a pointed object such as a ballpoint pen or a screwdriver, as this may scratch or damage the key. Use your only your fingers.
- ◆ Do not let objects fall on the pan in order to avoid damages to the weighing system.
- ◆ Never place into the balance an object heavier than the maximum weight of the balance's range, the sensor could be damaged.
- ◆ Never disconnect it by pulling the cable for security reasons.
- ◆ Place the unit in a location with adequate ventilation to prevent heat build-up the unit.
- ◆ Do not submerge the balance in liquids.
- ◆ The balance has a schuko plug. This plug must be connected to a switch connected to ground. This switch must be close in order to be able to disconnect it in case of emergency.
- ◆ In case of damage or doubt it can be solved getting in touch with the NAHITA Technical Service.

**Cleaning**

- ◆ For the cleaning of the inox steam parts use always a special cloth that does not give off fluff. Do not use any cloth that can produce damages to the equipment.
- ◆ For the cleaning of the different pieces of the equipment, we highly recommend to use neutral products that do not have a pH superior than 6 (soaped water without abrasives).

**ATTENTION!!**  
**IF AN EQUIPMENT IS NOT PROPERLY CLEAN AND DISINFECTED IT WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.**

Thank you for choosing this equipment. We sincerely wish that you enjoy your electronic balance NAHITA. We highly recommend looking after this equipment according to what is stated in this manual.

NAHITA develops its products being remarkable in the ergonomics and security of its user. The incorrect use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages. Please read the point of Maintenance, where we expose the security notes.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

**TO GET THE BEST RESULTS AND HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.**

Please bear in mind the following:

- ◆ This manual should be available for all the users of these equipments
- ◆ As an accuracy instrument, we should manipulate carefully avoiding the sudden movements, knocks, free fall of heavy objects, manipulation of sharp objects, and the pouring of liquids in the interior of the balance.
- ◆ Never dismantle the different pieces of the balance to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty.
- ◆ To prevent fire or shock hazard, do not expose the balance to rain or humid environments.
- ◆ If you have any doubt do not hesitate in contacting your distributor or our NAHITA Technical Assistant Department by email [asistencia@auxilab.es](mailto:asistencia@auxilab.es), by telephone 34 948310513 or by fax 34 948310500 or 948312071.
- ◆ Please do not forget filling the Warranty Certificate and send it before 15 days after the date of purchase.

**ATTENTION!!**  
**IF AN EQUIPMENT IS NOT PROPERLY CLEAN AND DISINFECTED IT WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.**

**LANGUAGE INDEX**

Spanish .....	1-10
English.....	11-19
French .....	20-28

**INDEX OF CONTENTS**

1. APPLICATIONS .....	12
2. DESCRIPTION .....	12
3. TECHNICAL SPECIFICATIONS .....	13
4. SETTING UP .....	13
5. FUNCTIONS OF THE BALANCE .....	15
6. MAINTENANCE AND CLEANING .....	17
7. TROUBLESHOOTING .....	19


## 1. APPLICATIONS

The Nahita analytical balance is an essential instrument in laboratory work when weighing tasks need to be carried out. It is a reliable and long-lasting instrument, which has been developed with the latest technology. Its simple, sturdy, ergonomic design makes this balance suitable for any industrial use, and its smooth surface makes the cleaning easier.


## 2. DESCRIPTION

Picture 1


- 1.1. Front panel
- 1.2. Level
- 1.3. Plate
- 1.4. Glass windbreak protection

Picture 2 - Front panel

- 2.1. ON/OFF
- 2.2. Print (out of use)
- 2.3. Tare
- 2.4. Calibration
- 2.5. Digital screen

Picture 3 - Rear panel

- 3.6. Plug socket
- 3.7. Fuse

- ◆ Open the port by removing the plug and you will find the underhook.
- ◆ Hang a fine thread through the hole.
- ◆ Place the balance on a table with a hole in it or a stand specially designed for below weighing, and let the thread through the hole
- ◆ Suspend a plate. For best results it is necessary to recalibrate the balance, erasing the information on the screen and you will see 0, 00000 g.

### Density determination

- ◆ Place the material with an unknown density in the plate when value 0, 00000 g is shown on the screen.
- ◆ Read the value of the material's weight suspended in the air.
- ◆ Press *Tare* key in order to put zero on the display.
- ◆ Immerse the container (with a known value now) in the water at 10°C.
- ◆ Read out the difference of weights between the weight of the material in air and the weight of the material in water.
- ◆ Calculate the density of the water density from the following board:

0°C	0,99984 g/cm <sup>3</sup>
4°C	0,99997 g/cm <sup>3</sup>
10°C	0,99970 g/cm <sup>3</sup>
15°C	0,99910 g/cm <sup>3</sup>
20°C	0,99821 g/cm <sup>3</sup>
25°C	0,99705 g/cm <sup>3</sup>
30°C	0,99565 g/cm <sup>3</sup>

- ◆ Calculate the density of the material:

$$\text{Weight in air (g)} \div \text{water density (cm}^3\text{)} = \text{g/cm}^3$$

## 6. MAINTENANCE AND CLEANING

*To get the best results and a higher duration of this equipment it is necessary to follow the processes of use.*

*Note:* All the processes of use mentioned before will not have any value unless you keep a continued and careful maintenance.

### General recommendations

- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of these equipments.
- ◆ Always use original components and supplies. Other dispositives can be similar but they can damage the equipment.
- ◆ Preserve the balance from heavy shocks and look for a dry place, away from dust, direct sun beams or air flows. Handle the balance gently, since it is an accurate instrument.

Check using an external calibration mass

- ◆ A 50 g standard mass should be used for checking. The mass must be accurate and the exact value should be known.
- ◆ Zero must appear on the screen; if not, then press *Tare* key (2.3) making sure that the pan is empty.
- ◆ Place the mass on the pan. If the difference between the value displayed and exact value of the mass is less than  $\pm 0,0001\text{g}$  your balance is calibrated. If the difference is more than  $\pm 0,0001\text{g}$ , please repeat the calibration procedure with the internal mass.

Note: If "CAL no" is shown in the screen during calibration, this means the balance can not be calibrated because the calibration weight used exceeded the limit of the calibration weight value (not more than 10%).

- ◆ If "CAL no" message remains in spite of using the appropriate weights, please ask NAHITA authorized Technical Service.


Simple weighing

- ◆ Take in mind the balance's temperature if you have moved it. In this case, wait for at least an hour to obtain the plain temperature's stabilization. Calibrate the balance if you consider it is necessary based on the *Calibration* point exposed before.
- ◆ Turn on the display by pressing *On/Off* key (1.7) with the pan empty
- ◆ When the value *0.0000 g* is shown in the screen with stabilized indicator (1.5) the balance is ready to take readings.
- ◆ You can tare the balance if you are going to use a recipient by putting it on the balance, then pressing *Tare* key (2.3) and waiting until the screen shows *0.0000 g*. After the value is stabilized, you can weight.

Note:

1. The value of subtractive tare is kept in the balance's memory until *Tare* key (2.3) is pressed once again.
  2. If you want to weight more than one substance in the same container and you require having their separated weights, press *Tare* key (2.3) again each time you reach your wanted weight and *0.0000g* will be shown on the screen, so you will obtain the independent/separate weights.
- ◆ If you need to weight a substance that is in the interior of a recipient, put the full recipient on the balance, then when the weight gets stabilized, press *Tare* key (2.3) and you will read value *0.0000 g* on the screen. With the help of a spatula or a similar tool, move the substance you want to weight to another recipient. When you finish moving the substance to the other recipient you will obtain its weight with negative sign.

Picture 6

Below weighing

5034 NAHITA balances are provided with a system to measure density and magnetic material's weighing.

Preparation

- ◆ Place the balance on a table with the front portion of the balance out of the edge of the table.
- ◆ A round plastic plug will be found under the balance (6.10).
- ◆ Take off the plug.

**3. TECHNICAL SPECIFICATIONS**

Among its various characteristics we would like to emphasize the following:

- ◆ Stable and rapid weighing.
- ◆ Easiness when using, with a simple and light line.
- ◆ Subtractive tare.
- ◆ Weighing mode: grams.
- ◆ In order to level off the scale it is provided with adjustable feet and centering level ( a spirit level).
- ◆ Automatic calibration via a built in calibration mass.
- ◆ Display of 8 digits.
- ◆ Adapter 220/240 V 50-60 Hz 9VDC.

<i>Code</i>	55034060	55034120	55034200
<i>Capacity</i>	200 g	120 g	60 g
<i>Calibration mass</i>	200 g	100 g	60 g
<i>Readability</i>	0.0001 g		
<i>Repeatability</i>	$\pm 0.0001$ g		
<i>Linearity</i>	$\pm 0.0002$ g		
<i>Operating temperature</i>	5°C ~ 40°C		
<i>Pan's diameter</i>	90 mm		
<i>Chamber size</i>	185 x 225 x 155 mm		
<i>Dimensions</i>	195 x 275 x 406 mm		
<i>Weight</i>	11 Kg		
<i>Power</i>	Adapter AC 220V $\pm 10\%$ AC 50-60Hz		

**4. SETTING UP**Preliminary inspection

- ◆ Unwrap the balance, take off the involving plastic and take off the poliespan protection in which it comes fitted. Please make sure that it does not present any damage because of the shipment. If it does please communicate to your transport agent or suministrator immediately, so that they can make the claims in the correct time limit.
- ◆ Please keep the original wrapping; you will always need it for returns enclosed with all the accessories suministrated.
- ◆ Please check that all the accessories are enclosed with the equipment:
  - Pan.
  - AC power adapter.
  - Wind breaker
  - 2 fuse 0.5 A.
  - Instruction manual.
  - Warranty.

**We will not accept any equipment for return unless it comes in its original wrapping.**

### Installation

- ◆ Before using this instrument, it is convenient for you to familiarize with its components and basic essentials.


**PLEASE READ THOROUGHLY THE INSTRUCTIONS BEFORE CONNECTING AND OPERATING WITH THIS EQUIPMENT.**

- ◆ Please put the balance on top of a horizontal, plane, stable table making a free space at least at 30 cm per side of the balance. Do not put the balance near any warm supply (burners, blowlamp...), do not expose directly to the sun, etc. During its functioning, do not use the equipment near fluids that can give off steams or produce explosive or inflammable mixtures.
- ◆ The balance is suminstrated with a Schuko standard wire.
- ◆ Please insert the wire that feeds the A/C in the base of current 220V  $\pm$ 10% 50-60Hz provided with earth wire.

*The manufacturer and the distributor will not assume any responsibility for the damages produced to the equipment in its installing or damages to persons suffered by the wrong use of the electric connection. The tension should be 220V  $\pm$ 10% 50-60Hz.*

### Setting up

- ◆ Please position the inox pan in the sensor of the balance.
- ◆ Once the balance is installed press the *On/Off* key. It will go on and the screen will show "-.-", once the weight is stabilized you will see on the screen "0.0000 g".
- ◆ For future movements follow the same steps as before this way you will protect the sensor.
- ◆ Position the pieces on the balance as shown in the following picture: dust guard, draft ring, pan support and pan.
- ◆ If your balance is not used for a long period, please turn the display off, then unplug the AC adapter. This way you will save electricity and your equipment will have a longer life, also you will avoid accidents and an unnecessary use of energy.

Picture 4

### Security

- ◆ The balance has to be used by qualified personal so it gets its correct use.
- ◆ You should put the balance in a horizontal plane stable table, having a free space of at least 30 cm per side.
- ◆ Do not put the balance near any warm supply (burners, blowlamp...), do not expose directly to the sun, avoid vibrations and dust.
- ◆ During its functioning, do not use the equipment near fluids that can give off steams or produce explosive or inflammable mixtures.
- ◆ If you are not going to use the equipment during a long period of time please make sure that it is unplugged from the net
- ◆ Never dismantle the different pieces of the balance to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product's warranty
- ◆ Made under European regulations.

## 5. FUNCTIONS OF THE BALANCE

### Weighing


- ◆ To increase your balance's accuracy while weighing, please check that the place where you are going to take your weights has an average temperature of  $\pm$  20 °C and a relative humidity of 50%. Before using your balance, let it 10 or 15 minutes for its stabilizing.
- ◆ If your balance is used almost every day, just turn the display off by pressing *On/Off* (2.1) key with power connected. In long periods of no usage, keep the balance unplugged.

### Calibration

- ◆ The balance should be calibrated the first time it is used. It is necessary to recalibrate your balance when it is moved to another place or after a period of use (about 30 days).
- ◆ This balance has an internal weight for its calibration.
- ◆ Make sure that the level of the weight is on the position of not weighing (just as indicated in the picture 5.8.b).
- ◆ Turn the balance on by pressing *On/Off* (2.1) key with the pan empty.
- ◆ When you turn the balance on it will show 8.8.8.8.8.8.8.8. and once it is stabilized it will show 0.0000 g.

*Note:* If the display does not show 0.0000 g press *Tare* (2.3) until you can see in the screen the fixed value: 0.0000 g and the arrow shows that the balance is stabilized.

- ◆ Press "CAL" key (2.4) for the calibration.
- ◆ "CAL in" will be displayed on the screen.
- ◆ After that "CAL dn" appears.
- ◆ Pull the handle located on the right side of the balance towards you gently until it is stopped (5.8.a). While "CAL ..." is displayed, do nothing but wait for next display. It shows that you should take the calibration mass off the weighing system when "CAL up" appears.
- ◆ Push the same handle backwards gently until it is stopped (5.8.b).
- ◆ "CALL..." will be displayed and then seconds after "CALL End" appears.
- ◆ After some seconds 0.0000 g will be displayed. It indicates that the balance came into weighing mode and you can perform weighing now. The calibration has ended.
- ◆ It is necessary to check the result by using the internal calibration mass or an external calibration mass after calibration.

Picture 5

### Check using the internal calibration mass

- ◆ First you should learn the exact value of the internal calibration mass from the name plate in the rear of the balance.
- ◆ When the display shows 0.0000g, put the mass on the weighing pan and check that the difference between the value displayed and the value of the internal calibration mass is  $\pm$ 0.0001 g (maximum). If the result differs from  $\pm$ 0.0001 g, please repeat the calibration procedure.