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Thank you for choosing this equipment. We sincerely wish that you enjoy your Nahita mini centrifuge. We highly recommend looking after this equipment according to what is stated in this manual.

Nahita develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

**TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT**

Please bear in mind the following:

- ◆ This manual is inseparable from the Nahita mini centrifuge, so it should be available for all the users of this equipment.
- ◆ Follow all required applicable laws about safety and health (electrical security devices, inflammable, toxic and/or pathological products, etc.).
- ◆ You should carefully handle the centrifuge avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- ◆ Never dismantle the different pieces of the centrifuge to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the equipment.
- ◆ To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting Nahita Technical Assistance Department by email to (asistencia@auxilab.es)
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the equipment warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the centrifuge warranty.
- ◆ Accessories (including their loss), are not covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Nahita Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.

**ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.**

**INDEX OF LANGUAGES**

Spanish ..... 2-8  
 English ..... 9-15  
 French ..... 16-22

**INDEX OF CONTENTS**

1. USES OF THE INSTRUMENT ..... 10  
 2. DESCRIPTION ..... 10  
 3. TECHNICAL SPECIFICATIONS ..... 11  
 4. THEORETICAL FOUNDATIONS ..... 11  
 5. INSTALLATION / SETTING UP ..... 12  
 6. ACCESORIES..... 14  
 7. MAINTENANCE AND CLEANING ..... 15  
 CE CERTIFICATE ..... 23

**1. USES OF THE INSTRUMENT**

The Nahita Mini centrifuge is a small, silent and desktop centrifuge very easy-to-use and with a great stability. The model 2507/24 is only for use in the car even for external or field applications.

It is very versatile since it can be used with Eppendorf tubes of 2, 1.5, 0.5 and 0.2 mL, thus being ideal for using in techniques such as PCR, HPLC or microfiltration.

**2. DESCRIPTION**

- 1.1 Lid
- 1.2 Rotor
- 1.3 Lid opening button
- 1.4 Power on/off
- 1.5 Speed selecting button (Model 2507/14 and 2507/24)

emergency.

- ◆ Never try to repair the centrifuge by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the centrifuge (burns, hurts...) or damages in nearby equipments.
- ◆ In the event of breakdown please contact your distributor to overhaul through Nahita Technical Assistance Department.

Cleaning

- ◆ For the cleaning of the different parts of the equipment, we recommend the use of neutral products with a pH value not higher than 6.
- ◆ Never use scourers or substances that can grate for cleaning metallic parts such as stainless steel, aluminium, coatings, etc. as they damage the centrifuge and produce an early ageing of the equipment.
- ◆ Use a fluff-free cloth dampened with soaped water that does not contain abrasives.

**7. ACCESORIES**

Microtubes for centrifuges

<i>Code</i>	<i>Description</i>
40420020	Centrifuge microtube PP, 0.2 ML, B/1000
40420050	Centrifuge microtube PP, 0.5 ML, B/1000
40420150	Centrifuge microtube PP, 1.5 ML, B/500

**INSTRUCTIONS ON ENVIRONMENT PROTECTION**

At the end of its life cycle, please, do not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

**Security**


- ◆ The centrifuge must be used by previously qualified staff that knows how the equipment works thanks to the user manual.
- ◆ Nahita centrifuges are not suitable for treatment of flammable liquids, explosive samples or samples capable of forming explosive mixtures.
- ◆ When working with toxic or pathologic substances close hermetically the tubes and strictly observe the safety regulations corresponding to such substances.
- ◆ When you are not using the centrifuge for a long period of time please make sure it is unplugged in order to avoid possible accidents.
- ◆ It is essential to have the equipment switched off and unplugged from the net before cleaning, checking components or replacing any piece.
- ◆ Never try to repair the centrifuge by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the equipment (burns, hurts...).
- ◆ Try not to spill any liquid inside the centrifuge though it is properly insulated. In case you have any doubts do please immediately unplug the centrifuge (see Maintenance).
- ◆ Made under the European regulations for electrical security, electromagnetic compatibility and security on machines.

**6. MAINTENANCE AND CLEANING**

To get the best results and a higher duration of this equipment it is essential to follow the processes of use.

*Note: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.*

- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of this equipment.
- ◆ Use the appropriate tubes resistant to the R.C.F. used.
- ◆ Balance the centrifuge by using tubes with the same weight, shape and size and put them in geometrically symmetrical positions; if it is necessary, use tubes filled with water. This is absolutely necessary because an upset of the rotor may produce vibrations that debilitate the centrifuge.
- ◆ Never force the stoppage of the centrifuge, it is a risk of accident and debilitation and it may produce as well a mixing of the separated parts.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ The centrifuge is supplied with a Schuko standard wire. It has to be plugged to an earth connection and the socket should be handy and ready to unplug the equipment in case of


**3. TECHNICAL SPECIFICATIONS**

<b>Code</b>	52507001	52507014	52507024
<b>Model</b>	2507/1	2507/14	2507/24
<b>Capacity</b>			
<i>Angle rotor</i>		6x1.5/2.0 mL 6x0.5 mL* 6x0.2 mL*	
<i>Strip rotor</i>		2 tiras de 8 x 0.2 mL	
<b>Speed</b>	10000 rpm	4000/6000/10000 rpm	
<b>FCR</b>	5000 g	900/2000/5000 g	
<b>Decibels</b>		≤ 52 dB	
<b>Power supply</b>	CA 100-240 V, 50/60 Hz		DC12V (coche)
<b>Dimensions</b>		17x15x13 cm.	
<b>Weight</b>		~0.8 Kg	

\* Using the included adapters

Centrifugation is a separation technique based on the movement of the particles that are displaced to the distal end of the rotation axis depending on their shape and mass.

The Relative Centrifugal Force (R.C.F.) is the required force to separate particles. It is given as some value times that of gravity (g).

The R.C.F. is calculated with the following formula:

$$F.C.R.=1,118 \cdot 10^{-5} \cdot r \cdot n^2$$

$1.118 \cdot 10^{-5}$  is a constant.

$r$  = radius, the horizontal distance in centimetres from the rotation axis to the bottom of the sample tube.

$n$  = rotation speed expressed in revolutions per minute (r.p.m.).

## 5. INSTALLATION / SETTING UP

### Preliminary inspection

Unwrap the centrifuge; take off the involving plastic and take off the polyspan protection in which it comes fitted. Check manually that the rotor spins freely without any obstacle and make sure the equipment does not present any damage because of the shipment. In case it presents any damage, tell it immediately to your transport agent or dealer so that they can make the claims in the correct time limit.

Please keep the original wrapping; you will always need it for returns enclosed with all the accessories supplied.

Please check that all the accessories are enclosed with the equipment:

- 6 Adapters for Eppendorf-like tubes of 0.5 mL (inner Ø of the adapter 8 mm)
- 6 Adapters for Eppendorf-like tubes of 0.2 mL (inner Ø of the adapter 6 mm)
- 1 Strip rotor for microtubes of 0.2 mL
- Schuko wire
- User's manual
- Warranty certificate

**We will only accept any equipment return within 15 days after delivery and provided it comes in its original wrapping.**

### Installation

- ◆ Before using this instrument, it is convenient for you to familiarize with its components and basic essentials, as well as with its control functions.
- ◆ Please put the centrifuge on top of a horizontal, plane and stable table making a free space of at least 30 cm per side. During operation, people and dangerous material such as flammable or infectious liquids must be out of this area. Do not put the equipment near any warm supply (burners, blowlamps...), nor expose it directly to the sun, etc.
- ◆ The centrifuge is supplied with a Schuko standard wire 220-240 V, 50-60 Hz. Please insert the wire that feeds the AC electric current in the base of current provided with

earth.

**Neither the manufacturer nor the distributor will assume any responsibility for the damages produced to the equipment during its installation or damages to people suffered by the improper use of the electric connection. The tension should be 220 V, 50 Hz  $\pm$  10%.**

If you are not using the centrifuge for a long period of time please make sure it is disconnected from the net and protected from dust (this way you will avoid accidents and will extend its working-life).


**PLEASE READ THOROUGHLY THE INSTRUCTIONS BEFORE CONNECTING AND OPERATING WITH THIS EQUIPMENT WITH THE AIM OF ACHIEVING THE HIGHEST FEATURES AND THE MAXIMUM DURABILITY OF THE INSTRUMENT.**

### Setting up


- ◆ Open the lid by pressing the corresponding button (1.3)

In case of using tubes of 0.5 or 0.2 mL put the corresponding adapters into the holes of the rotor.


- ◆ Bear in mind that the tubes must be perfectly balanced. Never centrifuge a unique tube; in case of having only one sample, fill another tube with an equivalent volume of water to balance the rotor.
- ◆ Make sure that all the tubes are well closed and insert them completely into the holes. The tubes must be placed opposite as shown in the picture.


2 tubes

3 tubes

4 tubes

6 tubes

- ◆ Select the centrifugation speed by the speed selecting button (1.5) (Model 2507/14 and 2507/24)
- ◆ Close the lid and press power on/off (1.4).

**NOTE: If the lid is not completely closed, the centrifuge will not start operating.**

- ◆ Once the process is finished, turn the centrifuge off and wait until the rotor completely stops to open the lid and take out the tubes.