## ANNEX I: CE CERTIFICATE

AUXILAB S.L.

**DECLARATION OF CONFORMITY**  
**PORTABLE pH METER 902/4 of Auxilab, S.L for the**  
**Directive of Machines (89/392/CEE modified) and the regulations**  
**adopted for their transposition**

NAME OF THE MANUFACTURER / IMPORTER:  
**AUXILAB, S.L.**

ADDRESS: **Polígono Morea Norte, 8**  
**31191 Beriáin (Navarra)**

WE STATE THAT:

**ph metro NAHITA 902/4**  
**Model 50902040**

Are designed and manufactured according to:

- ♦ Directive 89/392/CEE, including the modifications and national regulations that transpose them.
- ♦ Directive 89/336/CEE modified over electromagnetic compatibility
- ♦ And that the following harmonized rules have applied (or part of them):  
 UNE 292-1/-2/-2/A1, UNE-EN 1050, UNE-EN 614-1, UNE-EN 1037, UNE-EN 1088, UNE-EN 547, UNE-EN 953, UNE-EN 294, UNE-EN 418, UNE-EN 894-1, UNE-EN 894-2, UNE-EN 954-1, UNE-EN 60204-1, UNE 61010-1/A2, UNE-EN61010-2-051..

BERIAIN 16th July 2007

Signed by: **ALFONSO AINCIBURU SANZ**  
**DIRECTOR/MANAGER**

Polígono Morea Norte, 8 31191 Beriain (Navarra) - Spain. Tel. 948 310 513 Fax 948 312 071  
Internet: [www.auxilab.es](http://www.auxilab.es) · Email: [correo@auxilab.es](mailto:correo@auxilab.es)

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*This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.*

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### Glass body electrode

Code	Working temp	Diaphragm	Membrane	Electrolyte
90900651	5-60 °C	Ceramic	Bulb	Liquid, refillable
90900652	5-90 °C	Ceramic	Bulb	Liquid, refillable
90900953	5-60 °C	Ground glass	Bulb	Gel, sealed
90900955	5-90 °C	Ground glass	Bulb	Liquid, refillable

## 5. MAINTENANCE AND CLEANING

To get the best results and a higher duration of the pH meter it is essential to follow the processes of use.

*Note: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.*

- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of this equipment.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ Never try to repair the pH meter by yourself, since you will lose the warranty and may provoke damages to the general operating system, as well as injuries to the people that usually handle the pH meter (burns, hurts...).
- ◆ Maintenance of the electrode. Rinse the pH electrode in clean water after each measurement; for a thorough cleaning, clean the electrode with a cloth and a solution of HCl 0.1 M. After cleaning the electrode, calibrate the equipment.
- ◆ Do not use the light more than necessary to extend the life of batteries.
- ◆ If LOBAT is displayed, change the batteries.
- ◆ For long periods of no use, remove the batteries.
- ◆ Do not expose the equipment directly to the sunlight, since the LCD could suffer damages and produce an early ageing of the equipment.
- ◆ In the event of breakdown please contact your distributor to overhaul through Nahita Technical Assistance Department.

### Cleaning

- ◆ Never use scourers or substances that can grate for cleaning the case of the equipment as they damage the pH meter and produce an early ageing of the equipment.
- ◆ Use a fluff-free cloth dampened with soaped water that does not contain abrasives.

### INSTRUCTIONS ON ENVIRONMENT PROTECTION

*At the end of its life cycle, please, do not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.*

*The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.*

*Please inquire at the community administration for the authorized disposal location.*

- ◆ When the displayed value is stable press ENTER (3).
- ◆ Press MODE (2) and STD2 will be displayed at the lower part of the screen. This is the second calibration point.

◆ Clean the electrode with distilled water and submerge it in the buffer solution.

*NOTE: For the second calibration point use the solution with pH 4.01 if working in an acidic pH range or the solution with 10.1 if working in a basic pH range.*

- ◆ When the displayed value is stable press ENTER (3). The equipment is calibrated.
- ◆ To go to the H measuring range press MODE (2); MEAS will be displayed. The pH meter is ready to operate.

**pH and mV measurements**

*NOTE: Make sure the equipment is calibrated before performing any measurement.*

- ◆ Turn the ph meter on by pressing ON/OFF (1).
- ◆ Select the pH or mV measuring mode by pressing MODE (2). MEAS will be displayed at the lower part of the screen.
- ◆ Clean the electrode with distilled water and dry it with a filter or cellulose paper.
- ◆ Submerge the electrode 4 cm into the sample to be analyzed and stir it with the electrode. Wait until the reading is stable. pH or mV value will be displayed at the middle and below is shown the temperature value previously introduced.

*NOTE: If several samples are going to be analyzed successively it is recommended to clean the pH electrode with distilled water to avoid cross contamination.*

**Security**

- ◆ The pH meter must be used by previously qualified staff that knows how the equipment works thanks to the user manual.
- ◆ It is essential to have the equipment switched off before cleaning, checking components or replacing any piece.
- ◆ Never try to repair the pH meter by yourself, since you will lose the warranty and may provoke damages to the general operating system, as well as injuries to the people that usually handle the equipment (burns, hurts...).
- ◆ Try not to spill any liquid on the control panel, though it is properly insulated.
- ◆ Made under the European regulations for electrical security, electromagnetic compatibility and security on machines.

**Accessories**

We are provided with a wide range of combined ph electrodes to adapt to any application in chemistry pharmaceuticals, education, etc. All the electrodes present a BNC connector and a pH range from 0 to 14.

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◆ The Nahita pH meter 902/4 operates with 2 alkaline batteries type AA. Open the batteries compartment located at the back side of the equipment and put the batteries according to their polarity. Put again the lid of the compartment.

*NOTE: Remove the batteries if you are not using the instrument for a long period of time.*

◆ Connection of the pH electrode. Connect the pH electrode (8) to the connector (7) situated at the upper part of the equipment.

◆ Turn the pH meter on by pressing ON/OFF (1). After a few seconds the equipment will be ready to operate.

#### Working temperature adjustment

Nahita pH meter 902/4 is provided with a temperature compensation mechanism by entering data using the keypad. For this purpose:

- ◆ In the measuring mode press MODE (2) until the °C symbol starts blinking.
- ◆ Introduce the exact value of the temperature of the solution to be analyzed by using ▲ (4) and ▼ (5).
- ◆ Press ENTER (3) to save the temperature value.

#### Calibration of the pH meter

*Note: It is recommended to calibrate de pH meter before starting using the equipment and whenever observing variations in measurements, after changing the electrode and after a prolonged use with solutions of pH<2 or pH>12, organic solutions or fluorine containing solutions.*

#### Preparation of the standard buffer solutions:

The pH meter is supplied with the following standard buffer solutions: pH 4.01 (25 °C), pH 7.00 (25 °C), pH 10.01 (25 °C).

- ◆ Pour the content of one of the packet of buffer solution into a clean beaker and add 100 mL of distilled water.
- ◆ Dissolve well the powder; if it is necessary add more water.
- ◆ Once all the content is dissolved, pour the solution into a 250 mL volumetric flask and add distilled water until it reaches the neck of the flask.

#### Calibration:

- ◆ Connect the electrode to the equipment.
- ◆ Turn the pH meter on by pressing ON/OFF (1). The equipment will enter in the mV measuring mode.
- ◆ Press MODE (2) to make the temperature adjustment. When the symbol °C blinks introduce the temperature of the standard buffer solution by using ▲ (4) and ▼ (5). Save the value by pressing ENTER (3).
- ◆ Press MODE (2) again and STD1 will be displayed at the lower part of the screen. This is the first calibration point.

*NOTE: If the equipment does not need to be calibrated, press MODE (2) two more times and the equipment will enter the pH measuring mode.*

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*NOTE: For the first calibration point use the solution with pH 7.00 (neutral).*

### 3. TECHNICAL SPECIFICATIONS

<b>Code</b>	<b>50902040</b>
pH range	0.0 a 14.0
mV range	0 a $\pm 1400$ mV
pH resolution	0.01
mV resolution	1 mV
pH accuracy	$\pm 0.03$
mV accuracy	$\pm 0.2\%$
Temperature compensation	Entering data using the keypad
ATC range	0-60 °C
Calibration	1 or 2 points
Power	2 batteries 1.5 V, tipo AA

### 4. INSTALLATION / SETTING UP

#### Preliminary inspection

Unwrap the pH meter and, without turning the equipment on, make sure that it does not present any damage because of the shipment. In case it presents any damage, tell it immediately to your transport agent or dealer so that they can make the claims in the correct time limit.

Please keep the original wrapping; you will always need it for returns enclosed with all the accessories supplied.

Please check that all the accessories are enclosed with the equipment:

- pH electrode
- Standard buffer solutions pH 4.01, pH 7.00, pH 10.01
- 2 batteries 1.5 V, type AA
- Warranty certificate
- User's manual

**We will only accept any equipment return within 15 days after delivery and provided it comes in its original wrapping.**

#### Installation

Before using the pH meter, it is convenient for you to familiarize with its components and basic essentials, as well as with its control functions.

**PLEASE READ THOROUGHLY THE INSTRUCTIONS BEFORE CONNECTING AND OPERATING WITH THIS EQUIPMENT WITH THE AIM OF ACHIEVING THE HIGHEST FEATURES AND THE MAXIMUM DURABILITY OF THE INSTRUMENT.**

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
**1. USES OF THE INSTRUMENT**

The Nahita pH meter 902/4 is a light, compact and easy-to-use portable instrument ideal for in situ pH determinations. It presents an easy-to-read liquid crystal display that can be illuminated to facilitate readings in dark places and a damp-proofed membrane keypad with 5 functional keys to access to all the instrument possibilities.

The equipment, supplied with a glass body combined electrode, can perform pH and voltage readings of different solutions. To adapt to any application, the equipment can also be used with other type of electrodes provided they present a BNC connector. The Nahita pH meter 902/4 also presents a temperature compensation mechanism by entering data using the keypad.

**2. DESCRIPTION**

1. ON/OFF
2. MODE: selection of the measuring, calibration and temperature adjustment modes.
3. ENTER: parameter confirmation
4. ▲ : parameter increase
5. ▼ : parameter decrease
6. Light button
7. Electrode connector
8. pH electrode

Thank you for choosing this equipment. We sincerely wish that you enjoy your Nahita portable pH meter 902/4. We highly recommend looking after this equipment according to what is stated in this manual.

Nahita develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

**TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.**

Please bear in mind the following:

- ◆ This manual is inseparable from the Nahita portable pH meter 902/4, so it should be available for all the users of this equipment.
- ◆ You should carefully handle the equipment avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- ◆ Never dismantle the different pieces of the equipment to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the pH meter.
- ◆ To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting Nahita Technical Assistance Department by email to asistencia@auxilab.es or by telephone: +34 807 117 040 (0.30 Euros/min).
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the equipment warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the warranty.
- ◆ Accessories, including their loss are not covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Nahita Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.

**ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.**

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Referencia	Tª trabajo	Diafragma	Membrana	Electrolito
90900218	5-60 °C	Fibra	Esférica	Gel, no rellenable
90900219	5-60 °C	Fibra	Esférica	Líquido, rellenable
90900900	5-60 °C	Anillo PTFE	Esférica	Gel, no rellenable
90900901	5-60 °C	Anillo PTFE	Plana	Gel, no rellenable

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