# REFRACTÓMETRO ABBE ABBE REFRACTOMETER

**Zuzi**

**Modelo / Model 315**

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*This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.*

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**ANNEXE III: CE CERTIFICATE**

**AUXILAB S.L.**

**CE DECLARATION OF CONFORMITY**  
**ZUZI REFRACTOMETER SERIES 315 of AUXILAB, S.L. for the Directive of Machines (89/392/CEE modified) and the regulations adopted for their transposition**

NAME OF THE MANUFACTURER/IMPORTER:  
**AUXILAB, S.L.**

ADDRESS: **Polígono Morea Norte, 8  
 31191 Beriáin (Navarra)**

WE STATE THAT:  
**ZUZI REFRACTOMETER SERIES 315  
 MODELS 50315010 AND 50315000**

Are designed and manufactured according to:

- ◆ Directive 89/392/CEE, including the modifications and the national regulations that transpose them.
  - ◆ Directive 73/23/CEE modified over the electric security
  - ◆ Directive 89/336/CEE modified over the electromagnetic compatibility
- And that the following harmonized rules have been applied (or part of them):

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BERIAIN 23rd June 2004

Signed by: ALFONSO AINCIBURU SANZ  
 DIRECTOR/MANAGER

Polígono Morea Norte, 8 31191 Beriain (Navarra) - Spain. Tel. 948 310 513 Fax 948 312 071  
Internet: www.auxilab.es · Email: correo@auxilab.es

**ANNEXE II**

**Theoretical basis. Refractometry**


When a light ray gets to the boundary between different media, part of the light is reflected and part is refracted (it goes in the second medium), spreading the latter with an angle and speed different from the ones of the first medium. These two magnitudes can be distinguished by the refractive index, which is a characteristic value for each substance.

The refractive index is related to the mass, the charge and the number of particles of the substance through which the luminous radiation is transmitted. There are different types of refractometers, but most of them follow the same principle.

The first principle of a refractometer is agreement with the following Law of Refraction: being  $n_1$  and  $n_2$  the refractive indexes of the medium at both sides of the boundary between different media, where  $\alpha$ ' is the incident angle and  $\alpha''$  is the angle of refraction.

When a light ray goes from an optically dense medium to another medium optically lighter, the incident ray is narrower than the refraction angle. When we change the first so that the second rises  $90^\circ$ , the incident angle is called *critical angle*. The principle of this refractometer to set the refractive index is based on knowing the critical angle.

When the light ray passes over the surface  $AB$  from any angle (picture 4), its refraction angle is always bigger than  $i$ . If we observe the emerging beam through a telescope, the field of vision appears divided in two halves, one dark and another light. If we place the delimiting line of the two halves just in the middle we obtain the critical angle.  $ABCD$  is the refractive prism whose refractive index is  $n_2$ . Over the surface  $AB$  we place the substance to be measured (liquid or transparent solid) whose refractive index is  $n_1$ .

Picture 4

$$\text{Law of Refraction: } n_1 \times \text{sen} \alpha' = n_2 \times \text{sen} \alpha''$$

Thanks to the Law of Refraction we obtain the following equation:

$$\begin{aligned} n_1 \times \text{sen} 90^\circ &= n_2 \times \text{sen} \alpha \\ n_2 \times \text{sen} \beta &= \text{sen} i \\ \phi &= \alpha + \beta \\ \alpha &= \phi - \beta \end{aligned}$$

Transferring this to the first equation we obtain that:

$$\begin{aligned} n_1 &= n_2 \times \text{sen} (\alpha - \beta) = (\text{sen} \phi \times \text{cos} \beta - \text{cos} \phi \times \text{sen} \beta) \\ n_2^2 \times \text{sen}^2 \beta &= \text{sen}^2 i \\ n_2^2 \times (1 - \text{cos}^2 \beta) &= \text{sen}^2 i \\ n_2^2 - n_2^2 \times \text{cos}^2 \beta &= \text{sen}^2 i \end{aligned}$$

$$\text{cos} \beta = \sqrt{\frac{n_2^2 - \text{sen}^2 i}{n_2^2}}$$

$$n_1 = \text{sen} \phi \times \sqrt{n_2^2 - \text{sen}^2 i} - \text{cos} \phi - \text{sen} i$$

The refractive angle  $\phi$  and the refractive index  $n_2$  of the prism are known, so as we know the critical angle  $i$ , we can obtain the refractive index of the sample.

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

**ANNEXE I**

Correction chart of the refractive index with the temperature of the distillate water

Temperature °	Refractive Index
10	1.33369
11	1.33364
12	1.33358
13	1.33352
14	1.33346
15	1.33339
16	1.33331
17	1.33324
18	1.33316
19	1.33307
20	1.33299
21	1.33290
22	1.33280
23	1.33271
24	1.33261
25	1.33250
26	1.33240
27	1.33229
28	1.33217
29	1.33206
30	1.33194
31	1.33182
32	1.33170
33	1.33157
34	1.33144
35	1.33131
36	1.33117
37	1.33104
38	1.33090
39	1.33075
40	1.33061

- ◆ Avoid spilling water on the control panel.
- ◆ During its functioning keep the equipment away from dangerous materials such as inflammable material or liquids.
- ◆ Do not put the refractometer near any warm supply (burners, blowlamp...), do not expose directly to the sun.
- ◆ If you find out that there is degradation of the cement that joins the inferior prism to the carcase, disconnect the equipment immediately from the net and get in touch with the ZuziTechnical Service.
- ◆ If you are not going to use the equipment during a long period of time please make sure that it is unplugged from the net.
- ◆ Never dismantle the different pieces of the refractometer to repair yourself, since it could produce a defective use of the whole equipment and a loss of the product's warranty. For the overhaul of the equipment please contact the Zuzi Technical Service.
- ◆ Fabricated under the European regulations.

## 5. MAINTENANCE AND CLEANING

*To get the best results and a higher duration of this equipment it is necessary to follow the following recommendations.*

*Note:* All the processes of use mentioned before will not have any value unless you keep a continued and careful maintenance.

- ◆ Please follow the processes of use of this manual
- ◆ This manual should be available for all users of this equipment
- ◆ Do always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ After every use the refraction prisms should be cleaned with a cloth that does not give off fluff, damped with water or with neutral products (water with soap without abrasives).
- ◆ Solid impurities should not exist in the liquid samples and the solid samples should have polished and plane faces to avoid scratching the prism.
- ◆ If the lamp gets blown, change it by taking off the condensed lateral lens and then putting the lens on again.
- ◆ When you are not going to use the refractometer for a long period keep it in its box with a desiccant substance (Silica Gel).
- ◆ Avoid damages and rubbing between the prisms for an optimum functioning of the refractometer; it can be damaged by samples that contain crystalline or big fragments.
- ◆ If there is any damage do not try to repair the equipment yourself, since it could produce a defective use of the whole equipment and a loss of the product's warranty. It can also cause damages to the electric system.
- ◆ For all your doubts please contact Zuzi's official Service.

### Cleaning

- ◆ For the cleaning of the different pieces of the equipment, we recommend the use of neutral products that do not contain abrasive agents, e.g. some drops of cleaning glass soap reducing it in water.
- ◆ Never pour the dissolution over the equipment, always damp a soft cloth in the soap and water.

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◆ When you press the button *Read*, the "00000" will disappear from the screen, you will see "Waiting" and we will obtain the refractive index of the sample. You can make the choice of the reading mode once the first reading is done immediately after turning on the equipment or once it is made.

◆ When the refractive index appears on screen, it can be changed.

*Note:* Once the *Read* button is pressed, wait until the measure appears in screen; if you press continuously this button it will cause an obstruction in the system leaving the equipment useless.

◆ In the screen all the measure units are shown but, if we want to make a reading in Brix<sup>o</sup> with adjustment to temperature or we want to know the temperature at which the measures are being made, we will press the correct button situated in the control panel (see picture 2). The digital screen will show the principal form of the unit we are working with. When we want to work in Brix<sup>o</sup>, press the "BRIX" button and, as the principal measure, you will see Brix<sup>o</sup>. The same procedure is followed to work with the refractive index "nD" and in brix with the temperature adjustment "BX-TC", if you wish to see it as the principal unit.

*Note:* Each time we change the sample we have to press *Read* (2.3); this way the measurement will take place.

◆ If we want to measure the temperature of the sample we will press *Temp*. This measure can be made in every moment except when in the screen we can read "Waiting", after the button *Read* has been pressed.

◆ If you keep the equipment in a different place of that where you make the measures, bear in mind that it has to adapt to the environment's temperature of the working place; if not, it will be giving us the temperature at which the equipment is, depending where it has been kept.

◆ After making the appropriate measurements dissolve the sample with water (when the sample is a sweet dissolution) or use the suitable dissolvent for each sample; after that clean the prism carefully with a special cloth that does not give off fluff and does not scratch, damped with neutral products.

### Calibration

◆ The refractometers should be calibrated periodically, for this we use distillate water for the *test piece*.

◆ If you proceed in the correct way, as stated in the previous chapter, you will obtain the refractive index. If it does not coincide with the known value of the distillate water at 20° C (1'33299) or with the test piece at 20° C (engraved in the upper part), we will insert the screwdriver that is supplied in the small hole located in the dispersion correction handle.

◆ By turning it we will make the crossed lines go up or down until we make the reading fit with the known values. You will have to check when moving the crossed lines that the value read by the equipment is the value of reference

### Security

◆ The refractometers should be used by qualified staff, that is aware of the use of the equipment by the instruction manual.

◆ Please put the refractometer on top of a horizontal, plane and stable table making a free space of at least 30 cm per side.

**Installation**

- ◆ Before using this instrument, it is convenient for you to familiarize with its components and basic essentials.

**PLEASE READ THOROUGHLY THE INSTRUCTIONS BEFORE CONNECTING AND OPERATING WITH THIS EQUIPMENT.**


- ◆ Please put the refractometer on top of a horizontal, plane, stable table making a free space at least at 30 cm per side. Do not put the refractometer near any warm supply (burners, blowlamp...), do not expose directly to the sun, etc. During its functioning, do not use the equipment near fluids that can give off steams or produce explosive or inflammable mixtures.
- ◆ The refractometer is sumministrated with a Schuko standard wire.
- ◆ Please insert the wire that feeds the A/C in the base of current 220V  $\pm$ 10% 50-60Hz provided with earth wire..

*The manufacturer and the distributor will not assume any responsibility for the damages produced to the equipment in its installing or damages to persons suffered by the wrong use of the electric connection. The tension should be 220V  $\pm$ 10% 50-60Hz.*

- ◆ If you are not using the refractometer during a long period of time please make sure it is disconnected from the net.

**Setting up**

- ◆ Connect the refractometer to the current net.
- ◆ Press *Power* (2.2); it will light up (1.8) at the same time as the digital screen (2.1) shows "Welcome"; immediately after you will read "Waiting". The screen will show all the units that can be made, but by defect, as the principal until it will show en *nD*.
- ◆ Open the unit (1.6) where the prism are located and clean them carefully with a special cloth that will not harm or damage the prism -free from dust, it has to be impregnated with water and with neutral products-
- ◆ Put the sample on the surface of the refractive prism.
- ◆ If the sample is liquid, put two drops over the working surface in the lower prism, positionate the superior prism and now you can make the measures.
- ◆ If the sample is solid, it should have a plane and polished face over which you will, once it is clean, put some drops of Brome Naphthalene. We will put the sample over the prism making sure that there is a good contact of both plane surfaces and a good distribution of the brome naphthalene. It is not necessary to close with the superior prism to make measures of solids.
- ◆ Turn the illuminator arm (2.7) until you get a good illumination, where we can see more contrast between the clear and dark zones.
- ◆ Turn the selector handle (1.5) until you make the divisor line coincide clear/dark in the crossed lines. With the handle of the dispersion correction (1.2) we will make the most possibly sharp/clear the divisor line.
- ◆ Confirm that this line coincides exactly with the centre of the crossed lines and press the button *Read* (2.3) (see picture 3).

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**3. TECHNICAL SPECIFICATIONS**

<b>Measuring range</b>	
Refractive Index	1'3000 to 1'7000
<sup>o</sup> Brix with T <sup>e</sup> correction	0 to 95%
<sup>o</sup> Brix	0 to 95%
<b>Measure's precision</b>	
Refractive Index	± 0'0002
Brix with T <sup>e</sup> correction	± 0'1%
<sup>o</sup> Brix	± 0'1%
<b>Temperature</b>	
Display range	0 to 50 °C
Brix correction versus T <sup>e</sup>	15-45°C
<b>Exterior measures</b>	
	330 x 180 x 380 mm
<b>Weight</b>	
	10 Kg
<b>Electric supply</b>	
	220V±10% 50-60Hz
<b>Fuse</b>	
	250V 1A

*Note:* RS232 parameters: Bauds 2400, Data 8, Stop 1. The flow control is hardware and without parity.

**4. INSTALLATION / SETTING UP*****Preliminary inspection***

- ◆ Unwrap the Refractometer, take of the involving plastic and take of the poliespan protection in which it comes fitted, please make sure that it does not present any damage because of the shipment. If it does please communicate to your transport agent or sumministrador immediately, so that they can make the claims in the correct time limit.
- ◆ Please keep the original wrapping; you will always need it for returns enclosed with all the accessories sumministrated.
- ◆ Please check that all the accessories are enclosed with the equipment:
  - Refractometer
  - The following articles are enclosed in the case:
    - a) Test piece (1,51685).
    - b) Wrapping of mono brome naphthalene.
    - c) Two lamps 6,3V and 2,5W.
    - d) Screwdriver.
    - e) Standard Schuko cable.
    - f) Manual of use.
    - g) Warranty Certificate.

***We will not accept any return of equipments in warranty unless they are in their original wrapping.***

**1. APLICATIONS OF THE EQUIPMENT**

The *Refractometer* is used to measure the refractive index (nD) of transparent or translucent liquids and solids, being more used the transparent liquids.

The Zuzi Refractometer 315 RS232 has a circuit in its point of view for the thermostatzation of the samples at the temperature we need for working. For the correction of the refractive index with the temperature please go through the enclosed tables.


The *refractive index* is an important optic constant in the substances, and it can be used to determine the optic behaviour, purity, concentration, dispersion, etc.

It is possible to calculate the sugar percentage contained it the sacarose dissolutions in a range between 0 and 95 Brix°, what matches with a refractive index of 1.333 to 1.531.

This equipment has an extensive range of practical applications, being an essential tool in the petrol industry, of oil and grease, paintings, pharmacy, food, cosmetic, chemistry, and sugar, also in universities and centres of investigation and teaching.


**2. DESCRIPTION**

Picture 1 - Parts of refractometer

- 1.1. Eyepiece
- 1.2. Correction and dispersion control
- 1.3. Calibration screw
- 1.4. Membrane and screen
- 1.5. Regulation control
- 1.6. Thermostatization prism.
- 1.7. Articulated illuminator with condensing lens.
- 1.8. Power point (Schuko base).
- 1.9. Output RS 232.

Picture 2 - Membrane keyboard

- 2.1. Digital screen.
- 2.2. Power. On of interruptor of the equipment.
- 2.3. Read. Lecture button.
- 2.4. Temp. It shows the temperature.
- 2.5. BX-TC. It shows the Brix° with the temperature correction.
- 2.6. nD. It shows the reading of the Refractive Index.
- 2.7. BRIX. It shows the reading in Brix°.

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Thank you for choosing this equipment. We sincerely wish that you enjoy working with your Zuzi Refractometer model 315 RS232; we highly recommend looking after this equipment according to what is stated in this manual.

Zuzidevelops its products being remarkable in the ergonomics and security of its user. The incorrect use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, please read the point of maintenance, where we expose the security notes.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

**TO GET THE BEST RESULTS AND HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.**

Please bear in mind the following:

- ♦ This manual should be available for all the users of these equipments
- ♦ We should manipulate carefully avoiding the sudden movements, knocks, free fall of heavy objects, and manipulation of sharp objects.
- ♦ Never dismantle the different pieces of the Refractometer to repair yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty.
- ♦ To avoid fire o electric discharges, avoid humid environments.
- ♦ If you have any doubt about setting up, installation, functioning do not hesitate in contacting your distributor or our ZUZI Technical Assistant Department by email asistencia@auxilab.es
- ♦ The manufacturer reserves the right to modify or improve the manual or equipment.
- ♦ Please do not forget filling the Warranty Certificate and send it before 15 days after the date of purchase.

**ATTENTION!!**

**IF AN EQUIPMENT IS NOT PROPERLY CLEAN AND DISINFECTED IT WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.**

**INDEX OF LANGUAGES**

Spanish .....	1-10
English.....	11-19

**INDEX OF CONTENTS**

1. APLICACIONES OF THE EQUIPMENT .....	12
2. DESCRIPTION .....	12
3. TECHNICAL SPECIFICATIONS .....	13
4. INSTALLATION/SETTING UP .....	13
5. MANTEINANCE AND CLEANING.....	16
ANNEXE I: CORRECTION CHART OF REFRACTIVE INDEX WITH T <sup>a</sup> .....	17
ANNEXE II: THEORETIC FOUNDATIONS.....	18
ANNEXE III: CE CERTIFICATE.....	19