AUXILAB S.L.

Material de laboratorio  
Laboratory supplies

DECLARACIÓN DE CONFORMIDAD CE  
CE DECLARATION OF CONFORMITY  
DÉCLARATION DE CONFORMITÉ CE

El fabricante | The manufacturer | Le fournisseur:

AUXILAB, S.L.

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Fdo: Alfonso Ainciburu Sanz  
DIRECTOR | GERENTE

A handwritten signature in black ink, appearing to be 'Alfonso Ainciburu Sanz', written over a white background.

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The logo for Nahita, featuring the brand name in a stylized, bold, black font.

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*This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.*

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Thank you for choosing this equipment. We sincerely wish you enjoy the Rotary Evaporator Model 9200/9300. We highly recommend looking after this equipment according to what is stated in this manual.

Nahita develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its users.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

**TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.**

Please bear in mind the following:

- ◆ This manual should be inseparable from the rotary evaporator model 9200/9300. This manual should be available for all the users of this equipments
- ◆ We should manipulate carefully avoiding the sudden movements, knocks, free fall of heavy objects, manipulation of sharp objects, please avoid spilling liquids inside the equipment.
- ◆ Never dismantle the different pieces of the rotary evaporator to repair yourself, since it could produce a defective use of the whole equipment, damage to the users and a loss of the product warranty.
- ◆ To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting Nahita Technical Assistance Department by email to asistencia@auxilab.
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the equipment warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the rotary evaporator's warranty.
- ◆ Neither fuses (1A) nor accessories (including their loss), are covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Nahita Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Please do not forget filling the warranty certificate and send it before 15 days after the date of purchase.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.

## 6. RESOLUTION PROBLEMS

<i>PROBLEM</i>	<i>CAUSE</i>	<i>SOLUTION</i>
No work	No power supply	1. Check the power supply base 2. Check connetions 3. Check fuses of 5A 4. Press switch on button
Motor not turn	insufficient voltage (flak excess weight)	Reduce rpm
Error "LLL" in screen SV	Do not switch off heater	1. Check connection T sensor 2. Check cables
Work properly but water bath is not heating	No llega tensión a la resistencia calefactora Fusible de 10A roto	1. Checks the cable connection 2. Change fuse 10A

## 7. ACCESSORIES

<i>Code</i>	<i>Description</i>
99200012	Oil condenser. 29/32
99000112	Flask 1 L, 29/32
99000122	Flask 2 L, 29/32
99000132	Flask 3 L, 29/32
99000030	Pipette with stopcock 24/29
99200010	Vapour ducting piece 29/32
63226229	Clip 29/32
59001010	Vacuum pump 10L/min
59001025	Vacuum pump 25L/min

### INSTRUCTIONS ON ENVIRONMENT PROTECTION

-At the end of its life cycle, please, do not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.

-The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

-Please inquire at the community administration for the authorized disposal location.

### Security

The rotary evaporator must be used by previously qualified staffs that know how the equipment works thanks to the user manual.

- ◆ You should put the equipment in a horizontal plane stable table, having a safety area of at least 30 cm per side.
- ◆ Do not place the equipment near any warm supply (burners, blowlamps, etc), nor expose it directly to the sun. Avoid vibrations, dust and dry environments.
- ◆ During its functioning dangerous materials such as flammable or pathological substances must be out of the safety area.
- ◆ When you are not using the rotary evaporator for a long period of time please make sure it is unplugged in order to avoid possible accidents.
- ◆ It is essential to have the equipment switched off and unplugged from the net before cleaning, checking components or replacing any piece (e.g. replacement of a fuse).
- ◆ Never try to repair the equipment by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the equipment (burns, hurts...).
- ◆ Try not to spill any liquid on the control panel, though it is properly insulated. In case you have any doubts do please immediately unplug the evaporator (see Maintenance).
- ◆ Made under the European regulations for electrical security, electromagnetic compatibility and security on machines.

## 5. MAINTENANCE AND CLEANING

To get the best results and a higher duration of this equipment it is essential to follow the processes of use.

*Note: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.*

- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of this equipment.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ In case of breakdown please contact your distributor to overhaul through Nahita Technical Assistance Department.

### Cleaning

- ◆ Never use scourers or substances that can grate for cleaning metallic parts such as stainless steel, aluminium, coatings, etc. as they damage the evaporator and produce an early ageing of the equipment.
- ◆ Use a fluff-free cloth dampened with soaped water that does not contain abrasives.
- ◆ When you finish, we recommend to evaporate distill water for cleaning the equipment.

**ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.**

## INDEX OF LANGUAGES

Spanish .....	2-13
English .....	14-23
French.....	24-34

## INDEX OF CONTENTS

1. USES OF THE INSTRUMENT .....	15
2. DESCRIPTION .....	16
3. TECHNICAL SPECIFICATIONS .....	17
4. INSTALLATION/SETTING UP .....	18
5. MAINTENANCE AND CLEANING .....	22
6. RESOLUTION PROBLEM .....	23
7. ACCESSORIES.....	23
CE CERTIFICATE .....	36

## 1. USES OF THE INSTRUMENT

Nahita rotary evaporators are essential equipments in fields such as Biochemistry, Organic Chemistry, Pharmacy, etc. for all the tasks where you need either a quick distillation or a sample concentration.


The increase of the evaporation surface with the application of heat on its rotating flask by water bath and an additional vacuum system (not included), achieve a significant increase in evaporation rate comparatively with traditional methods. Its rotating flask generates, a high heat transfer, prevents the localized overheating and mixing content simultaneously.


## 2. DESCRIPTION

- 1.1 Main body
- 1.2 Digital display
- 1.3 Water Bath
- 1.4 inclination angle Button
- 1.5 height adjustment Button

- 1.6 Evaporating flask
- 1.7 Coil condenser
- 1.8 Pipette flow sample
- 1.9 Receiving flask

- 2.1 Power supply (Input 220V)
- 2.2 Power supply (Output 220V)
- 2.3 Sensor connection T
- 2.4 Bath connection (Input 220V)
- 2.5 Fuses

- ◆ Place the evaporating flask and screw the device subsection (5.4) onto the vapor ducting piece
- ◆ Select evaporation flask and introduce the sample to be processed.
- ◆ Insert the vapour ducting piece flask (5.3) and attach it firmly by the flask device subsection (5.4).

*\* Note: Before installing the evaporator and receiver flasks, put Vaseline on the frosted glass to avoid break when it is to dismount.*

### Setting up

- ◆ By the height adjustment button (1.5), place the lifting mechanism up.
- ◆ The rotary evaporator has the possibility to put the glass set in an angle of 45 degrees. To do this, take the turned body and loose the bottom (1.4), place it manually at the proper angle and re-adjust.

*Attention!: Turn the body carefully to avoid knock and breakage of glass set.*

- ◆ Fill the thermostatic bath with water at least with 2/3 of its capacity.

*Note: Check regularly the level of the bath to avoid the water run out during its works.*

- ◆ Put down the lifting mechanism, with the button height adjustment (1.5), till the height you consider.
- ◆ Connect the water system as it is described in section “coil condenser Installation”.
- ◆ Connect the sample power supply circuit (optional) when you considerer necessary.
- ◆ If they are not using the equipment, keep the key closed to allow the application of the vacuum system.
- ◆ Connect the vacuum circuit (optional) when you considers necessary.
- ◆ T Programming.
  - pull the HEATING switch (3.2) [On display (SV (3.6) (programming T) (PV (3.7) real) T)].
  - By temperature adjusting button (▲▼) (3.5), program, or change at any time, the T of the bath.
  - Pilots (3.4) indicate the heating state (OUT) and overheating (ALM) from the water bath (when the real T exceeds in X ° C the programming T).

- ◆ Used the rotary speed adjustment (3.3) smoothly to get speed wanted for the evaporation flask (Screen speed display (3.8) (Model 9200 / 1)).
- ◆ When you finish work, push switch off button (3.1), disconnect all circuits and pick up the processed products.


- ◆ Pull the general switch (3.1) and, through the height adjustment button (1.5), place the lifting mechanism up.
- ◆ Loosen the screw connection in the main body, remove the rubber seal and place it perfectly fitted in the rim union (4.1) of the coil condenser.
- ◆ Then enter the rim of the coil condenser (4.1) in this connection and adjust it in the thread firmly.

**Note: Make sure it is securely attached and can not be removed!**

- ◆ Install the water circuit at the input and output of the coil condenser (4.4). We recommended to connect the water input at the top and output at the bottom.
- ◆ Place the input pipette (1.8) in the connection (4.2) for a continuous supply of the sample (Optional).
- ◆ Install the vacuum system at the corresponding output (4.5) (optional).
- ◆ Select the flask receptor you want to use, place it in the connection (4.3) \* and fit it with the clip.


Evaporation system installation. (Figure 5)

- 5.1 Central body thread
- 5.2 Turn lock hole
- 5.3 Vapour ducting piece
- 5.4 Flask device subjection

- ◆ Through the height adjustment button (1.5), place the lifting mechanism up.
- ◆ Loosen the screw connection in the main body (5.1)
- ◆ In one of the four turn lock holes (5.2) insert screwdriver supplied with the equipment, to block the movement.
- ◆ Insert the vapour ducting piece (5.3) in the connection and fix it, adjusting firmly in the screw.
- ◆ Remove the screwdriver.

- 3.1 ON/OFF: General Switch on/off
- 3.2 HEATING: heating bath Button
- 3.3 Rotary speed adjustment.
- 3.4 Works OUT (heating) / alarm ALM (overheating)
- 3.5 Temperature adjustment button (▼▲)
- 3.6 SV Screen (Programming T)
- 3.7 PV Screen (Real T)
- 3.8 speed Screen display (Model 9200/1)

**3. TECHNICAL SPECIFICATION**

<b>Code</b>	59200320	59200321	59300330
<b>Model</b>	9200	9200/1	9300
<b>Bath capacity</b>	3 L	3 L	5 L
<b>Cover Elevator</b>	130 mm	130 mm	150 mm
<b>Inclination Angle</b>	45°		
<b>Velocity range</b>	20-200 rpm		
<b>T Range</b>	Ambiente -100°C		
<b>T Resolution</b>	1°C		
<b>T Precision</b>	±1°C		
<b>Consumption</b>	1300W		
<b>Power supply</b>	220-240V 50-60Hz +/-10%		
<b>Measure</b>	580x450x790 mm		
<b>Weight</b>	~30 Kg	~30 Kg	~38 Kg

**Main technical features:**

- ◆ Motorized lifting system: by pressing a single button and without user's effort, the glass set can be lifted or taken down softly along a distance of 150m.
- ◆ Adjusting of the glass set angle: by manual adjustment to move the glass set in the most appropriate position and inclination.
- ◆ Thermostatic bath: with digital temperature control and Teflon cover to assure a better cleaning and higher duration.
- ◆ Dual digital display: for temperature setting and visualization of the current temperature of the bath. (Model 9200 / 1 also includes display of speed)
- ◆ Refrigerant in double coil: with higher condensing surface.
- ◆ Vacuum connection (optional / not included): to accelerate the process of evaporation by decreasing the T sample.
- ◆ Continuous sample supply (optional / not included): through pipette which allow the supply of the sample to the evaporation flask.

**4. INSTALLATION/SETTING UP****Preliminary inspection**

Unwrap the rotary evaporator, take off the involving plastic and take off the polispán protection in which it comes fitted.

Make sure that it does not present any damage because of the shipment. In case the equipment presents any damage tell it immediately to your transport agent or dealer so that they can make the claims in the correct time limit.

Please check that all the accessories are enclosed with the equipment:

- 1 Thermostatic bath
- 1 Coil Condenser 29/32
- 1 Vapour ducting piece
- 1 Pipetta flow sample (Plastic tuve included)
- Flasks 29/32:
  - 2 units of 1L
  - 1 unit of 2L
  - 2 unit of 3L (9300 model)
- 1 Clip 29/32
- 2 Schuko cable
- Fuses (2x5A / 1x10A)
- Screwdriver (Ø 4 mm)
- Warranty
- Instruction Manual

**Installation**

Before using this instrument, it is convenient for you to familiarize with its components and basic essentials.

**TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.**

- ◆ Please put the rotary evaporator on top of a horizontal, plane and stable table making a free space at least at 30 cm per side.
- ◆ Fit the bath to the main body, fitting it on the rail available for this purpose.
- ◆ Remove the protective material clean it carefully.

**Electrical installation assembly.**

- ◆ Insert the power current alternate cord (AC) to the base of 220V 50Hz  $\pm$  10% supplied with ground and the other end to the connection (2.1) of the equipment.

*Neither the manufacturer nor the distributor assume any responsibility for damage of the equipment, installations or injuries to people due to incorrect installations process. The voltage has to be 220 V, 50 Hz  $\pm$  10%.*

- ◆ Connect the Schuko cable supplied by one of its ends to the connection (2.2) of the central body and the other to the water bath (2-4).
- ◆ Finally, connect the sensor temperature cable that comes from the water bath in the connection (2.3) of the central body.
- ◆ When not to use the equipment for long periods of time, make sure it is disconnected from the network to prevent accidents.

**Coil condenser installation. (Figure 4)**

- |                               |                               |
|-------------------------------|-------------------------------|
| 4.1 Rim union central body    | 4.4 Input / Output coil water |
| 4.2 Pipette connection        | 4.5 Vacuum connection         |
| 4.3 Connection flask receiver |                               |