AUXILAB S.L.

Material de laboratorio
Laboratory supplies

DECLARACIÓN DE CONFORMIDAD CE
CE DECLARATION OF CONFORMITY
DÉCLARATION DE CONFORMITÉ CE

El fabricante | The manufacturer | Le fournisseur:

AUXILAB, S.L.

Declara que el equipo | Declare that the equipment | Declare que l'appareil:

VISCOSÍMETRO | VISCOSIMETER | VISCOSIMÈTRE

Código | Code | Code: 50801000.

Modelo | Model | Modèle: 801.

Cumple las siguientes directivas | Meet the following directives | Accomplit les directives suivantes:

73/23/CE | Directiva de seguridad eléctrica
Directive for electrical safety
Directive de sécurité électrique

89/336/CE | Directiva de Compatibilidad electromagnética (CEM)
Directive for electromagnetic compatibility y (EMC)
Directive de compatibilité électromagnétique (CEM)

Cumple las siguientes Normas: | Meet the following Standards | Accomplit les normes suivantes:

EN 61326 | Material eléctrico para medida control y uso en laboratorio. Requisitos de compatibilidad electromagnética (CEM).
Electrical equipment for measurement, control and laboratory use. EMC requirements.
Matériel électrique de mesure, de commande et de laboratoire. Prescriptions relatives à la CEM.

EN 61010-1 | Requisitos de seguridad de equipos eléctricos de medida, control y uso en laboratorio. Parte 1: Requisitos generales.
Safety requirements for electrical equipments for measurement, control and laboratory use. Part 1: General requirements.
Règles de sécurité pour appareils électriques de mesure, de régulation et de laboratoire. Partie 1: Prescriptions générales.

Fdo: Alfonso Ainciburu Sanz
DIRECTOR | GERENTE

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VISCOSÍMETRO ROTACIONAL ROTATIONAL VISCOSIMETER VISCOSIMÈTRE ROTATIONNEL

Modelo / Model / Modèle 801

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◆ Never try to repair the viscosimeter by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the viscosimeter (burns, hurts...) or damages in nearby equipments.

◆ In the event of breakdown please contact your distributor to overhaul through Nahita Technical Assistance Department.

Cleaning

◆ Never use scourers or substances that can grate for cleaning metallic parts such as stainless steel, aluminium, coatings, etc. as they damage the viscosimeter and produce an early ageing of the equipment.

◆ Use a fluff-free cloth dampened with soaped water that does not contain abrasives.

◆ For cleaning spindles, use neutral products with pH not higher than 6.

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ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

INSTRUCTIONS ON ENVIRONMENT PROTECTION

At the end of its life cycle, please, do not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

Thank you for choosing this equipment. We sincerely wish that you enjoy your Nahita viscosimeter 801. We highly recommend looking after this equipment according to what is stated in this manual.

Nahita develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

Please bear in mind the following:

◆ This manual is inseparable from the Nahita viscosimeter 801, so it should be available for all the users of this equipment.

◆ Follow all required applicable laws about safety and health (electrical security devices, inflammable, toxic and/or pathological products, etc.).

◆ You should carefully handle the equipment avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.

◆ Never dismantle the different pieces of the equipment to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the viscosimeter.

◆ To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.

◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting Nahita Technical Assistance Department by email to asistencia@auxilab.es.

◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).

◆ Overhaul is not covered by the equipment warranty.

◆ Operations made by non-qualified staff will automatically produce a loss of the warranty.

◆ Accessories (including their loss), are not covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.

◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Nahita Technical Assistance Department you should enclose the original invoice or a copy as guarantee.

◆ Manufacturer reserves the right to modify or improve the manual or equipment.

ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

INDEX OF LANGUAGES

Spanish 2-10
 English..... 11-18
 French 19-27

INDEX OF CONTENTS

1. USES OF THE INSTRUMENT..... 12
 2. DESCRIPTION..... 12
 3. TECHNICAL SPECIFICATIONS..... 13
 4. PRINCIPLE OF OPERATION..... 14
 5. INSTALLATION / SETTING UP..... 15
 6. MAINTENANCE AND CLEANING..... 17
 7. ACCESSORIES..... 18
 ANNEX I: CE CERTIFICATE..... 28

1. USES OF THE INSTRUMENT

Nahita viscosimeter is used to measure resistance of fluids due to their viscosity. The equipment presents a wide field of applications since it allows determination of viscosity of very different fluids such as fats, paints, plastics, foods, pharmaceutical products, etc...

2. DESCRIPTION

- 1.1. Scale indicator
- 1.2. Power on/off
- 1.3. Bubble level
- 1.4. Speed adjustment knob
- 1.5. Indicator braking lever
- 1.6. Spindle holder
- 1.7. Adjustable legs
- 1.8. Holding screw, depth adjustment
- 1.9. Height adjustment knob
- 1.10. Holding screw, height adjustment


CHART 2. Table of coefficients according to speed and spindle

Spindle	Speed (r.p.m.)			
	60	30	12	6
1	1	2	5	10
2	5	10	25	50
3	20	40	100	200
4	100	200	500	1000

Security

- ◆ The viscosimeter must be used by previously qualified staff that knows how the equipment works thanks to the user manual.
- ◆ Put the viscosimeter on top of a horizontal, plane and stable table making a free space at least at 30 cm per side.
- ◆ Do not put the viscosimeter near any warm supply (burners, blowlamps...), nor expose it directly to the sun, etc. Avoid vibrations, dust and dry environments.
- ◆ During operation, dangerous material such as flammable liquids or pathological material must be out of the working area.
- ◆ When you are not using the viscosimeter for a long period of time please make sure it is unplugged in order to avoid possible accidents.
- ◆ It is essential to have the equipment switched off and unplugged from the net before cleaning, checking components or replacing any piece (e.g. replacement of a fuse).
- ◆ Never try to repair the viscosimeter by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the equipment (burns, hurts...).
- ◆ Try not to spill any liquid on the control panel, though it is properly insulated. In case you have any doubts do please immediately unplug the equipment (see Maintenance).
- ◆ Made under the European regulations for electrical security, electromagnetic compatibility and security on machines.

6. MAINTENANCE AND CLEANING

To get the best results and a higher duration of the viscosimeter it is essential to follow the processes of use.

Note: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.

- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of this equipment.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ The viscosimeter is supplied with a wire. It has to be plugged to an earth connection and the socket should be handy and ready to unplug the equipment in case of emergency.

base of current 220 V, 50 Hz \pm 10%

Neither the manufacturer nor the distributor will assume any responsibility for the damages produced to the equipment during its installation or damages to people suffered by the improper use of the electric connection. The tension should be 220V 50Hz \pm 10%.

Setting up

◆ Remove the protecting cap that covers the thread of the spindle holder (1.6) and screw the selected spindle by turning it clockwise (to unscrew, turn the spindle in anticlockwise direction).

NOTE: Hold the spindle holder (1.6) with one hand and with the other carefully screw the spindle.


- ◆ Prepare the fluid to be measured in a beaker or a cylindrical container of at least 70 mm diameter and control its temperature precisely.
- ◆ Slowly, raise the recipient with the fluid until the spindle is submerged up to the circular mark.
- ◆ Release the indicator braking lever (1.5) and switch the motor on by pressing power on/off (1.2); turn the speed adjustment knob (1.4) until the desired speed is upwards matching with the red mark. After several turns, 20 or 30seconds, indicator reading tends to be stable. Readings can also be made at specific times.
- ◆ If the rotating speed is fast (30 or 60 r.p.m.) press the indicator braking lever (1.5) to cause the reading to settle.
- ◆ When the indicated value is too high or too low, change spindle or speed; it is better to perform measurements at a speed range from 30-90 rpm.

Selection of range, coefficient, spindle and speed

- ◆ First of all, the approximate range of viscosity of the fluid must be estimated. Then, select the appropriate speed and spindle according to table of ranges (Chart 1).
- ◆ For instance, if the estimated viscosity is 3000 cP, one of the following combinations must be selected: spindle n° 2 at a speed of 6 r.p.m or spindle n° 3 at a speed of 30 r.p.m.
- ◆ When viscosity is difficult to estimate, a high viscosity must be supposed. Try the spindle from larger to smaller and speed from lower to higher, until finding the most adequate combination for the viscosity of the fluid.

CHART 1. According to speed and spindle

Spindle	Speed (r.p.m.)			
	60	30	12	6
1	100	200	500	1000
2	500	1000	2500	5000
3	2000	4000	10000	20000
4	10000	20000	50000	100000

3. TECHNICAL SPECIFICATIONS

Measuring range	10-100000 mPa·s (cP)
Precision	1%
Spindles	N° 1, 2, 3 and 4
Rotor speed	6, 12, 30 and 60 r.p.m.
Tolerance	+/- 5% (on Newtonians fluids)
Power	220 V, 50-60 Hz


4. PRINCIPLE OF OPERATION

1. As shown in *figure 2*, the synchronous motor which is connected to the graduated disk, rotates at a steady speed making the spindle rotate through the spring and the indicator. If no fluid resistance is found, the spindle rotates synchronously with the spring, the indicator and the graduated disk, being the reading "0" in this case.

Otherwise, if the fluid offers resistance due to its viscosity, this will produce spring torsion; when the spring is balanced against the fluid resistance, the indicator connected to the spring will make a particular reading (the torsion moment of the spring) on the graduated disk.

By multiplying this reading by the specific coefficient for each spindle and speed (as shown in *chart 2*, the fluid viscosity in mPa·s or cP will be obtained.

1. Synchronous motor
2. Graduated disc
3. Indicator
4. Hairspring
5. Fluid to be measured
6. Rotor


2. Depending on each fluid requirements, the rotating speed can be selected, by the corresponding knob, between four different values (6, 12, 30 and 60 rpm).

3. A total of four spindles (n° 1, 2, 3 and 4) are supplied together with the viscosimeter. The appropriate spindle and speed must be selected according to the viscosity of the fluid to be measured (see enclosed charts).

4. Viscosimeter is equipped with a device to fix the indicator when the rotating speed is too fast (30 or 60 r.p.m.) and reading is not possible.

5. The spindle protector stabilizes measurement and protects the spindle. If the protector is used, measurements will be more stable.

6. Viscosimeter is supplied with a support stand, it is advisable to level the viscosimeter on the support stand and then, use a laboratory

jack to raise or lower the recipient with the sample to be measured.

5. INSTALLATION / SETTING UP

Preliminary inspection

◆ Unwrap the viscosimeter, take off all protections and without connecting the viscosimeter to the power supply, make sure that the equipment does not present any damage because of the shipment. In case it presents any damage, tell it immediately to your transport agent or dealer so that they can make the claims in the correct time limit.

◆ Please keep the original wrapping; you will always need it for returns enclosed with all the accessories supplied.

◆ Please check that all the accessories are enclosed with the equipment:

- Head
- Support stand
- Box with 4 spindles
- Accessory to work at high speed
- Instructions manual
- Warranty certificate

We will only accept any equipment return within 15 days after delivery and provided it comes in its original wrapping.

Installation

Before using the viscosimeter, it is convenient for you to familiarize with its components and basic essentials, as well as with its control functions.

PLEASE READ THOROUGHLY THE INSTRUCTIONS BEFORE CONNECTING AND OPERATING WITH THIS EQUIPMENT WITH THE AIM OF ACHIEVING THE HIGHEST FEATURES AND THE MAXIMUM DURABILITY OF THE INSTRUMENT.

◆ Please put the viscosimeter on top of a horizontal, plane and stable table making a free space at least at 30 cm per side. Do not put the viscosimeter near any warm supply (burners, blowlamps...), nor expose it directly to the sun, etc.

◆ During operation, dangerous material such as flammable or infectious liquids must be out of the working area.

Mounting the equipment

◆ Screw the rod on the support stand with the cogged surface at the front. Tighten the lower thread to avoid rod movements.

◆ Put the head on the guideline of the holding clamp, adjust depth according to necessities and fix the head on this position with the screw (1.8).

◆ Level the viscosimeter through the three adjustable legs (1.7) located at the base of the support stand, until the bubble (1.3) coincides with the inner circle of the level device.

◆ Plug the equipment to the net; insert the wire that feed the AC electric current to the