**CENTRÍFUGAS ANGULARES DIGITALES
ECO DIGITAL CENTRIFUGE
CENTRIFUGEUR DIGITAL ECO**

Nahita

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Thank you for choosing this equipment. We sincerely wish that you enjoy your Centrifuges Series 2640 Nahita. We highly recommend looking after this equipment according to what is stated in this manual.

Nahita develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT

Please bear in mind the following:

- ◆ This manual is inseparable from the Nahita Centrifuge Series 2640 ECO, so it should be available for all the users of this equipment.
- ◆ Follow all required applicable laws about safety and health (electrical security devices, inflammable, toxic and/or pathological products, etc.).
- ◆ You should carefully handle the equipment avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- ◆ Never dismantle the different pieces of the equipment to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the equipment.
- ◆ To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting Nahita Technical Assistance Department by email to asistencia@auxilab.es.
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the equipment warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the equipment warranty.
- ◆ Neither fuses (1A) nor accessories (including their loss), are covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Nahita Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.

9. MAINTENANCE AND SERVICE

- ◆ The normal service life of the rotor of this machine is three years, make sure to replace it with a new one beyond this period.
- ◆ It is not allowed to impact the rotor with sharp and hard object, during the use, make sure to prevent cracks on the rotor due to scratch or external injury.
- ◆ Be sure to regularly check if there are corrosion stains, grooves, fine cracks etc. on the rotor, if any of the above is present, please immediately stop using this rotor and contact the manufacturer.
- ◆ If finding the sample under centrifugation during the use splashes out and drips on the rotor, immediately absorb the liquid and clean the rotor locally.
- ◆ When cleaning the rotor, in order to prevent destruction of surface, please use sponge or cotton cloth soaked with cleaning agent to clean the rotor, then wash off the cleaning agent with distilled water or wipe the rotor with 70% alcohol, after cleaning, it is allowed to hang it upside down and let it dry.
- ◆ Cut off the power supply when the centrifuge is out of service.

ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

Instructions on environment protection

At the end of its life cycle, please, do not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

8. ANALYSIS OF COMMON FAULTS AND ELIMINATION THEREOF

<i>Fault symptom</i>	<i>Cause analysis</i>	<i>Elimination methods</i>	<i>Remarks</i>
No display or disorderly display	Blown out fuse (5A) of primary circuit	Change for a fuse of same size	
	Dropout of wire connected to display board	Open the cabinet and re-weld the loose wire. Open the cabinet and re-weld the loose wire.	
	Misoperation of chip computer and disorderly operating procedure	Turn off power switch of the host machine and start it several minutes later	Operate according to normal operating procedure
Unstable rotating speed	Control circuit or brushless controller has fault		Contact Nahita Technical Service
Big noise	Damaged drive motor		
	Test tubes are placed unbalancedly	Place test tubes symmetrically	
	Loosened nut on the rotor	Tighten the nut on the rotor	
Display present, but centrifugation cannot be carried out	Door cover is not closed well	Close door cover again	
	Speed fails to rise due to insufficient power supply voltage	Contact Nahita Technical Service	
	Failed to open the door and change the sample after last time centrifugation finishes	Open the door, then close it securely	
	Pushbutton switch fails due to poor contact		Contact Nahita Technical Service
No display after operation, but the set speed value is displayed after a period of time	The program stops the machine for protection after failure to read the speed signal.		Contact Nahita Technical Service

ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

INDEX OF LANGUAGES

Spanish 2-11
 English 12-21
 French 22-31

INDEX OF CONTAINS

FOREWORD 14
 IMPORTANT HINT 14
 SAFETY PRECAUTIONS
 PRECAUTIONS FOR OPERATION 14
 PRECAUTIONS FOR SERVICE ENVIRONMENT 15
 TRANSPORT AND STORAGE CONDITIONS 15
 MAIN TECHNICAL PARAMETERS AND PERFORMANCE INDEXES
 MAIN TECHNICAL PARAMETERS 15
 PERFORMANCE INDEXES 16
 OPERATING PRINCIPLE 16
 STRUCTURAL DIAGRAM AND SYSTEM BLOCK DIAGRAM OF CENTRIFUGE 17
 USE AND OPERATION
 PREPARATION AND CHECKS PRIOR TO USE 18
 DESCRIPTION OF OPERATION OF THE APPARATUS CONTROL SYSTEM 19
 ANALYSIS OF COMMON FAULTS AND ELIMINATION THEREOF 20
 MAINTENANCE AND SERVICE 21

1. FOREWORD

Main features and application scope:

The Series ECO 2640 low-speed centrifuge (hereinafter referred to as centrifuge) adopts advanced technology such as its drive unit and microcomputer control, etc., which is featured with small volume, short startup time, stable running, low noise, convenient operation, high service efficiency, safety and reliability, etc., it can be widely used for separation of solution constituents in such fields as of biology, chemistry, medical science, etc

2. IMPORTANT HINT

- ◆ When the apparatus is working, the mass deviation of each centrifugal tube which contains centrifugation sample shall not exceed 3g, and the tubes must be placed in the rotor symmetrically.
- ◆ When the apparatus operates, make sure to ensure that its operating voltage falls within the specified range. When the apparatus cannot operate normally, first check if the power supply voltage is normal or not!
- ◆ Prior to startup of the apparatus each time, be sure to carefully confirm if the holddown nut of the rotor is fastened!
- ◆ When using the apparatus, be sure to use a single-phase 3-hole power socket and ensure that it is grounded reliably!
- ◆ In order to ensure the safety of the apparatus and the operator, when the apparatus runs, the set value of its rotational speed shall not exceed its highest speed!
- ◆ After power off of the apparatus, please wait for enough time (over 3 min), then switch on the power supply!

3. SAFETY PRECAUTIONS

Precautions for operation

1.1. Prior to power-on

- ◆ Each time when using this machine, be sure to carefully check if there are cracks on the rotor and centrifugal pipe used for this machine or if there is serious corrosion phenomenon, if any, change it immediately.
- ◆ Keep interior of centrifuge chamber clean and prevent water accumulation and prevent particular foreign matters from entering into the chamber.
- ◆ The power supply to be used must be identical with the input voltage specification of the centrifuge, that is, single-phase 220V/50Hz and ensure that the input end of the power supply has protective ground wire.

1.2 During the use

- ◆ During the apparatus acceleration or deceleration, presence of short-time vibration is a normal phenomenon, it is not necessary to turn off the power switch of the main machi-

Operation description

Turn on the power switch, the centrifuge will respectively display the set time and rotating speed reserved after the last operation. At this time, if you need to modify such parameters as "timing", "speed", "centrifugal force" etc., you can operate according to the following methods:

Modify time

- ◆ Timing setting can be divided into: continuous running and running to set minute.
- ◆ By pressing the "**min**" (2) key, the switchover will be made cyclically between two statuses, for continuous running status, the timing window displays $\square\square\square$; for minute status, the timing window displays $\square\square$ (1-99 min).
- ◆ By pressing the "timing setting" key for 1-2 times, the timing window will switch over to the corresponding status, at this time, the window will flicker, press "**▲**" "**▼**" to modify it to the required time value.

Modify rotating speed

- ◆ After pressing the "**speed**" (3) key once, use "**▲**" or "**▼**" to modify the rotating speed.

Modify centrifugal force

- ◆ Before startup, after pressing the "**RCF**" (4) key once, use "**▲**" or "**▼**" to modify the centrifugal force; after startup, by pressing the "centrifugal force" key once, it will display the centrifugal force corresponding to the current rotating speed. .

Note:

- *the speed and centrifugal force are set alternately, i.e. once one of them is set, the other is changed accordingly.*
- *The current modified parameters such as timing, speed or centrifugal force can be confirmed as effective values only after normal running, and stored at the same time.*

The centrifuge provides jog function, hold the "**Short spin**" (7) key, the centrifuge begins to operate at the set speed, if the "**Short spin**" (7) key is released midway, then the centrifuge begins to reduce its speed until it stops operation.

Fault alarm prompts:

<i>Alarm sound</i>	<i>Meaning of alarm sound</i>
Sound: tit-tit-tit...	Indicates the door is unexpectedly opened during the running.
Sound: tit, tit, tit, tit, tit	Indicates the centrifugation has finished.
Sound: tit-tit-tit-tit...	Buzzes once every minute to prompt that the sample has not been taken out after centrifugation.


7. USE AND OPERATION

Preparation and checks prior to use

Under the condition of apparatus power-off, gently press the door button on the right side of the apparatus with right hand, meanwhile gently pull the centrifuge door cover with left hand so that the door cover is opened upward.

Turn the rotor with hands, the rotation should be light and flexible. Tighten the nut on the rotor with a wrench, don't put it into the balanced test tube until there is no abnormal phenomenon.

Attention: When placing the test tubes, be sure to place them into the rotor symmetrically, if the number of centrifugal tubes is not even, use empty test tubes and fill water into them so as to make the number of the tubes to be even. Otherwise, machine vibration will take place during the centrifugation, which may damage the machine!

1 Display: Displays relevant setting status, when a status is enabled, the corresponding indication lamp turns on.

2 Min.: Used for setting the time required for sample centrifugation, the unit is minute or continuous running.

3 Speed: Used for setting the rotating speed required in the process of sample centrifugation

4 RCF: Displays the centrifugal force at current speed, and can be used for setting the centrifugal force required in the process of sample centrifugation

5 Tecla ▲: Used for increasing the set value

6 Tecla ▼: Used for decreasing the set number

7 Short spin: Used for jog running of the centrifuge

8 Start: Used to make centrifuge begin to run

9 Tecla Stop: Used to make centrifuge stop running

ne or "stop" button on the operating panel. In the event there happens the phenomenon like power failure midway or forced cutoff of power switch of the main machine etc., don't open the door immediately, it must be opened only after the motor has stopped rotation (for about 5~10 min).

◆ Each time after the centrifuging process finishes, the centrifuge will wait for the user to open the door and take the sample, if the centrifuge continues centrifuging without the door opening, it will not respond.

1.3 After finishing use, take care to protect the apparatus, particularly the rotor and centrifugal pipe so as to prevent it from corrosion due to contamination by acid or alkali liquid.

1.4 The once continuous running of the centrifuge should advisably not exceed 60 min.

Working conditions.

In order to ensure stable and reliable operation and normal functioning of the centrifuge, make sure to ensure it meets the following conditions:

◆ Ambient temperature: 5°C~40°C

◆ Relative humidity: ≤80%

◆ Supply voltage: 220V±22V, 50Hz±1Hz

◆ The service environment should be well ventilated, and no foreign matters like dust, flocks, metal chips etc. enter into the centrifuge.

◆ Resistant to corrosive gas and strong electromagnetic interference.

◆ When using it, place it on a firm and level operating platform which is resistant to vibration.

Transport and storage conditions

◆ According to the requirements of ordering contract, the centrifuge should be transported under the packed condition, during the transportation, the centrifuge should be prevented from fierce impact, rain and strong sunshine.

◆ Completely packed centrifuge should be stored in well-ventilated room with relative humidity not higher than 80% and without corrosive gas.

4. MAIN TECHNICAL PARAMETERS AND PERFORMANCE INDEXES

2.1 Main technical parameters

Model	2640/8	2640/12
Max. speed (rpm)	4000	4000
Max. FCR (×g)	2432	2432
Capacity	8 ×15ml	12 ×15ml

2.2 Performance indexes.

Noise: ≤ 60 dB
 Timer range: 1-99 min. / continuous / short-time centrifugation
 Temperature rise (°C) of chamber for 20min centrifugation: ≤10 °C
 Power supply voltage: AC220V±22V 50Hz±1Hz
 Fuse: Φ5×20 (5A)
 Power: 80 W
 External dimensions L×W×H: 380×305×290 mm
 Weight: 20 Kg

5. OPERATING PRINCIPLE

Upon centrifugal operation, place the centrifugal vessel (centrifugal bottle, centrifugal test tube etc.) containing equal-part test liquid symmetrically in the centrifugal holes around the rotor and separate the test liquid under the centrifugal force generated through high speed rotation of the rotor driven by the electric motor. The magnitude of its relative centrifugal force (RCF) depends on the rotary diameter r and speed n during the centrifugation of the sample, its calculating formula is as follows:

$$R.C.F.= 1.118 \cdot 10^{-5} \cdot r \cdot n^2 (\times g)$$

Where $1.118 \cdot 10^{-5}$ una constante.
 r= rotary diameter (cm)
 n= rotating speed (r/min)

The time T required for separation and settlement of particles in the mixed liquid is calculated by the formula:


$$T = \frac{27.4 \times (\ln R_{max} - \ln R_{min}) \mu}{n^2 r^2 (\sigma - \rho)} \text{ (min)}$$

Where:
 ρ – density of mixed liquid (g/cm³)
 μ – viscosity of mixed liquid (P)
 n – rotating speed (rpm)
 r – diameter of particles (cm)
 σ – density of particles (g/cm³)

Rmax – horizontal distance from bottom surface of centrifugal test liquid to shaft center (cm)

Rmin – horizontal distance from liquid surface of centrifugal test liquid to shaft center (cm)

Schematic diagram for conditions of test tube during centrifugation

6. STRUCTURAL DIAGRAM AND SYSTEM BLOCK DIAGRAM

Schematic structural diagram of centrifuge

The centrifuge mainly consists of the centrifuge body assembly, rotary assembly, shock damping system, control system etc., its structural diagram is shown as below:

1. Door cover assembly
2. Rotor system
3. Casing assembly
4. Motor
5. Shock damping system
6. Cushion block
7. Hinge