ANEX I: CE CERTIFICATE

AUXILAB S.L.

CE DECLARATION OF CONFORMITY
ZUZI MICROSCOPE SERIES 128 of AUXILAB, S.L. for the Directive of Machines
(89/392/CEE modified) and the regulations adopted for their transposition

NAME OF THE MANUFACTURER / IMPORTER:
AUXILAB, S.L.

ADDRESS: **Polígono Morea Norte, 8**
31191 Beriáin (Navarra)

WE STATE THAT:

ZUZI MICROSCOPE SERIES 128
CODES 50128000, 50128001, 50128002, 50128022

Are designed and manufactured according to:

- ◆ Directive 89/392/CEE, including the modifications and the national regulations that transpose them.
- ◆ Directive 73/23/CEE modified over the electric security.
- ◆ Directive 89/336/CEE modified over the electromagnetic compatibility.

And that the following harmonized rules have been applied (or part of them):

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BERIAIN 4th DECEMBER 2007

Signed by: **ALFONSO AINCIBURU SANZ**
DIRECTOR/MANAGER

Polígono Morea Norte, 8 31191 Beriain (Navarra) - Spain. Tel. 948 310 513 Fax 948 312 071
Internet: www.auxilab.es · Email: correo@auxilab.es

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This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.

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◆ Check the prepared slide is not the other way round and check also the slide's thickness, which has to have a standard thickness of 0.17 mm.

If the slide slips under the mechanical stage:

◆ Adjust the distance between the mechanical section and the stage with the screws located on the rear part of the mechanical stage (the screws have a standard size), making sure there is enough space to avoid the slide passing through and obtain a smooth movement without scratching the stage. Once it is properly adjusted notice you should not hold the microscope by the mechanical stage, as its weight could loosen the equipment. In case the problem persists send the equipment to Zuzi Technical Assistance Department through your distributor.

8. TROUBLESHOOTING

◆ As the head is rotary you should rather observe by the stage's frontal part, as it facilitates access to the mechanical control knobs.

◆ Use of the immersion objective:

- To use all the numerical aperture of the immersion objective (100x), add a drop of immersion oil in the space between the objective and the slide. Proceed as follows:

- Focus the sample with a lower magnifying objective (40x) and put it in the centre of the field of view. This is essential because it is very difficult to find structures with the immersion objective.

- Turn slightly the nosepiece so as no objective stays over the slide.

- Add a drop of immersion oil over the region of the sample to be observed.

- Turn completely the nosepiece and put the immersion objective over the slide (both the sample and the objective lens must be in contact with the immersion oil). Adjust focusing with the fine adjustment knob.

- Avoid the appearance of bubbles in the oil film, since it spoils the resulting image.

- After each observation, clean the objective lens with a special lens tissue to remove the rests of immersion oil. Never let the oil dry on the lens.

NOTE: never submerge the objective in any liquid nor use solvents such as xylol or toluene for cleaning, since it may provoke the detachment of the lens.

◆ Handling the stage: First, move the sample in consecutive steps, from the front to the back; each fraction is observed to detail by means of a lateral movement. Thus, the sample can be moved in a uniform way and so explore all its area.

INSTRUCTIONS ON ENVIRONMENTAL PROTECTION

Do not dispose of this product in the usual household garbage at the end of its life cycle; hand it over at a collection point for the recycling of electrical and electronic appliances. The symbol on the product, the instructions for use or the packing will inform about the methods for disposal.

The materials are recyclable as mentioned in its marking. By recycling, material recycling or other forms or re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

there is not any diffuser filter.

- ◆ Eyepieces are well fit and objectives well screwed.
- ◆ Eyepieces and objectives are perfectly cleaned and free of particles. To check cleaning of the whole optical system proceed as follows:
 - a) Make the eyepieces turn checking the little specks are moving; if so, clean them. Make the head turn. You should never disassemble the head, but you can clean it with delicacy by blowing the accessible surfaces of prisms with a plastic bulb.
 - b) Turn the objective and, in case the parasite images also turn, clean it with the help of a dry brush to remove the dust. Watch the front part of the microscope with a magnifying glass or an inverted eyepiece.
- ◆ There is enough immersion oil and it does not contain bubbles or impurities.
- ◆ Thickness of slides, cover glasses and means of assembly is not too much as to avoid focusing with medium or high magnification (there are standard sizes for both the slides and cover glasses). Slides and cover glasses should be clean. Check that the slide is properly located and there are not two cover glasses superimposed.
 - The combination objective-eyepiece is appropriate; the eyepiece can be too powerful for the selected objective.
 - The prepared slide is well done; check it with a "test" prepared slide.

8. TROUBLESHOOTING

Before sending the microscope back to Zuzi Technical Assistance Department, check the following:

If the lamp does not switch on:

Check it is properly connected to the right current intake.

Check both fuse and lamp are in good conditions. If not, replace the fuse or the lamp.

If the field of vision is cut:

◆ Check the nosepiece is properly fit. To do that, turn it slightly back and forth until it fits. Check it is well centred and if not, send it to Zuzi Technical Assistance Department through your distributor.

If there is dust or dirt in the field of view:

◆ Check there is dust on the precondenser lens, upper lens of the condenser or the eyepiece, as well as the cleanliness of the sample. Once you find out where the dirt is, clean it as explained before.

If an area in the field of view is out of focus:

◆ Check the objective is properly located on the luminous beam's path. If not, turn it until it is in an appropriate position.

◆ Check the sample is properly located on the stage, making sure it is held by the stage clips. Then, if it does not focus revise the coarse adjustment knob.

If the frontal lens of the objective touches the prepared slide when focusing:

◆ Check the prepared slide is not the other way round (being the slide over the cover glass) and placing it properly in case it is appropriate to do so.

If the frontal lens of the more powerful objective touches the prepared slide when passing from a low-magnification objective to a superior one:

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
Dependiendo del modelo:

voke malfunction.

- ◆ Please use the plastic cover provided to keep the microscope away from dust lying on the optical parts when it is not used for a long period of time.
- ◆ Please keep the original packaging either to transport the equipment, when it is not being used for a long time or when you send it for an overhaul.

Cleaning

NOTE: Zuzi recommends using the microscope cleaner kit (code 89000001).

- ◆ Never use scourers or substances that can grate for cleaning metallic parts such as stainless steel, aluminium, coatings, etc. as they damage the microscope and produce an early ageing of the equipment.
- ◆ Use a fluff-free cloth dampened with soaped water that does not contain abrasives.
- ◆ Lenses must not be disassembled by the user. Were there any dust or dirt to be cleaned, clean it with a natural horse hair brush or a smooth piece of cloth, fluff-free, dampened with a bit of xilol or toluene.
- ◆ To remove dust from lenses, blow them with a plastic bulb or clean them with a natural horsehair brush.
- ◆ When cleaning mechanical or plastic parts do not use abrasive solvents that could damage either paint or finish. Always use neutral detergents.

ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DEFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

6. CHOOSING OBJECTIVES AND EYEPIECES

Zuzi microscopes series 128 are supplied with WF10x eyepieces and with 4x, 10x, 40x and 100x achromatic objectives, but can be used with other eyepieces and objectives. To know our offer of eyepieces and objectives look up our catalogue or visit our web page (www.auxilab.com).

The observed image loses surface and sharpness as magnification increases. The mentioned increase should be done by changing objectives and putting each time a more powerful one and not by changing eyepieces to a higher magnification, as eyepieces only magnify the image obtained with the objective. Thus, the more magnification the eyepiece has, the higher will be the loss of sharpness, clarity and surface of the resulting image.

For routine observations use the eyepiece with lesser magnification with more powerful objectives. The eyepiece with the higher magnification should be kept back for particular occasions, bearing in mind that it decreases definition and does not increase resolution.

7. CAUSES OF A DEFECTIVE IMAGE

In case of a defective image check the following:

- ◆ Illumination is well done and luminous intensity is neither excessive nor too weak. Never adjust it with the condenser's diaphragm. Both the condenser and the lamp should be well centred. Check that between the field diaphragm and the aperture diaphragm

Security

- ◆ The microscope must be used by previously qualified staff that knows how the equipment works thanks to the user manual.
- ◆ Put the microscope on top of a horizontal, plane, stable table, having a safety area of at least 30 cm per side.
- ◆ Do not place the microscope near any warm supply (burners, blowlamps, etc), nor expose it directly to the sun. Avoid vibrations, dust and dry environments.
- ◆ During operation dangerous materials such as flammable or pathological substances must be out of the safety area.
- ◆ When you are not using the microscope for a long period of time please make sure it is unplugged in order to avoid possible accidents.
- ◆ It is essential to have the equipment switched off and unplugged from the net before cleaning, checking components or replacing any piece (e.g. replacement of a fuse).
- ◆ Never try to repair the microscope by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the equipment (burns, hurts...).
- ◆ Made under the European regulations for electrical security, electromagnetic compatibility and security on machines.

5. MAINTENANCE AND CLEANING

To get the best results and a higher duration of this microscope it is essential to follow the processes of use.

Note: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.

- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of this equipment.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ The microscope is supplied with a Schuko standard wire. It has to be plugged to an earth connection and the socket should be handy and ready to unplug the equipment in case of emergency.
- ◆ Never try to repair the microscope by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the equipment (burns, hurts...) or damages in nearby equipments.
- ◆ In the event of breakdown please contact your distributor to overhaul through Zuzi Technical Assistance Department.
- ◆ In case the lamp blows up, replace it for another one, making sure it is a 6 V 20 W Zuzi original one. Please take care and do not touch the lamp with bare hands. Do not use lamps with a higher power, as they could provoke over heating or any malfunction.

IMPORTANT: Before replacing the lamp make sure the microscope is disconnected from the electric intake. Check you are using the proper lamps, as other types can pro-

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◆ The stage must be properly centred; otherwise, when turning the stage the sample's object would be out of the field of view. To centre a prepared slide it is necessary to start with few magnifications and engrave, either with ink or a diamond, a cross or a circle with 1 mm diameter. Then lace the sample in the middle of the stage, fix it with the holding clips and then turn the stage.

◆ It is inevitable that the object gets out of the field of vision. As it must remain centred, you should move it with your hand. After sizing up it previously, you have to turn the stage a complete turn without the object being out of the field of vision. In this moment you have to turn the nosepiece to a higher magnification making sure you follow the same steps; doing it as explained you will keep the sample centred with any optical arrangement.

◆ Once the stage is centred and the turning sample is placed between the crossed polarizers, two cases can come up while observing:


- If the sample disappears 4 times and appears another 4 times in positions that separate 45° one from each other (which can be done and verified with the stage's graduated circle), we can state we have a birefringent material. Do not forget that with overly thin materials or slightly anisotropic ones, lighting must be extremely weak.

- If the field of vision remains dark during all the time we are turning the sample, we can consider three different possibilities:

- The structure is birefringent, but it is badly positioned. If its axis coincides with the microscope's, it would be acting as if an isotropic object. In this case, relocate the object or replace the prepared slide.

- The material is slightly birefringent, not enough to produce a luminous signal.

- Finally, we can consider we have an isotropic object.


3. Eyepiece video cameras

| Code | Resolution | Output |
|----------|------------------|-----------------|
| 59140052 | 640x480 pixels | Composite video |
| 59140060 | 640x480 pixels | USB 2.0 |
| 59150060 | 1280x1024 pixels | USB 2.0 |

4. Double head for demonstration (code 90119910)

This accessory allows two users to observe the sample at the same time.

5. Eyepieces and objectives

| <u>Eyepieces</u> | | <u>Achromatic objectives</u> | |
|------------------|---------------------------|------------------------------|---------------|
| Code | Description | Code | Magnification |
| 90100104 | Widefield eyepiece WF/5x | 90100203 | 20x |
| 90100126 | Widefield eyepiece WF/16x | 90100206 | 60x(R) |

Birefringence is a characteristic of certain substances that partially restore light when placed between two crossed polarizers, that is, they have double refraction (double refraction is known as the monochromatic light beam that makes differ two refracted beams when simultaneously entering through certain materials).

Uses of polarized light microscopes cover different fields such as the study of vitamins, starch or biological anisotropies; but studies of worn stones (crystal structures) until they constitute thin preparations are one of the most appealing applications. In these preparations constituent minerals are easily identified by its crystals' shape, as well as thanks to the bright colours acquired when seen through crossed polarizers.

There are some precautions to be taken when using this type of microscopes:

- ◆ As a consequence of absorption phenomena produced in both polarizers a loss of light is caused, so it is convenient to avoid any odd light when using polarized light. Thus, it is recommended to use a large sheet of black paperboard (large enough to hide the observer's head) as a protective screen.
- ◆ In the case of tasks that require quite an accurate and meticulous work with slightly birefringent structures we should make sure neither slides nor cover glasses present traces of birefringence; for this reason, we will previously observe slides and cover glasses with the stage and the analyzer.
- ◆ For the study of dried up crystalline precipitates it is important to examine it in two parts: a part being dry and the other prepared with Canadian balsam, as it stands out better the interference colours. In this case it is advisable to have a little amount of balsam for some peripheral crystals to be in the air.



To put the petrographic stage, proceed as follows:

- ◆ Remove the clips (model 128) or the attachable mechanical stage (models 128/1, 128/2 and 128/22) supplied with the fixed stage (*Figure 5*).
- ◆ Put the petrographic stage on the fixed stage matching the holes of the petrographic stage with the holes marked as "c" in the fixed stage. Hold firmly the petrographic stage with the screws supplied (*Figure 7*).
- ◆ Put the analyzer supplied on the eyepiece (*Figure 8*). The microscope is ready to be used

To make a birefringence study proceed as follows:

- ◆ Make the "0" of the rotary platform coincide with the mark on the stage so as to read the rotating angle.
- ◆ Put the analyzer on the eyepiece and, without placing any sample on the stage, turn the analyzer until observing a completely dark image (crossed polarizers). Please be careful and make sure you do not move the analyzer while observing, as the reading done would be entirely wrong.
- ◆ Place the sample on the stage.
- ◆ Focus the object.

Figure 7

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
ases the focusing depth. Diffraction limit must not be surpassed as an excuse to increase the focusing depth. Diaphragm must not be used to reduce luminous diffraction. In spite of going straight on, light simple harmonic motion disperses when passing through a very closed diaphragm, provoking fuzzy vision and limiting the magnifying capacity of the microscope.

Accessories

1. Attachable mechanical stage (code 90118590)

A mechanical stage can be attached on the stage of the Zuzi microscope model 128 to move the sample through the X and Y axes. The attachment present two lateral wheels and graduated axes to easily find structures of interest.

To put the attachable mechanical stage proceed as follows:

- ◆ By using a screwdriver, remove the clips supplied with the fixed stage (Figure 5).
- ◆ The attachable stage should fit in the holes of the stage marked as "b" (Figure 6). Put the attachable stage so as the screw of the attachment matches the small central hole at the upper part of the stage. The small metallic pieces at both sides of the screw must fit into the holes at both sides of the central hole of the stage (Figure 6).
- ◆ Once the attachment is perfectly and securely placed on the stage, put the slide and hold it with the clamp of the attachment. By turning the lateral wheels of the attachment, the slide is moved over the stage.



Figure 5

Figure 6


2. Petrographic stage (code 90118901)

Zuzi microscopes series 128 can be used with a petrographic stage to study the optical properties of different structures according to birefringence.

eyepieces until there is a total fusion of the two images.

- Dioptre compensation: a more or less marked difference in the focusing is appreciated by closing alternately one eye, then the other; the dioptre compensation corrects this effect. The microscope is provided with a dioptre compensation ring at the left eyepiece holder. First of all, adjust the dioptre ring to zero and looking through the right eyepiece, focus the sample specimen using the coarse and fine adjustment knobs. Once the sample is focused, look through the left eyepiece and adjust the dioptre compensation ring until obtaining a clear, sharp image.

Once the values of interpupillary distance and dioptre compensation are found it will be very useful to learn them by heart, mostly if the microscope is shared with other users. This is highly advisable to avoid repeating the adjustment each time the microscope is used.

- ◆ Lower the stage according to the appropriate working distance for the objective with lesser magnification before placing the sample on it.
- ◆ Put the slide on the stage and, by looking directly at it (not by the eyepiece), place the sample centred on the stage opening.
- ◆ Looking through the eyepiece using the objective with lesser magnification, raise the stage by using the coarse adjustment knob until the image of the sample appears. Then, move the fine adjustment knob until the image is focused.
- ◆ Once the sample is focused with the objective with lesser magnification, screw the focus limit screw (*Figure 4*), so as avoiding the sample to hit accidentally the objectives and producing damages on them and on the sample.
- ◆ **IMPORTANT:** operation should always start by using the objective with lesser magnification, as it facilitates focusing and makes impossible to ruin the slides or get the objectives dirty. Furthermore, it is essential for adjusting the lighting previous to using the higher magnification, and it helps to see the whole topographic structure of the sample so as later concentrate on the more appealing parts to observe them at higher magnifications.
- ◆ By turning the nosepiece, pass to an objective with higher magnification. As they are parfocal objectives it will only be necessary to adjust focusing with the fine adjustment knob.
- ◆ Adjust both the diaphragm aperture until getting the light intensity desired. Normally, when operating with an intense lighting the condenser's numerical aperture should be lesser than the objective's, in order to avoid glares on the field of view that may distort the image contrast.
- ◆ Iris diaphragm enables decreasing the numerical aperture of the condenser to a value similar to the objective's, thus avoiding marginal illumination that decreases contrast. Therefore, diaphragm must be more closed when using a lesser magnification objective with a small numerical aperture, and gradually open the diaphragm when user higher magnification objectives.
- ◆ The closing of the diaphragm by reducing the optical system's numerical aperture incre-

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◆ Please put the microscope on top of a horizontal, plane and stable table making a free space at least at 30 cm per side. Do not put the microscope near any warm supply (burners, blowlamps...), nor expose it directly to the sun, etc.

◆ **VERY IMPORTANT:** Never hold the microscope by the stage nor the tube, as all the equipment's weight would fall on the fine adjustment screw and mechanical and precision parts would be slowly eroded.

◆ Avoid inflammable or toxic substances in the working area.

◆ The operator has to assume a comfortable, upright position with the back straight. It is convenient to work in a dark table in order to avoid an awkward fatigue and parasite light, as it may dazzle the operator or decrease the image sharpness.

◆ Please use the plastic cover provided whenever the microscope is not used so as to avoid dust laying on the optics.

◆ If you are not using the microscope for a long period of time please make sure it is disconnected from the net and protected from dust (this way you will avoid accidents and will extend its working-life).

Setting up

◆ Lower the stage with the help of the coarse adjustment knob. Screw the objectives on the nosepiece so as when it is turned clockwise the objectives follow an ascendant way (4x-10x-40x) (**Figure 2**).

◆ The microscope is supplied with the eyepiece/s already placed. If changing the eyepieces is necessary, loosen the screw located at the ocular tube, change the eyepieces and fixed them again with the screw (**Figure 3**). The microscope head is normally placed facing the front of the microscope, but if it is necessary it can be placed in any other direction since it can be rotated in an angle of 360°.

◆ Zuzi microscopes series 128 are supplied with a Schuko standard wire. Please insert the wire that feeds the AC electric current in the base of current 220V 50Hz $\pm 10\%$ provided with earth wire. The tension should be 220 V, 50 Hz $\pm 10\%$.

Neither the manufacturer nor the distributor will assume any responsibility for the damages produced to the equipment during its installation or damages to persons suffered by the improper use of the electric connection. The tension should be 220V 50Hz $\pm 10\%$.

◆ Models provided with binocular observation system must have the eyepiece tubes properly adjusted. This will depend on:

- Interpupillary distance: Adjust the interpupillary distance by separating or joining the

Figure 2



Figure 3

- **Total magnification:** It is the result of multiplying the eyepiece's magnification by the objective's magnification.

| Eyepieces | Objectives | 4x | 10x | 20x* | 40x(R) | 60x(R)* | 100X(R)(I) |
|-----------|---------------------|------|------|------|--------|---------|------------|
| WF/10x | Total magnification | 40 | 100 | 200 | 400 | 600 | 1000 |
| | Field of view (mm) | 4.50 | 1.80 | 0.9 | 0.45 | 0.30 | 0.18 |
| WF/16x* | Total magnification | 64 | 160 | 320 | 640 | 960 | 1600 |
| | Field of view (mm) | 2.75 | 1.10 | 0.55 | 0.27 | 0.18 | 0.11 |

(R) Retractable (I) Immersion * Optional

- **Numerical aperture:** It determines the objective's properties. The bigger it is, the brighter and better resolved the image will be. The numerical aperture of the objective must always be lower than or the same as the numerical aperture of the condenser.

- **Working distance:** it is the distance, in mm, between the microscope slide and the objective's front lens when the image is focused.

- **Focal distance:** It is the distance, in mm, from the main image plane to its image focus.

- **Resolution:** It is the reciprocal value of the separating power, that is to say, the smallest distance between two points on a sample that can still be distinguished as two separate entities.

- **Field of view:** Dimensions of the real field that is being observed, expressed in mm.

4. INSTALLATION AND SETTING

Preliminary inspection

Unwrap the microscope, take off the involving plastic and take off the polyspan protection in which it comes fitted. Take off all the protective items and, without connecting the equipment to the net, make sure that it does not present any damage because of the shipment. In case the microscope presents any damage tell it immediately to your transport agent or dealer so that they can make the claims in the correct time limit.

Please keep the original wrapping; you will always need it for returns enclosed with all the accessories supplied.

Please check that all the accessories are enclosed with the equipment:

- | | |
|--|------------------------|
| - Microscope with a Schuko standard wire | - Fuse |
| - Achromatic objectives (depending on the model) | - Cleaner kit |
| - Eyepiece/s WF/10x | - Plastic cover |
| - Green filter | - User's manual |
| - Halogen lamp (6 V, 20 W) | - Warranty certificate |

We will only accept any equipment return within 15 days after delivery and provided it comes in its original wrapping.

Installation

◆ Before using this instrument, it is convenient for you to familiarize with its components and basic essentials.

PLEASE READ THOROUGHLY THE INSTRUCTIONS BEFORE CONNECTING AND OPERATING WITH THIS EQUIPMENT.

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(green filter of 32 mm Ø included)

- Abbe condenser (N.A. 1.25) height adjustable, with iris diaphragm and movable filter holder (green filter of 32 mm Ø included)

Iris diaphragm

Located under the stage and condenser, it controls the amount of light that enters the condenser.

Illumination

Low tension illuminator with halogen lamp 6 V 20 W, built-in feeder (220 V, 50 Hz ± 10%), switch and intensity adjustment potentiometer

Packaging


- Plastic protective cover
- Expanded polystyrene mould
- Cardboard box with protective items to send the equipment

Optical specifications

- **Achromatic objectives:** when talking about objectives, it is possible to distinguish between two types. Dry objectives are the ones where there is only air between the objective and the microscope slide, whereas immersion objectives are those that need immersion oil to be placed between the lens and the slide so as to obtain a higher luminosity.

Zuzi microscopes series 128 are provided with three dry objectives (4x, 10x and 40x) and one immersion objective (100x, only models 128/2 and 128/22). The 40x and 100x objectives are also retractile, which means that they present a built-in mechanism to avoid the breakage of the objective or the slide if the objective contacts the sample.

The features of objectives are inscribed on the barrel of each objective in the following way:

- 40:** Magnification
- 0.65:** Numerical aperture (N.A.)
- 160:** Mechanical tube length
- 0.17:** Cover glass thickness

| | 4x | 10x | 20x* | 40x (R) | 60x (R)* | 100x(R)(I) |
|------------------------------|-------|------|------|---------|----------|------------|
| Magnification | | | | | | |
| Working distance (mm) | 37.50 | 7.32 | 1.26 | 0.63 | 0.41 | 0.19 |
| Numerical aperture | 0.10 | 0.25 | 0.40 | 0.65 | 0.85 | 1.25 |

(R) Retractable (I) Immersion * Optional

- **Eyepieces:** They are composed by two lenses separated by a diaphragm. Eyepieces work in combination with microscope objectives to further magnify the intermediate image so that sample details can be observed.

| Type / magnification | WF/10x | WF/16x* |
|-------------------------------|----------|----------|
| Focal distance | 24.99 mm | 15.58 mm |
| Field of view diameter | 18 | 11 |

* Optional

3. TECHNICAL SPECIFICATIONS

| | | | | |
|---------------------|--|----------|---|----------|
| <i>Code</i> | 50128000 | 50128001 | 50128002 | 50128022 |
| <i>Model</i> | 128 | 128/1 | 128/2 | 128/22 |
| <i>Head</i> | Monocular | | Binocular | |
| <i>Eyepieces</i> | WF10x | | | |
| <i>Nosepiece</i> | Triple | | Quadruple | |
| <i>Objectives</i> | Achromatic 4x,10x,40x (R) | | Achromatic 4x,10x,40x (R), 100x (R)(I) | |
| <i>Focusing</i> | Independent coarse and fine adjustment knobs | | | |
| <i>Stage</i> | With clips | | Mechanical | |
| <i>Condenser</i> | Fixed, single lens N.A. 0.65 | | Abbe N.A. 1.25 | |
| <i>Diaphragm</i> | Iris | | | |
| <i>Illumination</i> | Halogen 6V, 20 W | | | |

Head

Depending on the model:

- Monocular, inclined 45° and rotary 360°
- Binocular, inclined 45° and rotary 360°. With interpupillary distance adjustment by free mechanism and dioptre compensation ring in one of the eyepieces holders.

Mechanical length of the tube

- 160 mm

Estateive

- Metallic, very stable, with independent coarse and fine adjustment knobs

Focusing travelling

- Coarse: 30 mm with stop limit
- Fine: 30 mm, 0-200 microns in 2 microns

Nosepiece

Depending on the model:

- Triple, clockwise and anticlockwise revolving nosepiece
- Quadruple, clockwise and anticlockwise revolving nosepiece

Stage

Depending on the model:

- Stage (120x120 mm) with clips
- Mechanical stage (120x120 mm). Movement range 60x30 mm.

Condenser

Located under the stage, it enables the light to be focused on the observed sample

Depending on the model:

- Single lens fixed condenser (N.A. 0.65), with iris diaphragm and movable filter holder

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 **ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.**

INDEX OF LANGUAGES

| | |
|---------------|-------|
| Spanish | 2-17 |
| English..... | 18-32 |

INDEX OF CONTENTS


| | |
|---|----|
| 1. USES OF THE INSTRUMENT | 19 |
| 2. DESCRIPTION..... | 19 |
| 3. TECHNICAL SPECIFICATIONS | 20 |
| 4. INSTALLATION / SETTING UP..... | 22 |
| 5. MAINTENANCE AND CLEANING..... | 22 |
| 6. CHOOSING OBJECTIVES AND EYEPIECES..... | 29 |
| 7. CAUSES OF A DEFECTIVE IMAGE | 29 |
| 8. TROUBLESHOOTING | 30 |
| 9. RECOMMENDATIONS | 31 |
| ANNEX I: CE CERTIFICATE | 32 |

1. USES OF THE INSTRUMENT

Thanks to their modern and developed design and to the modular configuration, Zuzi microscope series 128 becomes essential and ease-to-use equipment for laboratories in any education field.

2. DESCRIPTION

1. Eyepiece
2. Rotary head
3. Objectives
4. Estative
5. Focus limit screw
6. Coarse adjustment knob
7. Fine adjustment knob
8. Stage
9. Clips
10. Condenser with iris diaphragm
11. Precondenser
12. Switch
13. Light intensity adjustment knob

Thank you for choosing this equipment. We sincerely wish that you enjoy your Zuzi education microscope 128. We highly recommend looking after this equipment according to what is stated in this manual.

Zuzi develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

Please bear in mind the following:

- ◆ This manual is inseparable from the Zuzi education microscope model 128, so it should be available for all the users of this equipment.
- ◆ You should carefully handle the microscope avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- ◆ Never dismantle the different pieces of the microscope to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the microscope.
- ◆ To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting Zuzi Technical Assistance Department by email to asistencia@auxilab.es or by telephone: +34 807 117 040 (0.30 Euros/min).
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the microscope warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the microscope warranty.
- ◆ Neither fuses (1 A) nor accessories (including their loss), are covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Zuzi Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Please do not forget filling the warranty certificate and send it before 15 days after the date of purchase.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.

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