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◆ Never try to repair the centrifuge by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the centrifuge (burns, hurts...) or damages in nearby equipments.

◆ In the event of breakdown please contact your distributor to overhaul through Nahita Technical Assistance Department.

Cleaning

◆ Check and clean your centrifuge at regular intervals or immediately after spills occur to prevent premature wear and safety hazards

◆ Never use scourers or substances that can grate for cleaning metallic parts such as stainless steel, aluminium, coatings, etc. as they damage the centrifuge and produce an early ageing of the equipment.

◆ Wipe the rotor with a soft cloth moistened with a neutral detergent with extreme care and remove the detergent with a clean clamp cloth.

6.RESOLUCIÓN DE PROBLEMAS

| <i>FAILURE</i> | <i>PROBABLE CAUSE</i> | <i>TROUBLESHOOTING</i> |
|--|---|--|
| The centrifuge won't run | 1. The power cord is unplugged 2. Foreign objects wind around the rotor shaft 3. The motor is near the end of its expected life | 1. Plug the power cord in 2. Remove the foreign objects after unplugging the centrifuge 3. Call Nahita technical assistance department |
| Abnormal noise The centrifuge vibrates abnormally | 1. Imbalanced loads 2. Foreign objects are sticking on the rotor or chamber 3. Defective motor | 1. Balance the load properly 2. Remove the foreign articles after unplugging the centrifuge 3. Call Nahita technical assistance department |
| The rotor won't be installed properly | Damaged or deformed rotor or motor shaft boss | Call Nahita technical assistance department |

INSTRUCTIONS ON ENVIRONMENT PROTECTION

At the end of its life cycle, please, does not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

se it directly to the sun. Avoid vibrations, dust and dry environments.

- ◆ During its functioning dangerous materials such as flammable or pathological substances must be out of the safety area.
- ◆ This centrifuge is not equipped with an isolation system to contain aerosols in the event a tube leaks or breaks during a run.
- ◆ Never centrifuge explosive or flammable substances. These can ignite during the run and will cause the centrifuge to explode.
- ◆ Do not overfill the tubes. Tubes more than 70% full will cause the sample leaking during the run. This will expose you and the centrifuge to sample aerosols that may make you ill or may damage the centrifuge.
- ◆ Use caps for tubes. Tubes without caps will cause the sample leaking during the run, this will expose you and the centrifuge to sample.
- ◆ When you are not using the centrifuge for a long period of time please make sure it is unplugged in order to avoid possible accidents.
- ◆ It is essential to have the equipment switched off and unplugged from the net before cleaning, checking components or replacing any piece (e.g. replacement of a fuse).
- ◆ Never try to repair the centrifuge by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the equipment (burns, hurts...).
- ◆ Try not to spill any liquid on the control panel, though it is properly insulated. In case you have any doubts do please immediately unplug the centrifuge (see Maintenance).
- ◆ Made under the European regulations for electrical security, electromagnetic compatibility and security on machines.


5. MAINTENANCE AND CLEANING

To get the best results and a higher duration of this equipment it is essential to follow the processes of use.

Note: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.

- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of this equipment.
- ◆ Visually inspect the rotor for damage such as dents or chips from being dropped. If the rotor is damaged, do not use the centrifuge and call Nahita technical assistance department.
- ◆ Inspect the chamber for loose debris and excess water and remove as necessary.
- ◆ Inspect the lid is able to open and close smoothly, and inspect the lid is locked when pressing down.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ The centrifuge is supplied with a Schuko standard wire. It has to be plugged to an earth connection and the socket should be handy and ready to unplug the equipment in case of emergency.

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
will extend its working-life).

Operation

- ◆ To open the lid, press the [Open] button. Place the needed rotor and fix it to the shaft with the allen wrench.

Attention: Be sure to seat and secure the right rotor properly before operation. Improperly installed rotor may come off the motor shaft and may injure someone.

- ◆ Connect the centrifuge to the net and turn the power on.
- ◆ To load tubes on the rotor place them symmetrically into the rotor. Imbalance load will significantly reduce the life of the motor and centrifuge. In severe case, imbalances may damage the centrifuge and injure the operator.

- ◆ Once the tubes are load, close the lid firmly. Be sure that is completely closed before running the centrifuge.
- ◆ Use the [^] and [v] buttons to set the work speed. LED display will show the selected speed (1000x rpm).
- ◆ Press [] and the centrifuge will run.

Note: When the centrifuge is running, the speed can not be modified.

- ◆ To stop the centrifuge press [⏻] button, and the rotor will stop in some seconds. Pressing the [open] button, the rotor automatically will stop.
- ◆ When the rotor has stopped spinning, open the lid and remove the tubes.
- ◆ To turn the power off switch the button off and disconnect from the net.

Security

- ◆ The centrifuge must be used by previously qualified staff that knows how the equipment works thanks to the user manual.
- ◆ You should put the centrifuge in a horizontal plane stable table, having a safety area of at least 30 cm per side.
- ◆ Follow all required applicable laws about safety and health (electrical security devices, inflammable, toxic and/or pathological products, etc.).
- ◆ Do not place the centrifuge near any warm supply (burners, blowlamps, etc), nor expo-

4. INSTALLATION / SETTING UP

Preliminary inspection

- ◆ Unwrap the centrifuge, take off the involving plastic and take off the poliespan protection in which it comes fitted.
- ◆ Please make sure that it does not present any damage because of the shipment. If it does please communicate to your transport agent or provider immediately, so that they can make the claims in the correct time limit.

We will only accept any equipment return within 15 days after delivery and provided it comes in its original package and with all accessories and documents supplied

- ◆ Please check that all the accessories are enclosed with the equipment:
 - Hematocrit rotor
 - Tube-strip rotor
 - 6 tube holes angle rotor
 - 4 tube holes angle rotor
 - Tube adaptors
 - Allen wrench
 - Schuko wire
 - User manual
 - Microhematocrit reader

Installation

- ◆ Before using this instrument, it is convenient for you to familiarize with its components and basic essentials.

PLEASE READ THOROUGHLY THE INSTRUCTIONS BEFORE CONNECTING AND OPERATING WITH THIS EQUIPMENT

- ◆ Please put the centrifuge on top of a horizontal, plane and stable table making a free space at least at 30 cm per side. Do not put the centrifuge near any warm supply (burners, blowlamps...), nor expose it directly to the sun, etc.
- ◆ Avoid inflammable or toxic substances in the working area.
- ◆ The environmental conditions required are: room temperature: 5 to 35 °C, atmospheric pressure: 70 to 106 kPa and relative humidity: 30 to 85 %.

Please insert the wire that feeds the AC electric current in the base of current 100-220V 50-60Hz provided with earth wire and to the other end to the centrifuge connector.

Neither the manufacturer nor the distributor will assume any responsibility for the damages produced to the equipment during its installation or damages to persons suffered by the improper use of the electric connection. The tension should be 100-220V 50-60Hz.

- ◆ If you are not using the centrifuge for a long period of time please make sure it is disconnected from the net and protected from dust (this way you will avoid accidents and

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

*With adapters provided

ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

INDEX OF LANGUAGES

| | |
|---------------|-------|
| SPANISH | 2-9 |
| ENGLISH..... | 10-16 |

INDEX OF CONTENTS

| | |
|-----------------------------------|----|
| 1. USES OF THE INSTRUMENT..... | 11 |
| 2. DESCRIPTION | 12 |
| 3. TECHNICAL SPECIFICATIONS..... | 12 |
| 4. INSTALLATION / SETTING UP..... | 13 |
| 5. MAINTENANCE AND CLEANING..... | 14 |
| 6. TROUBLESHOOTING..... | 15 |

1. USES OF THE INSTRUMENT

This mini centrifuge Nahita Fugevet with an elegant appearance, compact design, make the product operation very smooth, quiet even under high speed. Particular cover switch design, pressing the opening button easily will open. It automatically stops the motor before open if it is running.

This centrifuge is equipped with three types of rotors: angular, strips and micro-hematocrit tubes and their applications are mainly oriented to analysis needs in veterinary fields

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
Thank you for choosing this equipment. We sincerely wish that you enjoy your mini centrifuge Nahita Fugevet. We highly recommend looking after this equipment according to what is stated in this manual.

Nahita develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

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Please bear in mind the following:

- ◆ This manual is inseparable from the mini centrifuge Nahita Fugevet, so it should be available for all the users of this equipment.
- ◆ Follow all required applicable laws about safety and health (electrical security devices, inflammable, toxic and/or pathological products, etc.).
- ◆ You should carefully handle the centrifuge avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- ◆ Never dismantle the different pieces of the centrifuge to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the centrifuge.
- ◆ To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting Nahita Technical Assistance Department by email to asistencia@auxilab.es.
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the centrifuge warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the centrifuge warranty.
- ◆ Neither fuses nor accessories (including their loss), are covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Nahita Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.