DECLARACIÓN DE CONFORMIDAD CE
CE DECLARATION OF CONFORMITY
DÉCLARATION DE CONFORMITÉ CE

El fabricante | The manufacturer | Le fournisseur:

AUXILAB, S.L.

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Thank you for choosing this equipment. We sincerely wish that you enjoy your Refrigerator Incubator NAHITA. We highly recommend looking after this equipment according to what is stated in this manual.

NAHITA develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

Please bear in mind the following:

- ◆ This manual is inseparable from the refrigerator incubator Nahita, so it should be available for all the users of this equipment.
- ◆ You should carefully handle the refrigerator incubator Nahita avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- ◆ Never dismantle the different pieces of the refrigerator incubator Nahita to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the Refrigerator Incubator Nahita.
- ◆ To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting marca Technical Assistance Department by email to asistencia@auxilab.es
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the Refrigerator Incubator Nahita warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the Refrigerator Incubator Nahita warranty.
- ◆ Neither fuses (1A) nor accessories (including their loss), are covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to NAHITA Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Please do not forget filling the warranty certificate and send it before 15 days after the date of purchase.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.

INSTRUCTIONS ON ENVIRONMENT PROTECTION

- At the end of its life cycle, please, do not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances.
- It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.
- The materials are recyclable as mentioned in its marking.
- By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.
- Please inquire at the community administration for the authorized disposal location.
- The R134A is Biodegradable.

6. RESOLUTION PROBLEMS

| <i>PROBLEM</i> | <i>CAUSE</i> | <i>SOLUTION</i> |
|-------------------------|--|--|
| Not work | 1 The plug is not connected 2 No power supply 3 The current limiter tripped breaker | 1 Make sure the cable is properly connected 2 Check the fuse and replace if necessary 3 Check the current limiter |
| The team vibrate | Not properly leveled | To level correctly |
| Unstable T | 1 Door uncorrectly close 2 Near heat sources 3 Cable Lid is not in its position | 1 Check door 2 Do not open the doors frequently 3 Remove of heat sources 4 Put lid |
| The door can not close | 1 Uncorrect position Handle 2 Inner Obstruction 3 Habitual uncorrect close | 1 Check handle 2 Replace inner elements 3 Close the door gently |
| Lamp is not working | 1 Uncorrectly place 2 Lamp has gone 3 No Switch connected | 1 Put the lamp correctly 2 Change lamp 3 Connect the switch |
| Turn on but not working | The schedule is uncorrect RUN/STOP | Schedule again and press Down/Setting |

ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

INDEX OF LANGUAGES

| | |
|---------------|-------|
| Spanish | 2-11 |
| English | 12-21 |
| French..... | 22-31 |

INDEX OF CONTENTS

| | |
|------------------------------------|----|
| 1. USES OF THE INSTRUMENT | 13 |
| 2. DESCRIPTION | 14 |
| 3. TECHNICAL SPECIFICATIONS..... | 15 |
| 4. INSTALLATION / SETTING UP | 16 |
| 5. MAINTENANCE AND CLEANING | 19 |
| 6. RESOLUTION PROBLEMS..... | 20 |
| ANNEX I: CE CERTIFICATE..... | 32 |

1. USES OF THE INSTRUMENT

Nahita refrigerated incubators series 639 are especially designed to be used in those processes that take place at moderate temperatures, between 5 and 60 °C, and that requires a good stability and precise control of the temperature.


Such processes are commonly developed in laboratories of cosmetics, microbiology, food analysis, biotechnology, botany or biochemistry and involve, among their different applications, BOD determination, measurement of the enzymatic activity, serum and plasma assays, conservation or incubation of samples. Applications and fields of use of Nahita incubators series 639 are thus, very different and extended.


2. DESCRIPTION

- 1.1 Control panel
- 1.2 General switch
- 1.3 Inner glass door
- 1.4 Recipient surplus liquid

- 1.5 Socket
- 1.7 Ventilation grille
- 1.8 Auxiliary Plug
- 1.9 Outlet cable for additional equipment

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
- 2.1 Digital display
- 2.2 Illuminating: Button on / off interior light
- 2.3 Setting: Programming Button

- 2.4 Add / Screen color: Button rise / light digital display
- 2.5 Down / Operating: Button reduction / Commissioning

- ◆ Never try to repair the refrigerator incubator Nahita by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the equipment (burns, hurts...).
- ◆ Try not to spill any liquid on the control panel, though it is properly insulated. In case you have any doubts do please immediately unplug the refrigerator incubator Nahita (see Maintenance).
- ◆ Made under the European regulations for electrical security, electromagnetic compatibility and security on machines.
- ◆ The gas R134A is inflammable so please have precautions if the equipments suffer damage pleased contact with your distributor.
- ◆ Pleased Turn off the equipment to clean or pieces replacement.
- ◆ The equipment have a Schuko cable, pleased have this connecting in an accessible place to remove in a emergency case.

5. MAINTENANCE AND CLEANING

To get the best results and a higher duration of this equipment it is essential to follow the processes of use.

Note: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.

- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of this equipment.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ The refrigerator incubator Nahita is supplied with a Schuko standard wire. It has to be plugged to an earth connection and the socket should be handy and ready to unplug the equipment in case of emergency.
- ◆ Never try to repair the refrigerator incubator Nahita by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the Refrigerator Incubator Nahita (burns, hurts...) or damages in nearby equipments.
- ◆ In the event of breakdown please contact your distributor to overhaul through Nahita Technical Assistance Department.

Cleaning

- ◆ Never use scourers or substances that can grate for cleaning metallic parts such as stainless steel, aluminium, coatings, etc. as they damage the Refrigerator Incubator NAHITA and produce an early ageing of the equipment.
- ◆ Use a fluff-free cloth dampened with soaped water that does not contain abrasives.

section TIME is showing the passed time, which begins upon reaching the T set.

Note: If you not set a specific time (Time Set 00:00), the incubators works in continuous (in this case PV shows the T real and SV (3.2) the T-programmed).

COOL(3.4) or HEAT (3.5) indicators mark the heating and cooling states respectively in terms of the T set.

If during the operation of the equipment the actual temperature deviates by ± 5 ° C of the T set, the ALM indicator will be activated and notified of the situation.

◆ An acoustic signal marking the end of scheduled time ("End" in TIME display section) and the equipment stops run (STOP ON). Pressing any key will stop the sound, otherwise it will after 30s.

◆ Pressing the Down / Operating (2.5) button for 3 to 5 seconds, the equipment will restart with the preselected program. To reset the equipment, proceed from the point 2 to 5.

Note: Your equipment can be reprogrammed at any time from the point 2 to 4. without stopping the machine.

◆ After finish run, turn off the equipment pressing general switch (1.2). This switch leaves without additional power supply plug (1.8)

The incubator is equipped with an auxiliary plug (1.8) outside grounded for connecting a stirrer or similar equipment through the hole (1.9). Do not connect any equipment with heating function.

To watch and control samples, open the outer door and turn on the light by pressing Illuminating (2.2) of the control panel. Keep the inner glass door closed to keep the temperature fluctuations.

Attention! When you press the general switch or in case of cut / power recovery, the equipment automatically starts and the previous programming starting from 00.00 up to time.

Security

◆ The refrigerator incubator Nahita must be used by previously qualified staff that knows how the equipment works thanks to the user manual.

◆ You should put the refrigerator incubator Nahita in a vertical position, with the grilles free from obstructions.

◆ Do not place the refrigerator incubator Nahita near any warm supply (burners, blowlamps, etc), nor expose it directly to the sun. Avoid vibrations, dust and dry environments.

◆ During its functioning dangerous materials such as flammable or pathological substances must be out of the safety area.

◆ When you are not using the refrigerator incubator Nahita for a long period of time please make sure it is unplugged in order to avoid possible accidents.

◆ It is essential to have the equipment switched off and unplugged from the net before cleaning, checking components or replacing any piece (e.g. replacement of a fuse).

Digital screen

3.1 PV: Real temperature

3.2. Temperature programmed

3.3 TIME: Passed time

3.4 ALM: Alarm indicator

3.5 COOL: Indicator cooling function

3.6:HEAT: Indicator warming function

3.7 STOP / RUN: Stop/ Indicator

3. TECHNICAL ESPECIFICATION

| | | | |
|----------------------|--------------------------------|---------------|--------------|
| Code | 50639070 | 50639150 | 50639250 |
| Capacity | 70 L | 150 L | 250 L |
| Timer | 99 h y 59' (pasos de 1 h / 1') | | |
| Temperature range | 5 °C a 60 °C | | |
| Resolution | 0.1 °C | | |
| Accuracy | ±0.5 °C | | |
| Homogeneity | ±1 °C | | |
| N ° heights | 7 (regulables cada 5 cm.) | | |
| Int. dimens. (HxLxW) | 42x35x50 cm. | 50x40x75 cm. | 50x50x95 cm |
| Ext. dimens. (HxLxW) | 57x56x120 cm. | 64x62x127 cm. | 66x70x147 cm |
| Power | 300 W | 500 W | 800 W |
| Weight | 59 Kg | 68 Kg | 71 Kg |
| Power | 220-240 V, 50-60 Hz ± 10% | | |

Main technical Benefit:

◆ The external metallic structure is covered with epoxy and the interior chamber is made of resistant to acids and corrosive agents and easy-to-clean stainless steel.

◆ Solid external door with aperture of 180° and inner transparent thermal glass door.

◆ Internal illumination by a fluorescent lamp with independent switch that allows controlling sample evolution through the glass door without variation of the temperature in the internal chamber.

◆ Provided with an external socket and a hole to introduce a cable, probe, etc., that allow the connection and operation of additional equipments inside the chamber, thus increa-

- sing the applications of the incubator.
- ◆ Supplied with height adjustable grill type shelves (2 units), which provides a high flexibility of use since it allows introducing flasks of different size.
 - ◆ Transverse flow forced air circulation system that provides a maximum uniformity and homogeneity of the temperature in the inner chamber.
 - ◆ Parameter control by microprocessor
 - Control panel at the upper side with tactile keyboard and LCD/LED digital display and back-lit general switch.
 - Height adjustable legs for an easy installation and levelling.
 - Protection system against overheating with visual and acoustic signal.
 - Hermetic and very silent compressor mounted on shock absorbers and evaporator with automatic defrost.
 - R134a biodegradable cooling gas, CFC and HCFC-free.

4. INSTALLATION / SETTING UP

Preliminary inspection

- ◆ Unwrap the Refrigerator Incubator NAHITA, take off the involving plastic and take off the polispán protection in which it comes fitted. Take off all the protective items and, without connecting the Refrigerator Incubator NAHITA to the net, make sure that it does not present any damage because of the shipment. In case the Refrigerator Incubator NAHITA presents any damage tell it immediately to your transport agent or dealer so that they can make the claims in the correct time limit.
- ◆ Please keep the original wrapping; you will always need it for returns enclosed with all the accessories supplied.
- ◆ Please check that all the accessories are enclosed with the equipment:

We will not accept any equipment in return period unless it comes in its original wrapping.

Accessories:

- 2 shelves grilles and 4 guides
- 1 Compact medium shelves (27x18 cm)
- Schuko standar Cable
- Fuse
- Warranty
- Instruction manual

Installation

Before using this instrument, it is convenient for you to familiarize with its components and basic essentials.

PLEASE READ THOROUGHLY THE INSTRUCTIONS BEFORE CONNECTING AND OPERATING WITH THIS EQUIPMENT

- ◆ Please put the refrigerator incubator Nahita in vertical position and with the grille (1.7) with any obstruction.
- ◆ Do not put the refrigerator incubator Nahita near any warm supply (burners, blowlamps...), nor expose it directly to the sun, etc.
- ◆ Avoid inflammable or toxic substances in the working area.
- ◆ The refrigerator incubator Nahita is supplied with a Schuko standard wire
- ◆ Pleased put the shelves in the wanted position
- ◆ Place the compact tray at the bottom of the equipment so that collects the defrost liquid through plastic tube.
- ◆ Please insert the wire that feeds the AC electric current in the base of current 220V 50Hz $\pm 10\%$ provided with earth wire and to the other end to the Refrigerator Incubator NAHITA connector (1.5).

Neither the manufacturer nor the distributor will assume any responsibility for the damages produced to the equipment during its installation or damages to persons suffered by the improper use of the electric connection. The tension should be 220V 50Hz $\pm 10\%$.

If you are not using the refrigerator incubator Nahita for a long period of time please make sure it is disconnected from the net and protected from dust (this way you will avoid accidents and will extend its working-life).

NOTE: *Pleased do not connect the refrigerator incubator Nahita with the other equipment in the same terminal can suffer a over warming .*

Implementation

Before first use or change the location, do not turn on the incubator until they have spent 2 h for oils and gases from the engine to be stabilized.

Adjust of the temperature and run time:

- ◆ General switch (1.2). Wait until the equipment make automatically the self-adjustment.
- ◆ Press Setting (2.3). On screen shows "Temp Set" and the command "Sp" for selecting the operating temperature by adjusting key parameters (\blacktriangle \blacktriangledown).
- ◆ Press Setting (2.3) again. The screen shows "Time Set" and the command "St" for selection of run time (range 00-99 selection of 00-h and 1 h 59'en steps 1 ') through the keys (\blacktriangle \blacktriangledown).
- Use the Setting button (2.3) to move from hours to minutes.
- ◆ Press again Setting (2.3) to record the parameters.
- ◆ Finally, press and hold Down / Operating for 3-5 s until the incubator start run (RUN ON).

On screen, the PV section shows the actual temperature inside the incubator, while the