ANNEX I: CE CERTIFICATE

AUXILAB S.L.

**DECLARATION OF CONFORMITY
CONDUCTIVITY METER 910/5 of Auxilab, S.L for the
Directive of Machines (89/392/CEE modified) and the regulations
adopted for their transposition**

NAME OF THE MANUFACTURER / IMPORTER:
AUXILAB, S.L.

ADDRESS: **Polígono Morea Norte, 8
31191 Beriáin (Navarra)**

WE STATE THAT:

**NAHITA CONDUCTIVITY METER 910/5
CODE 50910050**

Are designed and manufactured according to:

- ♦ Directive 89/392/CEE, including the modifications and national regulations that transpose them.
- ♦ Directive 89/336/CEE modified over electromagnetic compatibility.
- ♦ And that the following harmonized rules have applied (or part of them):
UNE 292-1/-2/-2/A1, UNE-EN 1050, UNE-EN 614-1, UNE-EN 1037, UNE-EN 1088, UNE-EN 547, UNE-EN 953, UNE-EN 294, UNE-EN 418, UNE-EN 894-1, UNE-EN 894-2, UNE-EN 954-1, UNE-EN 60204-1, UNE 61010-1/A2, UNE-EN61010-2-051.

BERIAIN 19th JUNE 2007

Signed by: **ALFONSO AINCIBURU SANZ**
DIRECTOR/MANAGER

Polígono Morea Norte, 8 31191 Beriain (Navarra) - Spain. Tel. 948 310 513 Fax 948 312 071
Internet: www.auxilab.es · Email: correo@auxilab.es

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This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.

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Accessories

Conductivity meter Nahita 910/5 is supplied without a conductivity cell. The following models are available in our catalogue.

Code	Constant	Range
90908010	$0.1 \pm 0.02 \text{ cm}^{-1}$	$0.1 \mu\text{S}-300 \text{ mS/cm}$
90908020	$1 \pm 0.2 \text{ cm}^{-1}$	$1 \mu\text{S}-100 \text{ mS/cm}$
90908030	$10 \pm 2 \text{ cm}^{-1}$	$300 \mu\text{S}-1000 \text{ mS/cm}$

5. MAINTENANCE AND CLEANING

To get the best results and a higher duration of the conductivity meter it is essential to follow the processes of use.

Note: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.

- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of this equipment.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ Never try to repair the conductivity meter by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the conductivity meter (burns, hurts...) or damages in nearby equipments.
- ◆ **Maintenance of the conductivity cell.** Rinse the conductivity cell in clean water after each measurement; for a thorough cleaning, clean the cell with a cloth or a non abrasive detergent. After cleaning the cell, calibrate the equipment.
- ◆ In the event of breakdown please contact your distributor to overhaul through Nahita Technical Assistance Department.

Cleaning

- ◆ Never use scourers or substances that can grate for cleaning the case of the equipment as they damage the conductivity meter and produce an early ageing of the equipment.
- ◆ Use a fluff-free cloth dampened with soaped water that does not contain abrasives.

INSTRUCTIONS ON ENVIRONMENT PROTECTION

At the end of its life cycle, please, do not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

- ◆ Press ENTER/PRINT (1.4) and all the data saved in that measuring mode will be printed.

Note: If there are no data saved in that measuring mode, "NULL" will be displayed and after a few seconds the instrument will return to the measuring mode.

5. Auto-off function

Nahita conductivity meter 910/5 presents the possibility to programme and auto-off function that will turn the equipment off after 2 hours of use. For this purpose:

- ◆ In any measuring mode, press MODE/MEAS (1.3).
- ◆ With ▲ /C/T/S (1.5) or ▼ /SAVE (1.6) select "ASD" on display.
- ◆ Press ENTER/PRINT (1.4).
- ◆ With ▲ /C/T/S (1.5) select ON (to activate the auto-off function) or OFF (to inactivate the auto-off function).

Note: By default, the equipment presents the auto-off function activated.

- ◆ Press ENTER/PRINT (1.4) to confirm the selection.
- ◆ Press MODE/MEAS (1.3) to return to the measuring mode.

Note: The equipment also presents a low battery alert. When the batteries are lack of energy, a sign of low battery will be displayed to remind the user to change the batteries.

If you are not using the conductivity meter for a long period of time please remove the batteries and protected the equipment from dust (this way you will avoid accidents and will extend its working-life).

Security

- ◆ The conductivity meter must be used by previously qualified staff that knows how the equipment works thanks to the user manual.
- ◆ Put the conductivity meter on top of a horizontal, plane and stable table making a free space at least at 30 cm per side.
- ◆ Do not put the conductivity meter near any warm supply (burners, blowlamps...), nor expose it directly to the sun, etc. Avoid vibrations, dust and dry environments.
- ◆ During operation, dangerous material such as flammable liquids or pathological material must be out of the working area.
- ◆ It is essential to have the equipment switched off before cleaning, checking components or replacing any piece.
- ◆ Never try to repair the conductivity meter by yourself, since you will lose the warranty and may provoke damages to the general operating system, as well as injuries to the people that usually handle the equipment (burns, hurts...).
- ◆ Try not to spill any liquid on the control panel, though it is properly insulated. In case you have any doubts do please immediately unplug the equipment (see Maintenance).
- ◆ Made under the European regulations for electrical security, electromagnetic compatibility and security on machines.

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

- ◆ Press MODE/MEAS (1.3).
- ◆ With ▲ /C/T/S (1.5) or ▼ /SAVE (1.6) select "VIEW" on display.
- ◆ Press ENTER/PRINT (1.4), so as the equipment enters the data visualization mode.
- ◆ The last result that has been saved will be displayed.

Note: If there are no saved data at that measuring mode, "NULL" will be displayed.

- ◆ Use ▲ /C/T/S (1.5) to ▼ /SAVE (1.6) to see in turn all the saved data. The position number of each datum will be shown at the lower left side of the screen.
- ◆ Press ENTER/PRINT (1.4) to exit the data visualization mode.
- ◆ Press MODE/MEAS (1.3) to return to the measuring mode.

3. Data deleting

This function allows deleting all the data saved in a particular measuring mode:

- ◆ Select the desired measuring mode by pressing ▲ /C/T/S (1.5).
- ◆ Press MODE/MEAS (1.3).
- ◆ With ▲ /C/T/S (1.5) or ▼ /SAVE (1.6) select "DEL" on display.
- ◆ Press ENTER/PRINT (1.4) so as the equipment enters the data deleting mode.
- ◆ "SURE" will be displayed
- ◆ Press ENTER/PRINT (1.4) and all the data saved in that measuring mode will be deleted.

Note: Only the data saved in that measuring mode will be deleted. To delete the data saved in any other measuring mode, firstly select that mode and then proceed as previously described.

- ◆ While data are deleted "DEL--" will be displayed.
- ◆ Once the process has finished "DEL" will be displayed.
- ◆ Press MODE/MEAS (1.3) to return to the measuring mode.

4. Data printing

The results obtained with the conductivity meter 910/5 can be printed with a TP-16 printer. For this purpose, connect the conductivity meter with the printer through the interface RS232.

To print the current result:

- ◆ In measuring mode used, press ENTER/PRINT (1.3).
- ◆ While the current result and the temperature of the solution analyzed are printed, "PRN" will be displayed at the lower left side.
- ◆ Once printing is finished, the instrument returns to the measuring mode.

To print the saved data:

- ◆ Select the measuring mode whose data are going to be printed.
- ◆ Press MODE/MEAS (1.3).
- ◆ With ▲ /C/T/S (1.5) or ▼ /SAVE (1.6) select "PRN" on display.
- ◆ Press ENTER/PRINT (1.4).
- ◆ "SURE" will be displayed.

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Referencia **50910050**

Conductividad

Rango	0.00-19.99 μ S/cm; 20.0-199.9 μ S/cm; 200-1999 μ S/cm; 2.00-19.99 mS/cm; 20.0-199.9 mS/cm
Resolución	0.01 μ S/cm; 0.1 μ S/cm; 1 μ S/cm; 0.01 mS/cm; 0.1 mS/cm

standard solution to be (25.0±0.1) °C.

- ◆ Press MODE/MEAS (1.3).
- ◆ With ▲ /C/T/S (1.5) or ▼ /SAVE (1.6) select "TCAL" on display.
- ◆ Press ENTER/PRINT (1.4) to enter the TDS conversion factor adjustment mode.
- ◆ When the reading on display is stable, press ENTER/PRINT (1.4). "CA-M" will be displayed.
- ◆ With ▲ /C/T/S (1.5) (increase) or ▼ /SAVE (1.6) (decrease) introduce the exact value of TDS according to the TDS standard value chart.
- ◆ Press ENTER/PRINT (1.4) to save the value; after a few seconds "TCAL" will be displayed.
- ◆ Press MODE/MEAS (1.3) to return to the measuring mode.

Conductivity, TDS or salinity determinations

NOTE: before performing measurements, make sure the equipment is calibrated.

- ◆ Turn the conductivity meter on by pressing ON/OFF (1.2); by default, the instrument starts at the last measuring mode used.
- ◆ Clean the conductivity cell with distilled water and dry it with a filter or cellulose paper.
- ◆ Submerge the cell and the temperature probe 4 cm into the sample to be analyzed and stir it with the cell. Wait until the reading is stable. Conductivity, TDS or salinity measurement is displayed at the middle and the temperature is shown below.

Note: Once the temperature probe is connected to the equipment, the measurements will be compensated to the temperature of reference according to the temperature coefficient established. If the temperature probe is not connected, the value displayed will be the original value without being compensated. The reference temperature is 25 °C for conductivity and TDS determinations and 18 °C for salinity determinations.

- ◆ To go to any other measuring mode press ▲ /C/T/S (1.5).

NOTE: If several samples are going to be analyzed successively it is recommended to clean the conductivity cell and the temperature probe with distilled water to avoid cross contamination.

Other functions

1. Data saving

Nahita conductivity 910/5 can store up to 250 readings of each measuring mode. For this purpose, at any measuring mode press ▼ /SAVE (1.6) and the result will be automatically saved. While saving the result, "SAVE" will be displayed.

2. Data visualization

To view the data saved in any measuring mode:

- ◆ Select the desired measuring mode by pressing ▲ /C/T/S (1.5).

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- ◆ Press ENTER/PRINT (1.4), to save the value; after a few seconds "CAL" will be displayed.
- ◆ Press MODE/MEAS (1.3) to return to the measuring mode.

2. Calibration of the TDS conversion factor

- ◆ Select the Standard solution according to the character of the solution to be measured and the measuring range.

Relationship between conductivity and TDS of the standard solution

Conductivity (μS/cm)	TDS standard value		
	KCl (mg/L)	Conductivity (μS/cm)	442* (mg/L)
23	11.6	10.7	14.74
84	40.38	38.04	50.5
447	225.6	215.5	300
1413	744.7	702.1	1000
1500	757.1	737.1	1050
2070	1045	1041	1500
2764	1382	1414.8	2062.7
8974	5101	4487	7608
12880	7447	7230	11367
15000	8759	8532	13455
80000	52168	48384	79688

*442: 40% Na₂SO₄, 40% NaHCO₃, 20% NaCl

All the values shown in the chart are at 25 °C

- ◆ Calibrate the conductivity cell constant (See Calibration of the conductivity cell).
- ◆ Select the TDS measuring mode
- ◆ Clean the conductivity cell with distilled water.
- ◆ Submerge the conductivity cell in the standard solution. Control the temperature of the

Calibration**Calibration of the conductivity cell**

Each conductivity cell presents a specific constant value; if user thinks that the value indicated by the factory is not correct, proceed as follows:

- ◆ Select the right standard solution depending on the constant value of the cell:

KCl standard solution to be used depending on the cell constant value

Conductivity cell constant (cm ⁻¹)	0.01	0.1	1.0	10.0
Approx. concentration of KCl solution (mol/L)	0.001	0.001	0.01 or 0.1	0.1 or 1

Chart for preparation of the KCl standard solutions.

Concentration (mol/L)	KCL concentration (g/L) at 20°C
1	74.2650
0.1	7.4365
0.01	0.7440
0.001	Dissolve 100 mL of the 0.01 mol/L solution in 1 L

◆ Connect the conductivity cell to the equipment and disconnect the temperature probe; this way, the conductivity value displayed will be the original value without temperature compensation.

- ◆ Turn the conductivity meter by pressing ON/OFF (1.2)
- ◆ Clean the conductivity cell with distilled water.
- ◆ Submerge the conductivity cell in the standard solution. Control the temperature of the standard solution to be (25.0±0.1) °C, (20.0±0.1) °C, (18.0±0.1) °C or (15.0±0.1) °C.
- ◆ Select the level of the constant value (0.01, 0.1, 1.0, 5.0 or 10.0) (see Parameter adjustment; conductivity cell constant adjustment) and press ENTER/PRINT (1.4).
- ◆ Press MODE/MEAS (1.3) to return to the measuring mode.
- ◆ Press again MODE/MEAS (1.3)
- ◆ With ▲ /C/T/S (1.5) or ▼ /SAVE (1.6) select "CAL" on display.
- ◆ Press ENTER/PRINT (1.4) to enter the calibration mode.
- ◆ When the reading on display is stable, press ENTER/PRINT (1.4). "CA-M" will be displayed.
- ◆ With ▲ /C/T/S (1.5) (increase) or ▼ /SAVE (1.6) (decrease) introduce the exact value of conductivity according to the following chart:

Exact value of conductivity regarding the approximate concentration of the KCl standard solution.

Rango de conductividad (μS/cm)	Rango de resistividad (Ω cm)	Constante recomendada de la célula de conductividad (cm ⁻¹)
0.05-2	20x10 ⁶ -500x10 ³	0.01, 0.1
2-200	500x10 ³ -5x10 ³	0.1, 1.0
200-2000	5x10 ³ -500	1.0
2000-20000	500-50	1.0, 10
2000-2x10 ⁵	500-5	10

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- ◆ Press ENTER/PRINT (1.4) to save the value.
- ◆ Press MODE/MEAS (1.3) to exit the constant adjustment mode and return to the measuring mode.

Note: To cancel the constant adjustment process press MODE/MEAS (1.3)

2. TDS conversion factor adjustment

Nahita conductivity meter 910/5 takes by default a value of 0.50 for the TDS conversion factor. The user can adjust this value within a range from 0.40 to 1.00. For this purpose, proceed as follows:

- ◆ In the TDS measuring mode press MODE/MEAS (1.3).
- ◆ With ▲ /C/T/S (1.5) or ▼ /SAVE (1.6) select "TCOE" on display.
- ◆ Press ENTER/PRINT (1.4), so as the equipment enters the TDS conversion factor adjustment mode.
- ◆ Press ▲ /C/T/S (1.5) (increase) or ▼ /SAVE (1.6) (decrease) to adjust the desired exact value of the TDS conversion factor.
- ◆ Press ENTER/PRINT (1.4) to save the value.
- ◆ Press MODE/MEAS (1.3) to exit the conversion factor adjustment mode and return to the measuring mode.

Note: TDS conversion factor can only be adjusted from the TDS measuring mode.

3. Temperature coefficient adjustment

Conductivity determination depends on temperature. The temperature variations can be compensated by introducing a compensation factor (temperature coefficient) that will be taken into account by the equipment to calculate the conductivity at the reference temperature.


By default, the equipment takes a value of 0.020 (2 %/°C) for the temperature coefficient; if this needs to be modified, proceed as follows:

- ◆ In the Conductivity or TDS measuring mode press MODE/MEAS (1.3).
- ◆ With ▲ /C/T/S (1.5) or ▼ /SAVE (1.6) select "COEF" on display.
- ◆ Press ENTER/PRINT (1.4), so as the equipment enters the temperature coefficient adjustment mode.
- ◆ Press ▲ /C/T/S (1.5) (increase) or ▼ /SAVE (1.6) (decrease) to adjust the desired exact value of the temperature coefficient.
- ◆ Press ENTER/PRINT (1.4) to save the value.
- ◆ Press MODE/MEAS (1.3) to exit the temperature coefficient adjustment mode and return to the measuring mode.


Note: Temperature coefficient can only be adjusted from the Conductivity and TDS measuring mode. The temperature coefficient for salinity determinations is already fixed in the equipment and can not be modified.

- When selecting TDS mode, "TDS" will be displayed

- When selecting Salinity mode, "SALT" will be displayed

Note: to assure precise results, it is necessary to adjust and establish some parameters (see Parameter adjustment) before performing any reading.

Parameter adjustment

Nahita conductivity meter 910/5 allows the adjustment of different parameters to adapt to any kind of application and assure precise results:

- ◆ Conductivity cell constant (for conductivity, TDS and salinity determinations)
- ◆ TDS conversion factor (for TDS determinations)
- ◆ Temperature coefficient (for conductivity and TDS determinations)

1. Cell conductivity constant adjustment

Each conductivity cell presents its own constant value; this value must be introduced in the equipment every time the conductivity cell is changed:

- ◆ In any measuring mode (Conductivity, TDS or Salinity), press MODE/MEAS (1.3)
- ◆ With ▲ /C/T/S (1.5) or ▼ /SAVE (1.6) select "CONT" on display.
- ◆ Press ENTER/PRINT (1.4), so as the equipment enters the cell constant adjustment mode.
- ◆ Press ▲ /C/T/S (1.5) to select the level that corresponds to the constant value of the cell to be used (0.01, 0.1, 1.0, 5.0 ó 10.0).
- ◆ Press ENTER/PRINT (1.4) to save the selected level. "ADJS" will be automatically displayed.
- ◆ Press ▲ /C/T/S (1.5) (increase) or ▼ /SAVE (1.6) (decrease) to adjust the displayed value to the exact value of the cell constant (supplied by the manufacturer).

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B) Selection of the conductivity cell according to the TDS working range

TDS range (mg/L)	Conductivity range (µ S/cm)	Recommended constant value of the conductivity cell (cm ⁻¹)
0-1000	0-2000	1.0
1000-10000	2000-20000	1.0, 10
1000-19990	2000-40000	10

C) For salinity determinations it is recommended to use always a conductivity cell with a constant value of 10 cm⁻¹.

3. Connection of the temperature probe

4. The equipment is supplied with a temperature probe (1.8), connect this probe to the corresponding connector (2.2) of the equipment.

Turn the conductivity meter on by pressing ON/OFF (1.2). After a few seconds the instrument will be ready to operate.


Selection of the measuring mode

Nahita conductivity meter 910/5 can perform determinations of:

- a. Conductivity (µS/cm)
- b. TDS, total dissolved solids (mg/L)
- c. Salinity (%)

By pressing ▲ /C/T/S (1.5) select the desired measuring mode:

- When selecting Conductivity mode, "CON" will be displayed

they can make the claims in the correct time limit.

Please keep the original wrapping; you will always need it for returns enclosed with all the accessories supplied.

Please check that all the accessories are enclosed with the equipment:

- Temperature probe
- RS232 cable
- Aluminium case
- User's manual
- Warranty certificate

We will only accept any equipment return within 15 days after delivery and provided it comes in its original wrapping.

Installation

Before using the conductivity meter, it is convenient for you to familiarize with its components and basic essentials, as well as with its control functions.

PLEASE READ THOROUGHLY THE INSTRUCTIONS BEFORE CONNECTING AND OPERATING WITH THIS EQUIPMENT WITH THE AIM OF ACHIEVING THE HIGHEST FEATURES AND THE MAXIMUM DURABILITY OF THE INSTRUMENT.

1. Placing the batteries

Nahita conductivity meter 910/5 operates with 4 alkaline batteries type AA:

- ◆ Open the battery compartment; for this purpose, unscrew the four screws placed at the back side of the instrument.
- ◆ Insert the batteries according to their polarity.
- ◆ Put the lid of the compartment and screw the four screws.

NOTE: Remove the batteries if you are not using the instrument for a long period of time.

2. Connection of the conductivity cell

- ◆ Connect the conductivity cell (1.7) (not included) into the corresponding connector (2.1) of the equipment.

Selection of the conductivity cell

To adapt to any application, Nahita conductivity meter 910/5 can be used with conductivity cells of different constant value (see Accessories) allowing to adjust the constant value of the cell to be used at any moment.

Different charts are shown below to indicate the adequate constant of the conductivity cell depending on the working range

- A) Selection of the conductivity cell according to the conductivity working range:

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
4. INSTALLATION / SETTING UP

Preliminary inspection

Unwrap the conductivity meter (it is supplied in an aluminium case for its transport) and make sure that the equipment does not present any damage because of the shipment. In case it presents any damage, tell it immediately to your transport agent or dealer so that

- 1.1 Liquid cristal display (LCD)
- 1.2 ON/OFF
- 1.3 MODE/MEAS
- 1.4 ENTER/PRINT
- 1.5 ▲ /C/T/S
- 1.6 ▼ /SAVE
- 1.7 Conductivity cell (not included)
- 1.8 Temperature probe

3. TECHNICAL SPECIFICATIONS

Code	50910050
Conductivity	
Range	0.00-19.99 μ S/cm; 20.0-199.9 μ S/cm; 200-1999 μ S/cm; 2.00-19.99 mS/cm; 20.0-199.9 mS/cm
Resolution	0.01 μ S/cm; 0.1 μ S/cm; 1 μ S/cm, 0.01 mS/cm; 0.1 mS/cm

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INDEX OF CONTENTS

1. USES OF THE INSTRUMENT	19
2. DESCRIPTION	19
3. TECHNICAL SPECIFICATIONS	20
4. INSTALLATION / SETTING UP	21
5. MAINTENANCE AND CLEANING	31
ANNEX I: CE CERTIFICATE	32

1. USES OF THE INSTRUMENT


The Nahita conductivity meter 910/5 is ideal for field determinations in the analysis of river, spring or cultivation waters. It is supplied in a sturdy case for an easy transport of all the accessories of the meter and moreover, presents a degree of protection IP65 (protection against dust penetration and water jets) and a membrane keypad with 5 button to enter all the functions of the equipment.

The conductivity meter 910/5 determines conductivity, total dissolved solids (TDS) and salinity at the same time and furthermore, is equipped with a temperature probe for determination and automatic compensation of this parameter. It can be used with conductivity cells with a constant value of 0.01, 0.1, 1, 5 and 10 cm⁻¹ according the specific necessities of each application. Several parameters such as conductivity cell constant, temperature coefficient and TDS conversion factor can also be modified depending on the type of the sample to be analyzed or the conductivity cell used.

2. DESCRIPTION

Connections

- 2.1 Conductivity cell connector
- 2.2 Temperature probe connector
- 2.3 RS232 interface

Thank you for choosing this equipment. We sincerely wish that you enjoy your Nahita portable conductivity meter 910/5. We highly recommend looking after this equipment according to what is stated in this manual.

Nahita develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

Please bear in mind the following:

- ◆ This manual is inseparable from the Nahita portable conductivity meter 910/5, so it should be available for all the users of this equipment.
- ◆ You should carefully handle the equipment avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- ◆ Never dismantle the different pieces of the equipment to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the conductivity meter.
- ◆ To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting Nahita Technical Assistance Department by email to asistencia@auxilab.es or by telephone: +34 807 117 040 (0.30 Euros/min).
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the equipment warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the warranty.
- ◆ Accessories, including their loss are not covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Nahita Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.

ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

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