

FILTRES RS extraction thimbles



Food industry



Atmospheric pollution analysis



Atmospheric pollution analysis

- 1 | Made of pure fat-free cellulose cotton in a unique piece.
- 2 | High mechanical resistance.
- 3 | For Soxhlet, Kurnagawa extraction apparatus and other extraction equipments.
- 4 | Application: food industry, atmospheric pollution analysis, asphalt analysis, etc. for solid-liquid extraction.

Code	Ø int. (mm)	Ø ext. (mm)	Heigh (mm)	Box of (uds.)
85001430	14	16	30	25
85001990	19	22	90	25
85001860	18	22	60	25
85002260	22	26	60	25
85002380	23	27	80	25
85026100	26	30	100	25
85002980	29	34	80	25
85003380	33	36	80	25
85003394	33	36	94	25
85033100	33	36	100	25
85034130	34	37	130	25
85037123	37	41	123	25
85038150	38	41	150	25
85041123	41	47	123	25
85046165	46	52	165	25
85054170	54	58	170	25
85058170	58	62	170	25
85065240	65	70	240	25
85075200	75	80	200	25
85075250	75	80	250	25
85092250	92	95	250	25

