

Centrífugas analógicas modelos 2650 y 2655

Analogical centrifuges models 2650 and 2655

Marca|Trademark: Nahita

Centrífugas analógicas de acero pintado al horno con tapa de metacrilato abatible y dispositivo de seguridad al abrir la tapa. Están equipadas con interruptor de encendido y testigo luminoso, regulación de velocidad por posiciones y temporizador regulable hasta 30 minutos para desconexión automática y conexión continua. Con 12 portatubos de plástico para tubos de 16x100 mm (15 mL).

El modelo 2655 a diferencia del 2650 tiene capacidad para 4 tubos de 50 mL tipo Falcon sin faldón. También incorpora 4 adaptadores para poder trabajar con tubos de 16x100 mm (15 mL).

Analogical centrifuges made of furnace-painted steel with recycling methacrylate lid and safety disconnection device when the lid is opened. They included on/off switch and light device, speed adjustment by steps and adjustable timer up to 30 minutes for automatic disconnection and continuous connection. Supplied with 12 plastic tube holders for tubes of 16x100 mm (15 mL). Model 2655 differs from 2650 in capacity, as the first can accommodate 4 Falcon tubes without self standing base. Apart from this, it has 4 adapters for working with tubes of 16x100 (15 mL).



Referencia / Code	52650000	52655000
Capacidad / Capacity	12x15 mL	4X50 mL/4x15 mL
FCR máx. / Max. RCF	2147 g	2147 g
Velocidad máx. / Max. speed	4000 rpm	4000 rpm
Temporizador / Timer	0-30 min	
Dimensiones / Dimensions	28x33x27 cm	28x33x27 cm
Peso / Weight	13 Kg	13. Kg