

Modelos 2690 y 2690/5 | Models 2690 and 2690/5

Marca|Trademark:Nahita



Centrífugas analógicas de diseño compacto y funcional con un nivel muy bajo de vibración y ruido gracias a la suspensión elástica. Dispone de sistemas de desconexión del motor al abrir la tapa, temporizador de 30 min regulables, regulación de la velocidad por posiciones, piloto de señalización de encendido y tapa de metacrilato transparente azul y abatible. Con botón de puesta en marcha y ocho portatubos de aluminio para tubos de 16 x 100 (15 mL).

El modelo 2690/5 se diferencia del 2690 en la capacidad, ya que permite la centrifugación de hasta 16 tubos de 12x75 mm (5 mL)

Analogical, compact and functional centrifuges that have adapted a new vibration absorber contributing to a noise and vibration reduction. They include a safety disconnection device to turn off the motor in case the lid is opened, adjustable timer up to 30 minutes, speed adjustable knob, pilot light and blue transparent methacrylate reclining lid. With power switch and eight alluminium tubes of 16x100 mm (15 mL) holders.

Model 2690/5 differs from 2690 in capacity, as it can accommodate up to 16 tubes of 12x75 mm (5 mL).

Referencia / Code	52690000	52690050
Capacidad / Capacity	8x15 mL	16x5 mL
Máxima FCR / FCR max.	3130 g	
Temporizador / Timer	0/30 min	
Velocidad / Speed	5000 rpm	
Potencia / Power	150 W	
Peso / Weight	15 Kg	
Dimensiones / Dimensions	350x310x270 mm	300x350x370 mm