

Centrífugas modelos 2698 y 2698/5
Centrifuges models 2698 and 2698/5


Marca|Trademark:Nahita

Centrífuga de diseño compacto y funcional con un nivel muy bajo de vibración y ruido gracias a la suspensión elástica:

- Cabezal angular de aluminio
- Regulación digital de la velocidad
- Funcionamiento durante periodos de tiempo concretos hasta un máximo de 99 min.
- Interruptor de seguridad que impide que la centrífuga entre en funcionamiento si la tapa está abierta y que para automáticamente la centrifugación en caso de apertura de la tapa.
- Motor de inducción libre de mantenimiento

Digital, compact and functional centrifuges that have adopted a new vibration absorber contributing to a noise and vibration reduction:

- Alluminium fixed- angle rotor
- Digital speed adjustment
- Operation during time periods up to 99 min.
- Safety switch that avoids centrifugation if the lid is open and that automatically stops operation in case of opening the lid.
- Maintenance-free induction motor



Referencia / Code	52698000	52698005
Capacidad / Capacity	8x5 mL	16x5 mL
FCR máx. / Max. RCF	2000 g	
Velocidad máx. / Max. speed	4000 rpm	
Temporizador / Timer	0-99 min	
Dimensiones / Dimensions	32x32x28 cm	
Peso / Weight	15 Kg	