

Balanzas 5151 | Balance 5151

Marca|Trademark:Nahita

Presentamos el nuevo modelo de Balanza Digital de Precisión Nahita con una capacidad de 600g y una precisión de 0.1g. Destacamos sus principales características:

- Pantalla retroiluminada con amplios dígitos para lecturas fáciles y cómodas
- Funciones de tara, recuento de piezas y calibración
- Indicador de estabilizado (o)
- Selección de unidades de pesada: gramos (g), granos (GN), quilates (ct), pennyweights (dwt), libras (lb), onzas (oz) y onza americana (ozt)
- Burbuja de nivel y patas roscadas para perfecta equilibración de balanza
- Suministrada con adaptador de corriente y posibilidad de uso mediante pilas AA (no incluidas)

We incorporate to our catalogue this new precision digital balance with a capacity of 600 g and a readability of 0.1 g. Its main features are :

- Back-light display with large digits for and easy visualization of results
- Tare, piece counting and calibration functions
- Stabilization indicator (o)
- Weighing units: grams (g), grains (GN), carats (ct), pennyweights (dwt), pounds (lb), ounces (oz) and troy ounces (ozt)
- Bubble level and screwed legs for a perfect levelling of the equipment
- Supplied with AC adapter and possibility of use with AA batteries (not included)



Referencia / Code	55151260
Capacidad / Capacity	600 g
Sensibilidad / Readability	0.1 g
Rango de Tara / Tare range	0-600 g
Diámetro / Diameter	116 mm
Dimensiones (aprox.) / Dimensions (approx.)	245x175x80 mm

*Nota: el fabricante se reserva el derecho de modificar las características y fotografías de los productos.

*Note: Manufacturer reserves the right to modify features and photographs of the goods.