

Serie 5041 | Series 5041

Marca|Trademark:Nahita

Balanza electrónica básica con un diseño funcional muy ligero y compacto que permite un uso práctico, rápido y sencillo en aquellas aplicaciones que no requieren una elevada precisión. Entre sus características cabe destacar:

- Muy económicas ideales para su uso en el hogar y en la enseñanza
- Fabricadas en ABS, con pantalla digital de fácil lectura y teclado hermético a prueba de salpicaduras de agua
- Fácil manejo mediante tres botones
- Pesada en gramos (g) y onzas (oz)
- Tara sustractiva en todo el rango de pesada
- Calibración automática con pesa de ajuste externo (no incluida)
- Mecanismo de desconexión automática tras 2 minutos de no uso.
- Alimentación mediante adaptador externo (incluido) o mediante pilas (no incluidas)

Basic electronic balances with a functional, compact and light design that allows a practical, quick and simple use in those applications where a high accuracy is not required.

Among their main features it is worth noting:

- Very economic prices, ideal to be used at home or in educational purposes
- Made of ABS, with easy-to-read display and hermetically sealed keyboard so as to avoid liquid splashes
- Easy-to-use by only three buttons
- Weighing in grams (g) and ounces (oz)
- Subtractive tare in the whole range
- Calibration by external weight (not included)
- Auto-off mechanism after 2 minutes of no use
- Operation with adapter (included) or with batteries (not included)



Referencia / Code	55041220	55041250	55041320	55041350
Capacidad / Capacity	200 g	500 g	2000 g	5000 g
Sensibilidad / Readability	0.1 g	0.1 g	1 g	2 g
Rango tara / Tare range	0-100 g	0-500 g	0-1000 g	0-2500 g
Ø plato / Pan Ø	148 mm			
Alimentación / Power	Adaptador CA 220 V / 4 pilas 1.5 V tipo AA / AC adapter 220 V / 4 batteries 1.5 V type AA			
Dimensiones / Dimensions	170x250x45 mm			