

Serie 5023 | Series 5023

Marca | Trademark: Nahita



Las balanzas de precisión Nahita serie 5023 se caracterizan por presentar un sensor electromagnético que proporciona una gran fiabilidad y exactitud. Con una precisión de 1 mg y una capacidad máxima de 500 g están especialmente indicadas para aquellas aplicaciones que requieran la máxima precisión en las medidas.

- Pantalla de fácil lectura
- Plato de pesada fabricado en acero inoxidable
- Urna de cristal con puertas laterales y superior deslizantes
- Plástico protector permanente para proteger el equipo frente a salpicaduras o derrames de líquido
- Patas roscadas y burbuja de nivel para perfecta nivelación de la balanza
- Sensor electromagnético de gran precisión y fiabilidad
- Indicador de estabilización en pantalla
- Unidades de pesada: gramos (g), quilates (ct) y onzas (oz)
- Gancho en la parte inferior para llevar a cabo pesadas por debajo
- Función de tara sustractiva en todo el rango de pesada
- Función de recuento de piezas con número de referencia seleccionable
- Calibración mediante pesa externa; pesa de calibración incluida
- Salida RS232 para transmisión de datos
- Selección del modo de pesada para adecuarse a las necesidades y ambiente de trabajo

Nahita precision balances series 5023 are characterized by presenting an electromagnetic sensor that provides high reliability and accuracy. With a readability of 1 mg and a maximum capacity of 500 g they are especially suited for those applications where a high precision in measurements is required.

- Easy-to-read display
- Weighing pan made of stainless steel
- Glass chamber with sliding lateral and upper doors
- Permanent protection plastic cover to protect the equipment against liquid splashes or spillages
- Threaded legs and bubble level for a proper leveling of the balance
- High precision and reliability electromagnetic sensor
- Stabilization indicator on display
- Weighin units: grams (g), carats (ct) and ounces (oz)
- Provided with a hook at the bottom to perform lower weighings
- Tare function in the whole weighing range
- Piece counting function with different reference quantities available
- Calibration by external weight; calibration weight included
- RS232 output for data transmission
- Selection of the most suitable weighing mode to adapt to necessities and environmental working conditions

*Nota: el fabricante se reserva el derecho de modificar las características y fotografías de los productos.

*Note: Manufacturer reserves the right to modify features and photographs of the goods.

Serie 5023 | [Series 5023](#)

Marca|[Trademark](#):Nahita

Referencia / Code	55023210	55023220	55023230	55023250
Capacidad / Capacity	100 g	200 g	300 g	500 g
Sensibilidad / Readability	0.001 g			
Repetibilidad / Repeatability	± 0.005 g (0-'3d m-'3d50) ±0.010 (50-'3dm-'3d200) ±0.015 (m '3'3d200)			
Unidades pesada / Weighing units	g, ct y oz			
Calibración / Calibration	Por pesa externa (incluida) / With external weight (included)			
Diámetro de plato / Pan diameter	80 mm	110 mm	110 mm	110 mm
Salida RS232 / RS232 output	Sí			
Alimentación / Power	220 V, 50 ± 1 Hz			
Dimensiones cámara (LxAxH) Chamber dimensions (LxWxH)	21x19x20 cm (aprox.)			
Dimensiones totales (LxAxH) Total dimensions (LxWxH)	32x20x29 cm (aprox.)			
Peso / Weight	7.5 Kg (aprox.)			

Accesorios | [Accessories](#)

Marca|[Trademark](#):Nahita

Bandejas de pesada

Prácticas bandejas de poliestireno, desechables, especialmente diseñadas para su uso en los procesos de pesada por su ligereza a la vez que resistencia, permitiendo un cómodo y práctico uso. Se presentan en bolsas de 500 unidades.

Weighing trays

Useful and practical weighing trays made of polystyrene; disposable, light and resistant, they are specially aimed to be used during weighing tasks. Available in bags of 500 units



Ref. / Code: 40350080



Ref. / Code: 40350045

Referencia Code	Dimensiones Dimensions	Forma Shape
40350045	4.5X4.5 cm	Romboidal / Rhombus
40350080	8X8 cm	Cuadrada / Square

*Nota: el fabricante se reserva el derecho de modificar las características y fotografías de los productos.

*Note: Manufacturer reserves the right to modify features and photographs of the goods.