

Balanzas de precisión serie 5153 | Precision balances serie 5153

Marca|Trademark:Nahita



Presentamos nuestra nueva serie de balanzas Nahita con una precisión de 1 mg y disponibles en distintas capacidades a un precio realmente económico. De construcción robusta y gran fiabilidad son especialmente adecuadas para aquellas determinaciones que requieran una gran precisión en las medidas.

- Carcasa fabricada en ABS y teclado hermético a prueba de salpicaduras
- Pantalla digital retroiluminada de fácil lectura
- Urna de cristal con puertas laterales y superior deslizantes
- Patas roscadas y burbuja de nivel para una correcta nivelación del equipo
- Sistema de protección del sensor para evitar que sufra daños durante el transporte del equipo
- Sensor de bandas extensiométricas
- Pesada en gramos (g), libras (lb), quilates (ct), onzas (oz), onzas americanas (ozt), granos (GN) y pennyweights (dwt)
- Función de tara en todo el rango de pesada
- Calibración mediante pesa externa; pesa de calibración incluida
- Función de recuento de piezas, con número de referencia seleccionable
- Salida RS232 para transmisión e impresión de datos
- Alimentación mediante adaptador externo (incluido) o mediante pilas (no incluidas)

We introduce our new series of Nahita balances with a readability of 1 mg and available in different capacities, at a really economic price. With a sturdy design and a high reliability they are specially suited to those processes that require a high precision in measurements.

- ABS plastic structure and hermetically sealed control panel to avoid liquid entrance
- Easy-to-read retro-illuminated digital display
- Glass chamber with sliding lateral and upper doors
- Threaded legs and bubble level for a proper leveling of the balance
- Protection system to avoid the sensor from suffering during transport
- Strain gauge type sensor
- Weighing in grams (g), pounds (lb), carats (ct), ounces (oz), troy ounces (ozt), grains (GN) and pennyweights (dwt)
- Tare function in the whole weighing range
- Calibration by external weight; calibration weight included
- Piece counting function, with different reference quantities available
- RS232 output for data transmission and printing
- Operation with AC adapter (included) or with batteries (not included)

\*Nota: el fabricante se reserva el derecho de modificar las características y fotografías de los productos.

\*Note: Manufacturer reserves the right to modify features and photographs of the goods.

Balanzas de precisión serie 5153 | [Precision balances serie 5153](#)

Referencia / Code	55153210	5153220	55153230
Capacidad / Capacity	100 g	200 g	300 g
Sensibilidad / Readability	0.001 g		
Repetibilidad / Repeatability	+/- 0.002 g		
Linealidad / Linearity	+/- 0.002 g		
Unidades pesada / Weighing units	g, lb, oz, ozt, GN y dwt		
Calibración / Calibration	Por pesa externa / With external weight		
Diámetro de plato / Pan diameter	80 mm		
Salida RS232 / RS232 output	Sí		
Alimentación / Power	Adaptador CA 220 V ó 8 pilas tipo AA / AC adapter 220 V or 8 batteries type AA		
Dimensiones cámara / Chamber dimensions	180x180x175 mm		
Dimensiones equipo / Total dimensions	270x190x265 mm		
Peso / Weight	2.6 Kg (aprox. / approx)		

Accesorios | [Accessories](#)

Marca | [Trademark](#): Nahita

Bandejas de pesada

Prácticas bandejas de poliestireno, desechables, especialmente diseñadas para su uso en los procesos de pesada por su ligereza a la vez que resistencia, permitiendo un cómodo y práctico uso.

Se presentan en bolsas de 500 unidades.

Weighing trays

Useful and practical weighing trays made of polystyrene; disposable, light and resistant, they are specially aimed to be used during weighing tasks.

Available in bags of 500 units



Ref. / Code: 40350080



Ref. / Code: 40350045

Referencia Code	Dimensiones Dimensions	Forma Shape
40350045	4.5X4.5 cm	Romboidal / Rhombus
40350080	8X8 cm	Cuadrada / Square

\*Nota: el fabricante se reserva el derecho de modificar las características y fotografías de los productos.

\*Note: Manufacturer reserves the right to modify features and photographs of the goods.