

Viscosímetros digitales serie 807 | Digital viscosimeters serie 807

Marca|Trademark: Nahita

Los viscosímetros rotacionales de regulación digital serie 807 de Nahita son instrumentos ideales para la determinación de la viscosidad media de diversos tipos de fluidos como aceites comestibles, lubricantes, pinturas, barnices y productos similares. Gracias a su diseño, son equipos estables, eficaces y cómodos de usar.

- Equipados con 4 husillos que permiten al usuario seleccionar las condiciones más adecuadas según la viscosidad del líquido a analizar
- Su panel de mandos nos permite seleccionar cómodamente los parámetros de funcionamiento (velocidad y husillo).
- Disponen de una sonda que sumergida en el líquido nos indica en pantalla su temperatura.
- La base presenta 2 patas roscadas para conseguir una perfecta nivelación del equipo mediante el nivel a burbuja situado en la parte superior.
- La columna, maciza con sistema de piñón y cremallera, y la pinza de sujeción del cabezal, con tornillo y guía, permiten colocar cómodamente el cabezal a la altura deseada asegurando una buena estabilidad del mismo.
- Suministrados en un robusto maletín metálico para mantenerlos siempre bien recogidos y en perfectas condiciones de uso.

Nahita rotary viscosimeters with digital adjustment series 807 are ideal to determine the viscosity of several type of fluids such as edible oil, lubricants, paints, varnishes and similar products. Thanks to their design, they are stable, efficient and easy-to-use equipments.

- Equipped with 4 rotors to allow user to select the most adequate conditions according to sample viscosity.
- The front panel allows selecting easily the operation parameters (speed and rotor).
- Provided with a probe that once is introduced in the liquid shows its temperature on display.
- The base presents 2 screwed legs to get a perfect levelling of the equipment by means of a bubble level placed at the upper side of the equipment.
- The column, solid and with pinion and rack system, and the head holding platform allows placing the head at the desired height assuring a good stability.
- Supplied in a sturdy metallic case to keep the equipments always in perfect conditions.



Referencia / Code	50807001	50807002
Rango medida / Measuring range	1-6x10 ⁵ mPa·s (o cP)	1-6x10 ⁶ mPa·s (cP)
Velocidad / Speed	1-60 rpm	0.1-99.9 rpm
Precisión / Accuracy	±2.0%	
Repetibilidad / Repeatability	0.5%	
Husillos / Rotors	nº 1, 2, 3, 4	