

Nahita le ofrece esta serie de morteros fabricados en porcelana e ideales para triturar, machacar y mezclar sustancias de poca dureza. De base ancha y estable y con pico para un cómodo vertido, presentan el exterior esmaltado y el interior y la base sin esmaltar. Cada mortero se suministra con su mano fabricada también en porcelana y con mango esmaltado y la superficie de molienda sin esmaltar. Las manos para morteros también están disponibles individualmente.

Nahita offers this series of mortars made of porcelain and ideal for grinding, crushing and mixing substances of low hardness. With a wide and stable base and a spout for an easy pouring, they are glazed outside except for the base and unglazed inside. Each mortar is supplied together with its porcelain pestle with the handle glazed and the grinding surface unglazed. Mortar pestles are also available individually.

Morteros de porcelana con mano | [Porcelain mortars with pestle](#)

Marca|Trademark: Nahita



Referencia Code	Capacidad Capacity	Diámetro Diameter	Altura Height	Long. Mano Pestle length	Ø Bola Ball Ø
20211060	25 mL	60 mm	35 mm	65 mm	19 mm
20211080	50 mL	80 mm	45 mm	90 mm	26 mm
20211100	100 mL	100 mm	45 mm	110 mm	27 mm
20211130	300 mL	130 mm	62 mm	125 mm	36 mm
20211160	500 mL	160 mm	74 mm	145 mm	40 mm
20211210	1000 mL	216 mm	90 mm	215 mm	46 mm
20211250	2000 mL	254 mm	110 mm	245 mm	53 mm

Manos para morteros | [Pestles for porcelain mortars](#)

Marca|Trademark: Nahita



Referencia Code	Longitud Length	Ø Mango Handle Ø	Ø Bola Ball Ø	Para mortero For mortar
20210060	65 mm	11 mm	19 mm	60 mm
20211080	90 mm	15 mm	26 mm	80 mm
20210100	110 mm	15 mm	27 mm	100 mm
20210130	125 mm	18 mm	36 mm	130 mm
20210160	145 mm	20 mm	40 mm	160 mm