

Evaporador rotativo modelo 500 | Rotary evaporator model 500

Marca|Trademark:Nahita

Este completo evaporador incluye baño calefactor en acero inoxidable con termostato y regulador digital de temperatura hasta 100 °C. Dispones de lectura digital de la temperatura del baño por sonda y dispositivo de vaciado, así como de un soporte elevador accionado por electromotor que ejecuta la subida (hasta 150 mm) y bajada automáticamente, simplemente pulsando un botón.

Con regulación electrónica de la velocidad de giro del motor, incluye juego de vidrio con matraz de 250 mL, aceptando matraces de hasta 1000 mL (esmerilado 29/32). El juego de vidrio gira 45° e incluye: refrigerante, pieza de unión de vidrio con 3 salidas, pipeta de entrada de vidrio y adaptador de unión esmerilada, además de los matraces de fondo redondo.

This complete evaporator includes a heating bath made of stainless steel provided with a thermostat and digital temperature adjustment up to 100° C. It can read temperature digitally by means of a probe and emptying device, as well as an electromotor-driven lifter that produces the up and down automatically, simply by pressing a button.

Among its features we stand out its electronic adjustment of the motor 's rotaring speed. It also includes a glass set with a flask of 250 mL, but it admits sizes up to 1000 mL (29/32 ground joint). The glass set can rotate 45° and is composed of: condenser , joint piece piece for pipette, pipette for adding and adapter wit tapered joint grinding, as well as flasks with round bottom.



Referencia / Code	59000500
Velocidad / Speed	10-90 rpm
Esmerilado / Ground joint	29/32
Temperatura / Temperature	100 °C