

Agitador rotativo multiposición

Rotator

9




Datos técnicos sin accesorios Technical datas without accessory	9
Potencia Power [W]	40
Movimiento de agitación Shaking movement	Rotatorio Rotating
Ajuste de la varilla rotatoria Range of rotating rod	0 - 90°
Velocidad Speed range [r.p.m.]	10 - 70
Pantalla velocidad Speed display	LCD
Temporizador Timer [mn]	1 - 1199
Tipo de funcionamiento Run type	Con tiempo / en continuo time / continuous
Dispositivo de frenado Braking device	Si Yes
Dimensiones Dimension (L x l x h) [mm]	510 x 190 x 260
Peso Weight [Kg]	8
Clase de protección según Protection class acc. to	IP 21 DIN EN60529

CARACTERÍSTICAS

- control digital de velocidad y pantalla LCD
- movimiento « end-over-end » o movimiento ondulante
- velocidad desde 0 hasta 70 r.p.m.
- regulación del ángulo de funcionamiento entre 0° y 90°
- modo de funcionamiento continuo o durante periodos programados de tiempo

FEATURES

- digital control speed and LCD display timer and speed
- end-over-end or rolling action, rotator gives gentle but effective mixing
- adjustable speed between 0 and 70 r.p.m.
- adjustable mixing angle between 0 and 90°
- timing and continuous operation modes



Modelo Model	Referencia
9 con accesorios tipo « asador » - with rotisserie accessory	50683090
9 con accesorios tipo « disco » - with disc accessory	50683091