

Los vasos de precipitados son un elemento indispensable en el trabajo diario de laboratorio para la preparación de mezclas, calentamiento de muestras o el trasvase de volúmenes poco precisos. Le ofrecemos una amplia gama de vasos fabricados en vidrio borosilicato Nahita, que cubren un extenso rango de volúmenes para adaptarse a las necesidades de cualquier usuario. Los vasos, de forma alta o baja, presentan pico y se encuentran graduados (excepto los vasos de 25 ml) para una mayor comodidad del usuario. Todos ellos cumplen la norma ISO 3819.

Beakers are an essential element in daily lab work for mixture preparation, sample heating and measuring not precise volumes. We offer a wide range of Nahita borosilicate glass beakers that cover a broad range of volumes to adapt to any user necessities. Tall and low form beakers present spout and are graduated (except 25 ml beakers) according to ISO 3819.

Vasos forma alta | Tall form beakers

Marca|Trademark: Nahita



Referencia Code	Capacidad Capacity	Diámetro Diameter	Altura Height
10153121	25 mL	30 mm	60 mm
10153151	50 mL	38 mm	70 mm
10153201	100 mL	48 mm	80 mm
10153211	150 mL	54 mm	95 mm
10153221	250 mL	60 mm	120 mm
10153241	400 mL	70 mm	130 mm
10153261	600 mL	80 mm	150 mm
10153311	1000 mL	95 mm	180 mm

Vasos forma baja | Low form beakers

Marca|Trademark: Nahita



Referencia Code	Capacidad Capacity	Diámetro Diameter	Altura Height
10155121	25 mL	35 mm	55 mm
10155151	50 mL	43 mm	60 mm
10155201	100 mL	51 mm	71 mm
10155211	150 mL	61 mm	86 mm
10155221	250 mL	70 mm	98 mm
10155241	400 mL	79 mm	115 mm
10155261	600 mL	90 mm	125 mm
10155311	1000 mL	105 mm	145 mm
10155321	2000 mL	130 mm	185 mm
10155331	3000 mL	150 mm	210 mm
10155351	5000 mL	170 mm	270 mm