

Tubos digestores | Digesting tubes

Marca|Trademark: Nahita

Los tubos digestores de Nahita están diseñados para su utilización en el proceso de determinación de nitrógeno y fósforo totales en digestores automáticos. Fabricados en vidrio borosilicato y de paredes gruesas y uniformes, son altamente resistentes a la acción de ácidos y bases fuertes y a las elevadas temperaturas requeridas en el proceso de digestión. Además, su diseño de fondo redondo y gran longitud permite una distribución uniforme del calor en la muestra y minimiza las pérdidas por salpicaduras.

Nahita digesting tubes are designed to be used in the process of total nitrogen and total phosphorus determination in automatic digestors. Made of borosilicate glass and with thick and uniform walls they are highly resistant to strong acids and bases and to high temperatures required in the process of digestion. Moreover, their round bottom and large length design allows a uniform heat distribution and avoids spilling.



Referencia Code	Diámetro Diameter	Longitud Length	Descripción Description	Digestor Digestor
13800010	42 mm	300 mm	Con estrangulamiento With constrictio	Buchi Buchi
13800020	42 mm	300 mm	Sin estrangulamiento Without constriction	Tecator y otros Tecator and others