


## DP 4000 water still 4 L/h VERTICAL

### CODE PDG007

The Nahita-Blue DP 4000 model distillation equipment is a compact distiller used for the production of distilled water in laboratories, institutes, the food industry, pharmaceuticals, and universities.

The glass parts are made of borosilicate glass to ensure better equipment function and safety.

Easy to operate and clean.



## Technical specifications

Reference	PDG007
Production capacity	4,0 L/h
Energy consumption	3 100 W
Water consumption	45 L / h
Power supply	230 V / 50 Hz
Conductivity at 20°C	1,5 - 2,5 $\mu$ S/cm
Width of equipment	190 mm
Length of equipment	190 mm
Height of equipment	660 mm