


Flat bottom flask



- 1 | Flat bottom spherical flask with frosted mouth.
- 2 | Complies with USP regulations.

Code	Cap.	Frosted	Height	Diameter
PBK001	50 mL	14/23	85 mm	51 mm
PBK002	50 mL	29/32	85 mm	51 mm
PBK018	50 mL	19/26	85 mm	51 mm
PBK019	50 mL	24/29	85 mm	51 mm
PBK003	100 mL	19/26	103 mm	64 mm
PBK004	100 mL	24/29	103 mm	50 mm
PBK005	100 mL	29/32	103 mm	64 mm
PBK006	250 mL	19/26	130 mm	85 mm
PBK007	250 mL	24/29	130 mm	85 mm
PBK008	250 mL	29/32	130 mm	85 mm
PBK009	500 mL	24/29	160 mm	103 mm
PBK010	500 mL	29/32	160 mm	160 mm
PBK011	1000 mL	24/29	187 mm	131 mm
PBK012	1000 mL	29/32	187 mm	131 mm
PBK013	2000 mL	29/32	230 mm	166 mm

