



## BOTTLES WIDE AND NARROW NECK BOTTLES

Code CHB001 - CHB008

Code CHG001 - CHG008

Endo Glassware offers a wide range of multipurpose bottles of different capacities and ideal for storing any type of solution. All bottles are available in transparent (white) glass or amber glass (for those photosensitive solutions or reagents).

With standard ground joint and polyethylene stopper.



### Technical specifications

Code	Capacity	Base Ø	Neck Ø	Glass
CHB001	100 mL	54 mm	95 mm	Transparent
CHB005	100 mL	54 mm	95 mm	Amber
CHB002	250 mL	74 mm	125 mm	Transparent
CHB006	250 mL	74 mm	125 mm	Amber
CHB003	500 mL	88 mm	165 mm	Transparent
CHB007	500 mL	88 mm	165 mm	Amber
CHB004	1000 mL	112 mm	200 mm	Transparent
CHB008	1000 mL	112 mm	200 mm	Amber



### Technical specifications

Code	Capacity	Base Ø	Neck Ø	Glass
CHG001	100 ml	54 mm	95 mm	Transparent
CHG005	100 ml	54 mm	95 mm	Amber
CHG002	250 ml	74 mm	125 mm	Transparent
CHG006	250 ml	74 mm	125 mm	Amber
CHG003	500 ml	88 mm	165 mm	Transparent
CHG007	500 ml	88 mm	165 mm	Amber
CHG004	1000 ml	112 mm	200 mm	Transparent
CHG008	1000 ml	112 mm	200 mm	Amber