

Matraces | Flasks

Marca|Trademark: Nahita

Los kitsatos son unos matraces robustos de paredes gruesas y fabricados en vidrio borosilicato de excelente resistencia física y química. De uso extendido en laboratorios de química y biología, poseen una oliva lateral para ser conectados a una bomba de vacío o trompa de agua y pueden ser utilizados con diversos tipos de embudos para la filtración rápida de soluciones. Presentan además graduación aproximada y amplia superficie de rotulación en blanco.

Filtering flasks are sturdily constructed with thick walls and are made of borosilicate glass of an excellent physical and chemical resistance. Widely used in chemistry and biology laboratories, they present a lateral hose connection to be connected with some sort of vacuum assembly and they can be used with any kind of funnel for fast solution filtration. They also feature a large marking area and approximate graduation.



Referencia	Capacidad	Ø Base	Ø Boca	Altura	Ø Oliva
Code	Capacity	Base Ø	Neck Ø	Height	Hose connect.
12418220	250 mL	85 mm	35 mm	140 mm	8-10 mm
12418250	500 mL	105 mm	38 mm	180 mm	9-11 mm
12418310	1000 mL	140 mm	44 mm	230 mm	10-12 mm