





Single rotor benchtop digital centrifuges

- 1 | Angular rotor (painted aluminum)
- 2 | Digital adjustment of parameters
- 3 | Furnace-painted steel casing
- 4 | Methacrylate lid with locking anchorage
- 5 | Suction cup legs
- 6 | Lid contact switch and electromagnetic lock
- 7 | Falcon type tubes, (16x100 mm). Using the adapter GBK003 Vacutainer tubes of 13x75 mm and 13x100 mm can be centrifuged (removing the plastic cushion for 13x100 mm)



[01]



Rotor: For Falcon type tubes, (16x100 mm). Vacutainer tubes of 13x75 mm (using adapter) and 13x100 mm (removing the plastic cushion)

[02]



Control pannel

Code	GBF011	GBF012
Model	2650-D	2650-D+
Motor	With brushes	Induction type (brushless)
Capacity	12 tubes of 15 mL	
Speed range	500-4000 rpm (increment 50 rpm)	
Max RCF	2110xg (increment 10xg)	
Timer	10 sec - 99 min 59 sec	
Screen	Backlit LCD	
Dimensions	330x 285x270 mm	
Weight	9,75 kg	
Power	220V/50Hz	