

Laboratory bottles, GL45, borosilicate glass 3.3

- Made of high-quality borosilicate glass 3.3.
- According to EN ISO standard 4796-1.
- The bottles have a retrace code, a durable volume scale as well as a GL45 thread (GL32 thread for 50ml bottles).
- The bottles as well as the closures are resistant against a wide range of chemicals.

1 Laboratory bottles, clear, with PP screw caps

With blue pouring rings and screw caps with integrated lip seal made of PP. Autoclavable at 121 °C.

Capacity	Ext. Ø	Height	Neck thread	PK	Cat. No.
ml	mm	mm	GL		
50	46	88	32	10	7.690 060
100	56	100	45	10	7.690 045
250	70	138	45	10	7.690 046
500	86	176	45	10	7.690 047
1000	101	225	45	10	7.690 048
2000	136	260	45	10	7.690 049
5000	181	330	45	1	7.690 050
10000	227	410	45	1	7.690 051

Height without screw cap