



GLASS WATER STILL BASIC PH4



- 1 | The glass water still BASIC PH4 is specially design for an easy and continuous distillation process with total safety and confidence thanks to its built-in safety thermostat.
- 2 | Provided with a 3 Kw nickel-chromium plated copper resistance, the water still BASIC PH4 has a production of 4 L/h of distilled water.
- 3 | The horizontal boiler is made of high quality borosilicate glass. CE certified.

Code	PDG001
Conductivity (with tape water)	3.0-4.0 $\mu\text{S}/\text{cm}$
Conductivity (with treated water)	1.5-2.0 $\mu\text{S}/\text{cm}$
pH	5.5-6.0
Dimensions	50x15x45 cm

Accessories

In case of damage, breakage, etc. of any of the components of the equipment; there are available the following pieces for replacement.

Code	Description
PDJ004	Glass boiler
PDJ005	Glass condenser
PDJ003	Hose kit with red conectors
PDJ002	Joint hit+screw+2 fixers for boiler
PDJ006	Metal support
PDJ007	2 fixers for boiler