



Pocket size digital pH meter



- 1 | Small, portable size for use anywhere, anytime.
- 2 | Low power consumption. High accuracy. Wide measuring range.
- 3 | It is not necessary to keep the electrode immersed in 3M KCl.
- 4 | Automatic temperature compensation.
- 5 | Buttons for:
 - ON/OFF switch
 - HOLD button
 - CAL button for calibration
- 6 | Acidic or basic pH range calibration. 3-point calibration.
- 7 | Includes carrying case, manual and 3 sachets of pH 4.01, 6.86 and 9.18 pH buffer solutions.

Code	KDB011
Display	LCD
Replaceable electrode	Code KDG046
Range	0.00 – 14.00 pH
Resolution	0.01 pH / 0.1 pH (<i>can be changed</i>) 0.1 °C / 0.1 °F
Accuracy	± 0.05 pH
Power supply	3 x 1.5V button cell batteries type LR44 (<i>included</i>)
Operating temperature	0°C~50°C (32~122°F)