

## Analogical stirrer with heating - model 690/1

Code **LBB001**

These small analogical magnetic stirrers with heating are ideal for those places with little space due to their compact design. For mixing up to 2 l of all types of solutions in which introducing a stirring magnet is possible.

- Metallic structure covered with epoxy.
- Speed and temperature regulation knob.
- On/Off lateral switch.
- Equipped with rod and clamp, supplied with a teflon covered magnet.



### Technical specifications

Code	LBB001
Speed range	100-200 rpm
Máx. temperature (plate)	180°C ± 10%
Stirrable volume	50-2000 ml
Plate dimensions	120x120 mm
Power	220V/50 Hz
Dimensions (LxWxH)	200x120x90 mm
Weight	2 kg