



Digital magnetic stirrer rslab-11

Code **LBG010**, **LBJ001**, **LBD005**, **LBD006**

- Temperature sensor Pt1000 controls (For RSLAB-11C model).
- Real-time display setting and actual temperature (For RSLAB-11C model).
- LED display speed (RSLAB-11NCD and RSLAB-11C models).
- Work plate material: plastic for 11NCA and RSLAB-11NCD model without heating and stainless steel work plate for RSLAB-11C model.
- Protection class according to DIN EN60529: IP42.

Without heating



Code	LBG010*	LBJ001*
Model	RSLAB-11NCA, analogical	RSLAB-11NCD, digital
Display	-	LED
Speed range	0-1500 rpm	200-1500 rpm
Max. volume	3 L	
Plate diameter	135 mm	
Plate material	Plastic	
Dimensions	260 x 150 x 80 mm	
Weight	1,8 Kg	

*Delivered without stirrer bar.

With heating



Code	LBD005	LBD006
Model	RSLAB-11C, digital	
Display	LED	
Speed range	200-1500 rpm	
Temperature range	RT-280°C	
Temperature display accuracy	1 °C	
Max. volume	3 L	
Plate diameter	135 mm	
Plate material	SS plate with ceramic-coated	Stainless steel
Dimensions	260 x 150 x 80 mm	
Weight	1,8 Kg	



Accessories for RSLAB-11C

COMPATIBLE ACCESSORIES NOT INCLUDED WITH THE EQUIPMENT



LBM002
Temperature sensor Pt1000



LBM005
Temperature sensor Pt1000
with glass coated



LBM004
Stand kit for sensor Pt1000
for RSLAB-11C



LBM006
Adapter for round bottom
flask 50 mL




LBM007
Adapter for round bottom
flask 100 mL



LBM013
Black poly-bloc with 4 holes and
diameter 28 mm, width 43 mm,
for 40 mL flask



LBM014
Green poly-bloc with 6 holes and
diameter 17,8 mm, width 26 mm,
for 8 mL flask



LBM015
Gold poly-bloc with 4 holes and diameter
21,6 mm, width 31,7 mm, for 16 mL flask



LBM008
Adapter for round bottom
flask 250 mL



LBM009
Adapter for round bottom
flask 500 mL



LBM010
Red poly-bloc with 11 holes and
diameter 12,5 mm, width 20 mm,
for 4 mL flask



LBM011
Purple poly-bloc with 4 holes and
diameter 28 mm, width 24 mm,
for 20 mL flask



LBM012
Blue poly-bloc with 4 holes and
diameter 28 mm, width 30 mm,
for 30 mL flask



LBM016
Blue support for 4 poly-bloc stands



LBM017
Blue ring for 4 poly-bloc stands

Check our Magnetic stirrer bar section in Catalogue