



## Magnetic stirrers RSLAB-2

Code LBD002, LBD008

- Temperature sensor, safety temperature can be set.
- Temperature sensor Pt1000 included (Stand kit not included).
- Warning display indicates presence of any residual heat when unit is switched off.
- Connection according to DIN 12878 for the use of contact thermometers.
- With watertight ring to protect the body of equipment.

### With heating RSLAB-2C



Code	LBD002	LBD008
Display	LCD	
Speed range	100-1500 rpm	
Maximum stirring volume	20 L	
Temperature range	RT - 340°C	
Safety range	360°C	
Plate material	Ceramic	Stainless steel
Plate dimensions	ø 135 mm	
Approx. dimensions	280 x 160 x 85 mm	
Approx. weight	2,8 Kg	



Accessories (not included)



Temperature probe PT1000  
Code LBM002



Temperature probe PT1000  
Code LBM005  
With glass coating



Support kit for probe PT1000  
Code LBM003 and LBM004



Adapter for round bottom  
Code LBM006  
Flasks 50 mL



Adapter for round bottom  
Code LBM007  
Flasks 100 mL



Poly-block black stand  
Code LBM013  
4 holes for 40 mL tubes



Poly-block green stand  
Code LBM014  
6 holes for 80 mL tubes



Poly block yellow stand  
Code LBM015  
4 holes with 16 mL tubes



Adapter for round bottom  
Code LBM008  
flasks 250 mL



Adapter for round bottom  
Code LBM009  
flasks 500 mL



Poly-block red stand  
Code LBM010  
11 holes for 4 mL tubes



Poly-block purple stand  
Code LBM011  
4 holes for 20 mL tubes



Poly-block blue stand  
Code LBM012  
4 holes for 30 mL tubes



Blue support  
Code LBM016  
4 poly-block stands



Blue ring  
Code LBM017  
4 poly-block stands



Protective cover  
Code LBM018  
Resistant up to 135°