



Model 700 Mini Vortex Mixer

Code LGB006

The mixer is manufactured using the design concept of compression type circular eccentric rotary oscillation to generate vortices. Its appearance is simple, with a low center of gravity structure that is compact and firm, and easy to operate. Its constant speed ensures rapid mixing of samples.

- Zinc alloy integrated base, low center of gravity balance, stable operation without displacement during use.
- External battery packs can be connected to meet on-site inspection needs outside of the laboratory.



Technical specifications

Reference	LGB006
Motor	DC12V / 1 A
Power (W)	10
Maximum load (kg)	0,6
Speed	4000
Oscillation mode	Circumference
Circumferential diameter	4,5 mm
Dimensions	95 x 95 x 80 mm
Weight (kg)	0,9
Noise level	≤ 40 dB (A)
For tubes	Up to Ø 30 mm
Protection level	IP20