



Manual pipette controllers

Code MBJ007 - MBJ008 - MBJ009

- Pipette pumps provide an effective and affordable solution for a comfortable and safe method of pipetting liquids of any nature.
- Its ergonomic design avoids awkward pipetting motions and allows drawing and dispensing liquids without changing hand position and with the only movement of the thumb.
- The thumbwheel can be rotated for precision either aspirating or dispensing, while the plunger may be pressed for quick emptying. Thanks to its flexible rubber chuck they can be used with several capacity pipettes made of glass or plastic.
- Pipette pumps present different colours depending on the maximum volume they can aspirate, are resistant to acids and alkalies, and easily disassembled for cleaning.



Technical specifications

Code	MBJ007	MBJ008	MBJ009
Volume	2 mL	10 mL	25 mL
Colour	Blue	Green	Red