

GRADUATED MEASURING CYLINDERS

Graduated measuring cylinders are an essential instrument in the daily lab work for the preparation of solutions or the measurement of any type of liquid.

Endo Glassware presents a complete range of tall form graduated measuring cylinders class A and B made of borosilicate glass; they are available with spout or with standard PP stopper.

All of them present an hexagonal base for a higher stability and have been calibrated at 20°C to contain (In); they are made according to ISO 4788.

Graduated measuring cylinders, class B (white graduation)



Reference	Capacity	Diameter	Height	Graduation	Tolerance
11636050	5 mL	11 mm	112 mm	0.1 ml	± 0.1 ml
11636110	10 mL	12 mm	135 mm	0.2 ml	± 0.2 ml
11636120	25 mL	20 mm	160 mm	0.5 ml	± 0.5 ml
11636150	50 mL	25 mm	195 mm	1.0 ml	± 1.0 ml
11636210	100 mL	30 mm	240 mm	1.0 ml	± 1.0 ml
11636220	250 mL	40 mm	320 mm	2.0 ml	± 2.0 ml
11636250	500 mL	50 mm	380 mm	5.0 ml	± 5.0 ml
11636310	1000 mL	65 mm	420 mm	10.0 ml	± 10.0 ml
11636320	2000 mL	80 mm	520 mm	20.0 ml	± 20.0 ml

Graduated measuring cylinders, class A (blue graduation)



Reference	Capacity	Diameter	Height	Graduation	Tolerance
11637050	5 mL	11 mm	112 mm	0.1 ml	± 0.05 ml
11637110	10 mL	12 mm	135 mm	0.2 ml	± 0.10 ml
11637120	25 mL	20 mm	160 mm	0.5 ml	± 0.25 ml
11637150	50 mL	25 mm	195 mm	1.0 ml	± 0.50 ml
11637210	100 mL	30 mm	240 mm	1.0 ml	± 0.50 ml
11637220	250 mL	40 mm	320 mm	2.0 ml	± 1.00 ml
11637250	500 mL	50 mm	380 mm	5.0 ml	± 2.50 ml
11637310	1000 mL	65 mm	420 mm	10.0 ml	± 5.00 ml
11637320	2000 mL	80 mm	520 mm	20.0 ml	± 10.00 ml

Graduated measuring cylinders class B with PP stopper (white graduation)



Reference	Capacity	Diameter	Height	Grad.	Tolerance	Mouth
11638050	5 mL	11 mm	112 mm	0.1 ml	± 0.1 ml	10/19
11638110	10 mL	12 mm	135 mm	0.2 ml	± 0.2 ml	10/19
11638120	25 mL	20 mm	160 mm	0.5 ml	± 0.5 ml	19/26
11638150	50 mL	25 mm	195 mm	1.0 ml	± 1.0 ml	19/26
11638210	100 mL	30 mm	240 mm	1.0 ml	± 1.0 ml	29/32
11638220	250 mL	40 mm	320 mm	2.0 ml	± 2.0 ml	29/32
11638250	500 mL	50 mm	380 mm	5.0 ml	± 5.0 ml	34/35
11638310	1000 mL	65 mm	420 mm	10.0 ml	± 10.0 ml	45/40
11638320	2000 mL	80 mm	520 mm	20.0 ml	± 20.0 ml	45/40

Graduated measuring cylinders class A with PP stopper (blue graduation)



Reference	Capacity	Diameter	Height	Grad.	Tolerance	Mouth
11639050	5 mL	11 mm	112 mm	0.1 ml	± 0.05 ml	10/19
11639110	10 mL	12 mm	135 mm	0.2 ml	± 0.10 ml	10/19
11639120	25 mL	20 mm	160 mm	0.5 ml	± 0.25 ml	19/26
11639150	50 mL	25 mm	195 mm	1.0 ml	± 0.50 ml	19/26
11639210	100 mL	30 mm	240 mm	1.0 ml	± 0.50 ml	29/32
11639220	250 mL	40 mm	320 mm	2.0 ml	± 1.00 ml	29/32
11639250	500 mL	50 mm	380 mm	5.0 ml	± 2.50 ml	34/35
11639310	1000 mL	65 mm	420 mm	10.0 ml	± 5.00 ml	45/40
11639320	2000 mL	80 mm	520 mm	20.0 ml	± 10.00 ml	45/40