



DIGIPETTE Junior fixed volume micropipettes

Code MGD095 - MGD101

Compact and light pipettes specially designed for small laboratories and school laboratories, with excellent performance and ergonomics.

The Junior line is composed of 7 references, covering the volume range from 5 μ L to 200 μ L.

Features:

- Compact working size of 130mm ensures fatigue free operation
- Two step plunger allows for "reverse pipetting" technique
- Calibrated as per EN ISO 8655-6 standard
- Can be sterilised in autoclave
- Stainless steel piston for reliable performance and years of durability
- Coloured push buttons for easy identification of volumes



Code	Nominal volume (μ L)	Test volume (μ L)	Inaccuracy ($\pm$)%	Imprecision ($\pm$)%
MGD095	5	5	1,3	1,2
MGD096	10	10	0,8	0,8
MGD097	20	20	0,6	0,5
MGD098	25	25	0,5	0,3
MGD099	50	50	0,5	0,3
MGD100	100	100	0,5	0,3
MGD101	200	200	0,4	0,2



Fixed volume
5 μ L



Fixed volume
10 μ L



Fixed volume
20 μ L



Fixed volume
25 μ L



Fixed volume
50 μ L



Fixed volume
100 μ L



Fixed volume
200 μ L