




**mini magnetic stirrer**  
**NAHITA-BLUE**  
models 680



The new range of Nahita-Blue presents an up-to-date and sturdy design and they make up a series of useful instruments for general applications in the laboratory.

	50692010	50692000	50682000	50680000
<b>Control</b>				
Digital	•			
Analogical		•	•	•
<b>Heating</b>				
With heating	•	•		
Without heating			•	•
<b>Temperature control</b>				
Solution	•			
Plate		•	-	-



**MINI MAGNETIC  
stirrer  
NAHITA-BLUE**

The new mini magnetic stirrer Nahita-Blue is a basic laboratory instrument that presents a very dynamic design and an economic price. It is characterized for being compact and lightweight, which makes this stirrer an essential instrument for laboratories with limited space. It also comes with a vehicle power adapter, allowing it to be ideal for field work, where its use is required but without having a nearby power outlet.



mini magnetic stirrer NAHITA-BLUE

## mini magnetic stirrer NAHITA-BLUE



## technical features

Code	50680000
Modelo	680
Speed range	100-2000rpm
Max. volume	1L
Plate dimensions	110mm
Power	AC adapter 220V/12V
Dimensions (LxWxH)	120x120x50mm
Weight	0.4Kg

## mini magnetic stirrer NAHITA-BLUE equipment parts

### mini magnetic stirrer NAHITA-BLUE features

It is a robust and durable product, as it is made of ABS plastic, allowing resisting minor splashing of almost any chemical. Simple to use, it has an analog control that allows agitation speed from 100 to 2000rpm

### equipment parts back side



## mini magnetic stirrer NAHITA-BLUE accessories

It is supplied with :

- support rod to use with probes for temperature, pH, ...
- stirring magnet.
- vehicle power adapter.

Creating a fully functional and complete product for almost any need that the user may have.



## AUXILAB, S.L.

Polígono Morea Norte Calle D Nº 6  
31191 Beriain - Navarra - España  
T: +34 948 310 513  
[comercial@auxilab.es](mailto:comercial@auxilab.es)  
[www.auxilab.es](http://www.auxilab.es)

## ROGO-SAMPAIC

23 Rue Ampère, ZI de Vilemilan  
91325 Wissous - Paris - France  
T: +33 1 69 53 67 67  
[rogosampaic@rogosampaic.com](mailto:rogosampaic@rogosampaic.com)  
[www.rogosampaic.com](http://www.rogosampaic.com)



**AUXILAB S.L.**  
Material de laboratorio  
Laboratory supplies