



high speed
multirotor
centrifuge
nahita-blue

high speed mult rotor centrifuge

Model 2624/2

This centrifuge that reaches 15000 rpm is specially suitable for use in diverse fields such as biology, veterinary science, chemistry, genetic engineering and immunology. This centrifuge offers its versatility due to its multiples rotors, making it an indispensable element in the routine laboratory work.

technical specifications

Code	52624002
Rotors	Fix angle (not included)
Max. speed	15000 rpm*
RCF max.	20375g*
Timer	1-99 minutes
Power	230V, 50Hz-60Hz
Dimensions	36X28X26 cm
Weight	15 Kg

* See rotor specifications



high speed multirotor centrifuge: Model 2624/2



[01] It is specially aimed to protocols that require centrifugation at high RCF.

[02] It is an equipment with a robust and resilient configuration, which has a furnace-painted steel casing and base.

[03] It has an elastic suspension system to guarantee a soft performance, silence and without vibrations.

[04] It is provided with an electromagnetic security closure that allows a safety work with the equipment. Moreover, it has a manual opening mechanism to use in case of failure of power or of the electric opening system.

[05] The lid is equipped with suspension mechanisms to avoid accidents during the centrifuge loaded or emptied.

[06] The cabinet is equipped with an auxiliary ventilation system that helps to prevent an excess overheating and the deterioration of the centrifuged samples.

[07] It's induction motor without brushes makes this equipment free of periodic maintenance.

[08] It has digital control of parameters with a double LED display.

high speed
multirotor
centrifuge
model 2624/2
FRONT PANEL



high speed multirotor centrifuge equipment parts



1
Lid

2
Front panel

4
robust and resilient configuration, which has a furnace-painted steel casing and base.



3
Power switch on/off

2
Socket

high speed multirotor centrifuge

high speed multirotor centrifuge: **different rotors**

In order to be suitable for different user needs, the equipment has 5 easy installing fix-angle rotors, which are supplied separately. Moreover, the centrifuge has a system that relates the rpm / g depending on the number of rotor that is going to be used.

Code	Description
92624001	n° 1: tubes 1.5/2.2mL x 24 Máx. 15000rpm, 20375g
92624002	n° 2: haematocrit capillary x 24Máx. 12000rpm, 15455g
92624003	n° 3: tubes bos 5mL x 10 Máx. 13500rpm, 12920g
92624004	n° 4: PCR plate (0.2x8x4cm) Máx. 14800rpm, 16200g
92624005	n° 5: tubes 0.5mL x 36 Máx. 13500rpm, 13250g

Code	92624001
Tubes	1.5/2.2mL x 24 Máx. 15000rpm, 20375g



Code	92624002
Haematocrit capillary	x 24Máx. 12000rpm, 15455g



AUXILAB, S.L.

Polígono Morea Norte Calle D Nº 6
31191 Beriain - Navarra - España
T: +34 948 310 513
comercial@auxilab.es
www.auxilab.es

ROGO-SAMPAIC

23 Rue Ampère, ZI de Vilemilan
91325 Wissous - Paris - France
T: +33 1 69 53 67 67
rogosampaic@rogosampaic.com
www.rogosampaic.com



AUXILAB S.L.
Material de laboratorio
Laboratory supplies