

▶ With spout, PP

Kartell

- 1| Excellent chemical resistance. Embossed graduation.
- 2| Can endure temperatures up to +120 °C.

| Code | Cap. | h. | Grad. | Tol. |
|--------|---------|--------|---------|----------|
| DBP001 | 10 mL | 88 mm | 0.25 mL | ±0.1 mL |
| DBP002 | 25 mL | 107 mm | 1 mL | ±0.25 mL |
| DBP003 | 50 mL | 143 mm | 2.5 mL | ± 0.5mL |
| DBP004 | 100 mL | 178 mm | 5 mL | ±1.0 mL |
| DBP005 | 250 mL | 264 mm | 10 mL | ±2.5 mL |
| DBP006 | 500 mL | 305 mm | 10 mL | ±5 mL |
| DBP007 | 1000 mL | 332 mm | 25 mL | ±10 mL |
| DBP008 | 2000 mL | 370 mm | 50 mL | ±20 mL |

