

Filtration equipment, 300 mL

Code PJB001

Filtration equipment, 300 mL - Endo Glassware's filtration equipment is ideal for vacuum filtration of any type of solutions for particle contamination analysis and also for solvent filtration for HPLC.

Technical specifications

- The glass parts are made of borosilicate glass with ground joints 40/35, and the clamping clip is made of anodized aluminum.
- The vacuum connection is integrated into the body with a porous plate, and the flask mouth is positioned above the filtrate outlet, preventing leaks of filtered liquid into the vacuum tube.
- The porous plate has a diameter of 40 mm, and the filtration funnel is graduated.
- The various parts that make up the equipment can be purchased separately for the replacement of damaged pieces.



Code	Description
PJB001	Complete filtration equipment, 300 mL

Composed of:

Code	Description
PBH030	Erlenmeyer flask ground joint 40/35, 1000 mL
PJD003	Base de vidrio con placa porosa
PJD002	Embudo de filtración, 300 mL
PJD006	Pinza de sujeción de aluminio anodizado