



## Stainless steel vacuum manifold kit - 1 position

Code. PJB021

The connection of the funnel to the filter holder is threaded and therefore does not require an aluminium clamp.

Vacuum manifold kit consisting of:

1. Vacuum manifold with 1 position
2. Stainless steel filter holder
3. Stainless steel funnel 500mL
4. Adapter with stainless steel valve
5. Filters PTFE 5 $\mu$ m

Code	Description
PJB021	Stainless steel vacuum manifold kit, 1 position



## Stainless steel vacuum manifold kit - 3 positions

Code. PJB022

The connection of the funnel to the filter holder is threaded and therefore does not require an aluminium clamp.

Vacuum manifold kit consisting of:

1. Vacuum manifold with 3 positions
2. 3 stainless steel filter holders
3. 3 stainless steel funnels 500mL
4. 3 adapters with stainless steel valve
5. Filters PTFE 5 $\mu$ m

Code	Description
PJB022	Stainless steel vacuum manifold kit, 3 positions



## Stainless steel vacuum manifold kit - 6 positions

Code. PJB023

The connection of the funnel to the filter holder is threaded and therefore does not require an aluminium clamp.

Vacuum manifold kit consisting of:

1. Vacuum manifold with 6 positions
2. 6 stainless steel filter holders
3. 6 stainless steel funnels 500mL
4. 6 adapters with stainless steel valve
5. Filters PTFE 5 $\mu$ m

Code	Description
PJB023	Stainless steel vacuum manifold kit, 6 positions