



Plastic Pasteur pipettes

Code **DLB010, DLB011, DLB022, DLB023, DLB024, DLB025, DLB026, DLB027, DLB028, DLB029, DLB030, DLB031**

- Made of transparent low-density polyethylene (LDPE)
- DNase/RNase, endotoxin and pyrogen free
- Good chemical stability
- Embossed graduation
- Accurate and with good repeatability
- Individually wrapped sterile pipettes. Sterilised with ethylene oxide.



Code	Packaging
DLB010	500 pcs/box
DLB011	500 pcs/box
DLB022	1000 pcs./bag
DLB023	1000 pcs./bag
DLB024	500 pcs./bag
DLB025	500 pcs./bag
DLB026	500 pcs./bag
DLB027	500 pcs./bag
DLB028	100 pcs./bag
DLB029	100 pcs./bag
DLB030	100 pcs./bag
DLB031	100 pcs./bag

Technical specifications

Code	DLB010	DLB011	DLB022	DLB023	DLB024	DLB025	DLB026	DLB027	DLB028	DLB029	DLB030	DLB031
Nominal volume (mL)	1	3	0,2	0,5	1	2	3	3 (long)	5	10	1	3
Total volume (mL)	5	7	0,9	1,8	3	5	6	8,2	12	23	5	7
Drop volume (μL)	45	48	20	35	42	48	48	48	40	55	45	48
Total length (mm)	145	155	65	126	145	155	155	180	210	300	145	155
Sterile	No	No	No	No	No	No	No	No	No	No	Yes	Yes